



V.A.C.®
Therapy

3M™ V.A.C.® Simplicity Therapy System

3M™ V.A.C.® Therapy in the Simplest Form

“One Touch” functionality provides the easiest
3M™ V.A.C.® Therapy System for patient operation



KCI -125 mmHg



Keeping it simple for both patients and clinicians

From the makers of 3M™ V.A.C.® Therapy, comes a pre-set system designed to deliver the right therapy, with fewer steps. 3M™ V.A.C.® Simplicity Therapy System offers the benefits and safety features of the most clinically proven NPWT on the market—with a single touch.

Easy-to-use and pre-loaded with the features you'll find only on 3M™ V.A.C.® Therapy

Only V.A.C.® Therapy features proprietary technology to reduce erroneous alarms, and if an alarm does occur, helps to resolve issues.

3M™ SensaT.R.A.C.™ Technology continually monitors and regulates -125 mm Hg negative pressure at the wound site

- Systems without regulated pressure feedback at the wound site are limited because they cannot sense or measure pressure delivered at the wound¹


3M™ Easyclear Purge™ Technology is proactive against blockages

- Helps to clear fluid in the lumen tubing to aid in the prevention and clearance of blockages¹

Intuitive alerts and alarms with minimal button pushes

- Low battery
- Blockage
- Leak

3M™ Easyclear Purge™
Technology Lumen



3M™ SensaT.R.A.C.™ Pad



Easy to Operate:

- Simple interface minimizes device training time
- Automatically delivers the most commonly prescribed NPWT setting of continuous -125 mm Hg¹
- Quick release canisters allow for efficient replacement
- 300ml canister minimizes canister changes
- Long battery life (8+ hours) with fast recharge time

Built for patient comfort:

- Lightweight and portable
- Carrying case allows discreet delivery of therapy on-the-go
- Alarms help provide patient safety
- Multi-position canister for mobility

Minimal training time and seamless transition

Easy: Simple and intuitive interface with no settings to learn

Compatible: 3M™ V.A.C.® Dressings can be used with either the 3M™ V.A.C.® Simplicity Therapy System or 3M™ ActiV.A.C.™ Therapy System

Efficient: Dressings changes up to 3 days

The right fit for every patient, from the 3M™ V.A.C.® Therapy you trust



Feature	3M™ V.A.C.® Simplicity Therapy System	3M™ ActiV.A.C.™ Therapy System
300 ml Canister Size	✓	✓
Customizable Dressings	✓	✓
Blockage Alarm	✓	✓
Leak Alarm	✓	✓
Battery Alarm	✓	✓
3M™ SensaT.R.A.C.™ Technology	✓	✓
Battery Life	8hrs.	14hrs.
6-Hour Charge Time	✓	✓
Light Weight and Portable	2.3lbs.	2.4lbs.
Therapy Option	Continuous	Continuous or Intermittent
Pressure Setting	-125 mm Hg	-25 mm Hg - 200 mm Hg

A world leader in skin and wound care right by your side.

As your partner, we're here to help you help your patients. When we combine our science with your expertise, amazing things happen.



Science-Based Solutions

3M products are trusted in more than 60,000 hospitals, homes and businesses worldwide. Our comprehensive portfolio of advanced wound care solutions is supported by clinical evidence across new and growing categories—including dressings, disposables, and digital technology and connectivity.



Ongoing Support

From product ordering to patient treatment and discharge, our support is seamless and efficient. Getting started is easy with online ordering and scheduled same-day delivery, and our team is available 24/7 to assist you with onsite clinical and technical support.



World-Class Education

We act as an extension of your team—empowering you with hands-on training and award-winning education. Access free medical education live and on-demand at: **3M.com/MedicalEducation**

Ready to take the next step?



Contact Your 3M
Account Representative



Visit [3M.com/KCI](https://www.3m.com/KCI)



Call **1-800-275-4524**

Available 24/7, including
Clinical and Technical support



3M Company
2510 Conway Ave
St. Paul, MN 55144 USA

Phone 1-800-275-4524 (NPWT products)
1-800-228-3957
Web [3M.com/medical](https://www.3m.com/medical)

Reference:

1. KCI. Cumulative NPWT Wounds - 10 million_ 2013-2015_Internal Report (v1.0). 2018

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.

© 2021 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.
70-2011-8246-9 PRA-PM-US-03106 (04/21)