

V.A.C.® SIMPLICITY™

THERAPY SYSTEM

V.A.C.® Therapy in the Simplest Form

“One Touch” functionality provides the easiest V.A.C.® Therapy System for patient operation



3M + **KCI**™



Keeping it simple for both patients and clinicians

From the makers of V.A.C.® Therapy, comes a pre-set system designed to deliver the right therapy, with fewer steps. V.A.C. SIMPLICITY™ Therapy System offers the benefits and safety features of the most clinically proven NPWT on the market—with a single touch.

Easy-to-use and pre-loaded with the features you'll find only on V.A.C.® Therapy

Only V.A.C.® Therapy features proprietary technology to reduce erroneous alarms, and if an alarm does occur, helps to resolve issues.

SENSAT.R.A.C.™ Technology continually monitors and regulates -125mmHg negative pressure at the wound site

- Systems without regulated pressure feedback at the wound site are limited because they cannot sense or measure pressure delivered at the wound¹

EASYCLEAR PURGE™ Technology is proactive against blockages

- Helps to clear fluid in the lumen tubing to aid in the prevention and clearance of blockages¹

Intuitive alerts and alarms with minimal button pushes

- Low battery
- Blockage
- Leak

EASYCLEAR PURGE™ TECHNOLOGY LUMEN



SENSAT.R.A.C.™ PAD



Easy to Operate:

- Simple interface minimizes device training time
- Automatically delivers the most commonly prescribed NPWT setting of continuous -125mmHg¹
- Quick release canisters allow for efficient replacement
- 300ml canister minimizes canister changes
- Long battery life (8+ hours) with fast recharge time

Built for patient comfort:

- Lightweight and portable
- Carrying case allows discreet delivery of therapy on-the-go
- Alarms help provide patient safety
- Multi-position canister for mobility

Minimal training time and seamless transition

Easy: Simple and intuitive interface with no settings to learn

Compatible: V.A.C.® Dressings can be used with either the V.A.C. SIMPLICITY™ System or ACTIV.A.C.™ System

Efficient: Dressings changes up to 3 days

The right fit for every patient, from the V.A.C.® Therapy you trust



Feature	V.A.C. SIMPLICITY™ System	ACTIV.A.C.™ Therapy System
300 ml Canister Size	✓	✓
Customizable Dressings	✓	✓
Blockage Alarm	✓	✓
Leak Alarm	✓	✓
Battery Alarm	✓	✓
SENSAT.R.A.C.™ Technology	✓	✓
Battery Life	8hrs.	14hrs.
6-Hour Charge Time	✓	✓
Light Weight and Portable	2.3lbs.	2.4lbs.
Therapy Option	Continuous	Continuous or Intermittent
Pressure Setting	-125mmHg	Range of Settings

Your partner, supporting the business of healing

10 Million

wounds worldwide Have been
treated with V.A.C.® Therapy¹

400

clinically trained representatives
support therapy outcomes



>90%

of published NPWT clinical evidence
is based on V.A.C.® Therapy²

65

dedicated team of nurses
monitor patient Progress

KCI National Contact Center

24/7 assistance with ordering, delivery,
billing, clinical, and technical support
800-275-4524

KCI Reimbursement Hotline

Specialists assist with insurance
coding, coverage guidelines, and other
reimbursement information.

E-mail: ReimbursementEducation@kci1.com

Fax: 800-311-9291

M-F: 8:00am–8:00pm EST

To learn more about the V.A.C. SIMPLICITY™ Therapy System,
call **800-275-4524** or visit mykci.com

Reference:

1. KCI. Cumulative NPWT Wounds. 2018

2. KCI. Percentage of V.A.C.(R) Therapy Articles vs Comp Articles. 9/20/2017

NOTE: Specific indications, contraindications, warnings, precautions, and safety information exist for these products and therapies. Please consult a physician and product instructions for use prior to application. Rx only.

© 2020 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. LIT#29-A-440 • PRA-PM-US-00150 (03/20)

