

Wound management at the touch of a button

Easy to operate

- Simple interface minimizes device training time
- Automatically delivers the most commonly prescribed NPWT setting of Continuous -125 mmHg
- SENSAT.R.A.C.™ Dressing/Technology continually monitors and regulates Management pressure at the wound site
- Quick Release canisters allow for efficient replacement
- Multi-orientational canister allows for patient mobility and product position flexibility
- 300ml canister minimizes canister changes
- Long battery life (8+ hours) with fast recharge time

Patient comfort

- Lightweight and portable
- Carrying case allows discreet delivery of therapy on-the-go
- Alarms help provide patient safety

Specifications:

V.A.C. SIMPLICITY™ Therapy unit: P/N 414063

Dimensions:7.6" W (19.3cm) x 6" H (15.2cm) x 2.5" D (6.4cm)

Weight:2.3 lbs. (1.04 kg)

Canister Volume:300ml

Battery Life:.....8+ hours (recharge time <6 hours)

Electrical Charger/ External Supply Input:100-240 VAC .72 A @ 115 VAC 50/60 Hz

Charger/External Supply Output:....12 V, 3.3 A

Maximum Electrical Leakage:.....<100 Microamps

IEC Classification: Type BF, Applied Part Class II IP22

Specifications are subject to change without notice.

Therapy options:

Negative Pressure:.....-125 mmHg

Mode of Operation:Continuous

Intensity:.....High

Storage conditions

Temperature Range:.....-13(o)F (-25(o)C) to 158(o)F (70(o)C)

Relative Humidity Range:0-93% non-condensing

Operating conditions

Temperature Range:.....41°F (5°C) to 104°F (40°C)

Relative Humidity Range:15-93% non-condensing



NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for KCI products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.