

3M Science.
Applied to Life.™

3M™ RelyX™ Universal
Resin Cement

Goodbye
complexity.



The two that can do it all.

3M™ RelyX™ Universal Resin Cement works both as a standalone, self-adhesive cement and as an adhesive cement when combined with 3M™ Scotchbond™ Universal Plus Adhesive.

Together they eliminate the hassle and confusion of multiple resin cements, primers and adhesives, **simplifying your direct and indirect restorative workflows.** With fewer products to stock and replenish, you save space, time and money.

So simple yet so powerful.

- Two components for virtually all indirect and direct bonding indications
- Outstanding self-adhesion of RelyX Universal Resin Cement
- Scotchbond Universal Plus Adhesive further enhances bond strength to all substrates
- Unique 3M automix syringe
- Easy excess clean-up
- High aesthetics



Game-changing syringe.

Proprietary 3M design takes automix cement applications to the next level. Developed with your needs in mind for easy dispensing, cleaning and storage. The small, low-waste 3.4 g syringe provides enough material for approximately 15 to 20 applications.*

Easy excess clean-up.

Efficient to the last step. Its unique chemistry ensures excess outflow stays put at the restoration margins, and can be easily removed after tack curing.



Self-adhesive cementation of a 3M™ Lava™ Esthetic Fluorescent Full-Contour Zirconia bridge; Image courtesy of Dr. Jan Gueth, Germany

Less stock, more savings

95% of dentists agreed the universal character of RelyX Universal Resin Cement helps to **reduce stock and save costs.**

Easier for operatory staff

97% of dentists agreed using one universal versus multiple resin cements is **easier for operatory staff.**

* Based upon the amount of cement used in a restoration from 60 µL (e.g. inlay, onlay, veneer) to 90 µL (e.g. crown). ** Per application compared to currently available standard automix systems. Source: M. Cowen, J.M. Powers, Dental Advisor, Biomaterials Research Report, Number of Automix Applications and Mixing Efficiency, No. 134 – June 23, 2020. *** Source: Field Evaluation EU/USA conducted by 3M.

80%

less cement waste** with Micro Mixing Tip



More ergonomic shape with **50%** less plastic waste**

Hygienic self-sealing syringe – for storage without the used mixing tip.



Thin, long elongation tip for easy root canal application



The strong partner for adhesive cases.



3M™ Scotchbond™ Universal Plus Adhesive is the first radiopaque, all-in-one universal adhesive for all direct and indirect bonding indications. And, it is a fully aligned system: the adhesive is cured by the cement, so no light-cure is needed.

Scotchbond Universal Plus Adhesive serves as universal primer for all materials including glass ceramics and further boosts the cement's excellent self-adhesive bond e.g. to zirconia.



Outstanding self-adhesion for maximised efficiency.


With its excellent self-adhesive bond strength to dentine, RelyX Universal Resin Cement reliably treats most cases with no need for an additional adhesive or primer. For light- or self-cure.



Virtually no post-operative sensitivity

99% of dentists reported that RelyX Universal Resin Cement shows virtually no post-operative sensitivity***

High aesthetics with four colour-stable, fluorescent shades



Translucent (TR) A1 A3 Opaque (A3O) White Opaque (WO)

Ordering information

3M™ RelyX™ Universal Resin Cement

Item no. Product information

Intro Kit

56968 4 x 3M™ RelyX™ Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60 x 3M™ RelyX™ Universal Micro Mixing Tips, 20 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (5 ml), 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

Trial Kits

Each trial kit contains: 1 x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15 x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (1.5 ml), 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides

56969 Shade TR **56970** Shade A1

Refills

Each refill contains: 1 x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15 x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips

56971 Shade TR **56973** Shade A3O

56972 Shade A1 **56974** Shade WO

Value Packs

Each value pack contains: 3 x 3M™ RelyX™ Universal Resin Cement syringes (3.4g each), 45 x 3M™ RelyX™ Universal Micro Mixing Tips, 15 Elongation Tips

56977 Shade TR **56978** Shade A1

Mixing Tips

56975 30 x 3M™ RelyX™ Universal Micro Mixing Tips

56976 15 x 3M™ RelyX™ Universal Micro Mixing Tips with Elongation Tips

3M™ Scotchbond™ Universal Plus Adhesive

Vial

41293 **Intro Kit:** 1 x 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml), 50 disposable brushes, 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 dispensing tips

Refill: 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml)

41294 1x 5 ml vial **41295** 3x 5 ml vial

Unit Dose

3M™ Scotchbond™ Universal Plus Adhesive Unit Dose L-Pops, 1 x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 dispensing tips

41297 **Intro Kit:** 50x L-Pop™ **41299** 200x L-Pop™

41298 100 x L-Pop™

3M™ Scotchbond™ Universal Etchant Etching Gel

41263 2 x 3M™ Scotchbond™ Universal Etchant (3 ml each), 50 dispensing tips



3M™ RelyX™ Universal Resin Cement Intro Kit (56968)



3M™ RelyX™ Universal Resin Cement Trial Kit TR (56969)



3M™ RelyX™ Universal Resin Cement Refill TR (56971)



3M™ Scotchbond™ Universal Plus Adhesive Refill (41294)



3M Australia Pty Limited
Building A, 1 Rivett Rd
North Ryde, NSW 2113
1300 363 454
3M.com.au/RelyXUniversal

3M New Zealand Limited
94 Apollo Drive
Rosedale, Auckland 0632
0800 808 182
3M.co.nz/RelyXUniversal

Scan to learn more



Always follow the Instructions for Use (IFU) and refer to IFU for full indications, precautions and warnings. Claims supported by 3M data on file.

3M, "3M Science. Applied to Life.", Lava, RelyX and Scotchbond are trademarks of 3M. All other trademarks are owned by other companies. © 3M 2020. All rights reserved.