





## 3M™ Attest™ Vaporized Hydrogen Peroxide Tri-Metric Chemical Indicator, Type 4

### Unprocessed



3M™ Attest™ Vaporized Hydrogen Peroxide Tri-Metric Chemical Indicators, Type 4 comply with ISO 11140-1:2014.

### Processed



ACCEPT



ACCEPT



ACCEPT



ACCEPT

Sample Comparison

**If a color change to pink has entered the ACCEPT zone it is considered a pass.**

### Processed



REJECT



REJECT



REJECT

Sample Comparison

**If a color change to pink has not entered the ACCEPT zone, it is considered a fail.**

Please reference the instructions for use for additional information.