



**Choose the right stethoscope to help you perform at your best.**

3M™ Littmann® Stethoscopes feature adjustable headsets with soft-sealing eartips. Select Littmann mechanical stethoscopes are designed with tunable technology to give you the ability to hear different frequency sounds by simply adjusting the pressure applied to the chestpiece. Get the acoustic performance, high quality and comfort you need with a Littmann stethoscope.



	Lightweight II S.E. Stethoscope	Classic II Paediatric and Infant Stethoscopes	Master Classic II™ Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	Electronic Stethoscope Model 3200
	Ideal for taking blood pressure readings and making limited patient assessments.	Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.			Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.		High-performance auscultation. Hear more of what you need to hear with amplification of 24x, ambient noise reduction and optional connectivity.
<b>Approximate Length</b>	71cm	71cm	69cm	69cm	69cm	69cm	69cm
<b>Approximate Weight</b>	118g	105g, 95g	160g	150g	177g	185g	185g
<b>Warranty</b>	2 year	3 year	3 year	5 year	7 year	7 year	2 year
<b>Diaphragm</b>	Tunable	Non-Tunable (Floating)	Tunable	Tunable (2)	Tunable (2)	Tunable	Electronic
<b>Approximate Chestpiece Diameter</b>	Double-sided 5.2 x 4.5cm 3.5 x 3.2cm	Double-sided Paediatric 3.4cm 2.7cm Infant 2.8cm 1.7cm	Single-sided 4.3cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm	Double-sided Adult side 4.3cm Paediatric side 3.4cm	Single-sided 4.8cm	Single-sided 5.1cm
<b>Features</b>	Teardrop-shaped chestpiece Fits under blood pressure cuff	Specially designed for the smallest patients	Designed for easier grip and maneuvering	Convertible open bell	Convertible open bell Dual-lumen tubing	Dual-lumen tubing Best acoustic response in Littmann mechanical line	Eliminates 85% (on average) of ambient noise Bluetooth® wireless technology Recording and playback



**The power to learn, in the palm of your hand.**

3M™ Littmann® Learning Institute App – Free with purchase of a 3M™ Littmann® Stethoscope!

**Need more help? Visit [Littmann.co.uk](http://Littmann.co.uk)**



# 3M™ Littmann® Stethoscope Ordering Information and Tube Colour Guide

**Components**  
(Chestpiece/Eartube/Stem)

Lightweight II SE Stethoscope	Classic II Infant Stethoscope	Classic II Paediatric Stethoscope	Master Classic II™ Stethoscope	Classic III™ Monitoring Stethoscope	Cardiology IV™ Diagnostic Stethoscope	Master Cardiology™ Stethoscope	Electronic Model 3200 Stethoscope
-------------------------------	-------------------------------	-----------------------------------	--------------------------------	-------------------------------------	---------------------------------------	--------------------------------	-----------------------------------

## Tube Color

All Standard Finish	2450 Black	2114 Black	2113 Black	2144L Black	5620 Black	5627 Burgundy	6152 Black	6156 Plum	2160 Black	3200BK27 Black
	2451 Burgundy	2114R Red	2113R Red	2146 Burgundy	5621 Gray	5633 Pearl Pink	6154 Navy	6158 Raspberry	2163 Burgundy	3200BU Burgundy
	2452 Caribbean	2124 Caribbean	2119 Caribbean	2147 Navy	5622 Navy	5831 Plum	6155 Hunter Green		2164 Navy	3200NB Navy
	2454 Ceil Blue		2122 Raspberry		5623 Caribbean	5832 Lavender			2167 Plum	
	2456 Pearl Pink				5630 Ceil Blue	5835 Turquoise			2178 Caribbean	
						5839 Lemon-Lime				
All Black				2141 Black	5803 Black	5868 Burgundy	6163 Black	6168 Navy	2161 Black	
All Brass Finish									2175 Black	
All Copper Finish					5809 Chocolate					
All Smoke					5811 Black	5871 Raspberry	6162 Black		2176 Black	
					5812 Dark Olive		6166 Plum	6171 Turquoise	2182 Dark Olive	
Rainbow/Brass/Brass	2157 Raspberry	2153 Caribbean			5806 Raspberry	5807 Caribbean				
Brass/Smoke/Smoke							6164 Black			
Rainbow/Black/Black					5870 Black		6165 Black			
Mirror/Smoke/Smoke					5862 Raspberry	5863 Navy				
					5959 Ceil Blue					
Mirror/Standard/Mirror							6177 Black	6170 Burgundy		
Champagne/Smoke/Smoke					5861 Black	5864 Burgundy	6179 Black			
Champagne/Black/Champagne							6176 Burgundy	6190 Caribbean		
High-Polish Rainbow/Black/Black							6242 Navy			
High-Polish Rainbow/Smoke/Smoke							6240 Black	6241 Raspberry		
<b>COLOURED STEM SERIES: (square inset represents stem colour)</b>										
Black/Black/Colour							6200 Black	6203 Black		
							6201 Black			
High-Polish Smoke/Black/Colour							6202 Navy	6204 Black		
							6232 Black	6234 Caribbean		
High-Polish Smoke/Smoke/Colour							6238 Gray			
Rainbow/Black/Colour							6205 Plum			
Champagne/Champagne/Colour							6206 Hunter			
Smoke/Smoke/Colour					5872 Turquoise	5875 Lime				
					5873 Gray					
Mirror/Smoke/Colour					5960 Plum	5962 Pearl Pink				
Mirror/Standard/Colour					5874 Caribbean					

**Components:**  
Chestpiece/Eartube/Stem

Standard Finish



Black



Brass Finish




Copper Finish




Smoke



Rainbow




Mirror




Champagne



High-Polish Smoke



High-Polish Rainbow



Colours shown are approximate.

[www.littmann.co.uk](http://www.littmann.co.uk) | Also available: Veterinary Stethoscope (1392), Teaching Stethoscope (2138)

3M, Littmann, the L Littmann logo, Classic III, Cardiology IV, Master Classic II and Master Cardiology are trademarks of 3M. Bluetooth is a registered trademark of Bluetooth SIG. © 3M 2020. All rights reserved.