

# Too much on your tray?

Check out 3M's universal products. They're all you need to perform virtually any restoration.



# One system.

## Virtually all bonding techniques.



Create positive patient experiences

What if you could cut down inventory to just the essentials?  
The universal series from 3M allows you to perform virtually any restoration with fewer products, helping you and your practice.



Simplify direct and indirect workflows



Save time



Worry less about cost and waste

## One family of products. A universe of possibilities.

### 3M™ RelyX™ Universal Resin Cement

Works both as a standalone, self-adhesive cement and as an adhesive cement when combined with 3M™ Scotchbond™ Universal Plus Adhesive.



### 3M™ Scotchbond™ Universal Plus Adhesive

A true one-bottle system<sup>1</sup> with gold standard bond strength, suitable for all direct and indirect indications, all etching techniques, and all dental substrates. It also serves as a universal primer for virtually all restoration materials.



### 3M™ Filtek™ Universal Restoratives

Simplify your single-shade restorations while providing natural-looking esthetics with top-of-the-line polish retention and wear resistance.



### 3M™ Filtek™ Supreme Ultra Universal Restorative

Exceptional esthetics and handling with excellent polish retention and outstanding strength for both anterior and posterior restorations.



<sup>1</sup>No separate primer, silane, or dual-cure activator required.

# A cement that can do it all. **Unparalleled simplicity.**

Eliminate the hassle and confusion of stocking multiple resin cements, primers and adhesives. **3M™ RelyX™ Universal Resin Cement** works both as a standalone, self adhesive cement and as an adhesive cement with 3M™ Scotchbond™ Universal Plus.

**99%**  
clinical  
rating\*



- Virtually all indirect indications with self-adhesive mode or as a two-component system with 3M™ Scotchbond™ Universal Plus Adhesive
- Outstanding self-adhesion for maximum efficiency
- 80% less cement waste with micro mixing tip<sup>2</sup>
- Virtually no post-operative sensitivity<sup>3</sup>
- Higher radiopacity than enamel
- Easy excess clean-up



Image courtesy of  
Dr. Jan Gueth,  
Germany

The **Micro Mixing Tip** is removed right after use enabling hygienic storage without a used mixing tip. The syringe is cleanly sealed by a unique valve mechanism.

\*Dental Advisor Editors' Choice - 99% clinical evaluation rating

<sup>2</sup>Per application compared to currently available standard automix systems. <sup>3</sup>Source: Field Evaluation EU/USA conducted by 3M.

# The universal adhesive that set the bar ...now raises it.

## 3M™ Scotchbond™ Universal Plus Adhesive

### First radiopaque universal adhesive



Non-radiopaque  
adhesive



3M™ Scotchbond™  
Universal Plus Adhesive

- 💧 Bonds and seals caries-affected dentin
- 💧 Gold standard adhesion to virtually all dental substrates, including glass ceramics
- 💧 Universally suitable for all direct and indirect indications and all etching techniques
- 💧 True one-bottle system with no need for a separate primer, silane or dual-cure activator
- 💧 Virtually no post-operative sensitivity



Also available in a unit dose for efficient hygiene management —  
reduces the risk of cross-contamination.



**100%**  
recommend to  
a colleague\*

\*Dental Advisor Editors' Choice - 100% recommend to a colleague




# Universal restoratives **for every procedure.**

Designed to cover every class of restoration, **3M™ Filtek™ Universal Restorative** and **3M™ Filtek™ Supreme Ultra Restorative** have you covered for all of your clinical needs. Discover why Filtek™ restoratives have been trusted in over 1 billion restorations worldwide.

## Cover dark stains

Versatile pink opaquer specifically designed to mask metal and dark tooth structure



**96%**  
clinical  
rating\*



**96%**  
overall  
rating\*\*



## 3M™ Filtek™ Supreme Ultra Universal Restorative

- 36 shades, 4 opacities
- Exceptional esthetics
- Natural looking fluorescence
- Outstanding strength for anterior and posterior
- 15 years of clinical success

## Both restoratives

- Can be used for every class of restoration (Classes I-V)
- TRUE nanotechnology
- Excellent wear resistance, polish retention and handling
- Minimally invasive
- Can be safely warmed

## 3M™ Filtek™ Universal Restorative

- 8 designer shades
- Universal opacity
- Extra white
- Pink opaquer
- Predictable, repeatable esthetics
- Low shrinkage stress/stress relief

# Ordering Information

## 3M™ RelyX™ Universal Resin Cement

Item No.	Product Information
<b>Intro Kit</b>	
<b>56968</b>	4x 3M™ RelyX™ Universal Resin Cement syringes shades TR, A1, A3O, WO (3.4g each), 60x 3M™ RelyX™ Universal Micro Mixing Tips, 20 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides
<b>Trial Kits</b>	
<b>56969</b>	<b>Shade TR</b>
<b>56970</b>	<b>Shade A1</b>
	Each trial kit contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips, 1 vial 3M™ Scotchbond™ Universal Plus Adhesive (1.5 ml), 1x 3M™ Scotchbond™ Universal Etchant (3 ml), 25 blue etchant tips, 50 orange microbrushes, step-by-step guides
<b>Refills</b>	
<b>56971</b>	<b>Shade TR</b>
<b>56972</b>	<b>Shade A1</b>
<b>56973</b>	<b>Shade A3O</b>
<b>56974</b>	<b>Shade WO</b>
	Each refill contains: 1x 3M™ RelyX™ Universal Resin Cement syringe (3.4g), 15x 3M™ RelyX™ Universal Micro Mixing Tips, 5 Elongation Tips
<b>Value Packs</b>	
<b>56977</b>	<b>Shade TR</b>
<b>56978</b>	<b>Shade A1</b>
	Each value pack contains: 3x 3M™ RelyX™ Universal Resin Cement syringes (3.4g each), 45x 3M™ RelyX™ Universal Micro Mixing Tips, 15 Elongation Tips



3M Oral Care  
2510 Conway Avenue  
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249  
Web 3M.com/dental

## 3M™ Scotchbond™ Universal Plus Adhesive

Item No.	Product Information
<b>Vial</b>	
<b>41293</b>	<b>Intro Kit:</b> 1x 3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml), 50 disposable brushes, 1x 3M™ Scotchbond™ Universal Etchant (3ml), 25 dispensing tips
<b>Vial Refills</b>	
<b>41294</b>	<b>1x 5 ml vial</b>
<b>41295</b>	<b>3x 5 ml vial</b>
	3M™ Scotchbond™ Universal Plus Adhesive vial (5 ml)
<b>Unit Dose</b>	
<b>41297</b>	<b>Intro Kit: 50x L-Pop</b>
<b>41298</b>	<b>100x L-Pop</b>
<b>41299</b>	<b>200x L-Pop</b>
<b>41304</b>	<b>400x L-Pop</b>
	3M™ Scotchbond™ Universal Plus Adhesive Unit Dose L-Pops, 1x 3M™ Scotchbond™ Universal Etchant (3ml), 25 dispensing tips

## 3M™ Scotchbond™ Universal Etchant Etching Gel

Item No.	Product Information
<b>41263</b>	2x 3M™ Scotchbond™ Universal Etchant (3ml each), 50 dispensing tips

## 3M™ Filtek™ Universal Restorative

**Shade Refills:**  
**Capsules** include 20 – 0.2g capsules; **Syringes** include 1 – 4g syringe  
**Pink Opaquer (PO) Refills:**  
**Capsules** include 10 – 0.2g capsules; **Syringes** include 1 – 4g syringe  
All refills include a Technique Guide and Instructions for Use

	Shade	Syringe	Capsule
<b>Refills</b>			
	<b>A1</b>	6555A1	6550A1
	<b>A2</b>	6555A2	6550A2
	<b>A3</b>	6555A3	6550A3
	<b>A3.5</b>	6555A3.5	6550A3.5
	<b>A4</b>	6555A4	6550A4
	<b>B1</b>	6555B1	6550B1
	<b>B2</b>	6555B2	6550B2
	<b>D3</b>	6555D3	6550D3
	<b>XW—Extra White</b>	6555XW	6550XW
	<b>PO—Pink Opaquer</b>	6555PO	6550PO

## 3M™ Filtek™ Supreme Ultra Universal Restorative

**Syringe Shade Refills**—include 1 – 4g syringe; Instructions for Use  
**Capsule Shade Refills**—Dentin, Enamel and Translucent Shade Refills include 10 – 0.2g capsules; Body Shade Refills include 20 – 0.2g capsules; Instructions for Use

	Shade	Syringe	Capsule
<b>Body</b>			
	<b>A1B</b>	6028A1B	6029A1B
	<b>A2B</b>	6028A2B	6029A2B
	<b>A3B</b>	6028A3B	6029A3B
	<b>A3.5B</b>	6028A3.5B	6029A3.5B
	<b>A4B</b>	6028A4B	6029A4B
	<b>A6B</b>	6028A6B	6029A6B
	<b>B1B</b>	6028B1B	6029B1B
	<b>B2B</b>	6028B2B	6029B2B
	<b>B3B</b>	6028B3B	6029B3B
	<b>B5B</b>	6028B5B	6029B5B
	<b>C1B</b>	6028C1B	6029C1B
	<b>C2B</b>	6028C2B	6029C2B
	<b>C3B</b>	6028C3B	6029C3B
	<b>D2B</b>	6028D2B	6029D2B
	<b>D3B</b>	6028D3B	6029D3B
	<b>WB</b>	6028WB	6029WB
	<b>XWB</b>	6028XWB	6029XWB
<b>Dentin</b>			
	<b>A1D</b>	6028A1D	6029A1D
	<b>A2D</b>	6028A2D	6029A2D
	<b>A3D</b>	6028A3D	6029A3D
	<b>A4D</b>	6028A4D	6029A4D
	<b>B3D</b>	6028B3D	6029B3D
	<b>C4D</b>	6028C4D	6029C4D
	<b>WD</b>	6028WD	6029WD
<b>Enamel</b>			
	<b>A1E</b>	6028A1E	6029A1E
	<b>A2E</b>	6028A2E	6029A2E
	<b>A3E</b>	6028A3E	6029A3E
	<b>B1E</b>	6028B1E	6029B1E
	<b>B2E</b>	6028B2E	6029B2E
	<b>D2E</b>	6028D2E	6029D2E
	<b>WE</b>	6028WE	6029WE
	<b>XWE</b>	6028XWE	6029XWE
<b>Translucent</b>			
	<b>Clear</b>	6028CT	6029CT
	<b>Blue</b>	6028BT	6029BT
	<b>Grey</b>	6028GT	6029GT
	<b>Amber</b>	6028AT	6029AT

3M, Filtek, RelyX and Scotchbond are trademarks of 3M or 3M Deutschland GmbH. Printed in USA.  
© 3M 2022. All rights reserved. 70-2013-7334-0 Rev Apr 2022