



Dermatac™

Drape

An innovative mix to support the next level of healing

Apply with ease: conforms to different anatomical locations and can be repositioned at initial placement

Seal in the heal: maintains a highly-effective seal for negative pressure

Remove with kindness: skin-friendly removal supports optimal healing and patient comfort

Specifications

Drape size:

- 24.9cm x 21cm
(overall dimensions with handle bars)
- 19.7cm x 21cm
(after handle bar removal)

Environmental/storing conditions

- Store at room temperature

Packaging dimensions

- 3M™ Dermatac Drape, 10-pack, 12.63" x 10.38" x 2"
(32.07cm x 26.35cm x 5.08cm)

3M™ Dermatac Drape and V.A.C.® Granufoam™ Dressing Kit

- 5-packs: 11.92" x 10.74" x 12.99"
(30.28cm x 27.28cm x 32.99cm)
- 10-packs: 11.92" x 12.99" x 21.10"
(30.28cm x 32.99cm x 53.59cm)



| Product Details | Drape(s) included - per kit | SKU |
|------------------------------------------------------------------------|-----------------------------|-----------|
| Drape only | | |
| 3M™ Dermatac™ Drape 10-Pk | 10 | DTAC10LDP |
| 3M™ V.A.C.® Granufoam™ Dressings | | |
| 3M™ Dermatac™ Drape and V.A.C.® Granufoam™ Dressing Kit – Small 5-Pk | 1 | DTGF05PKS |
| 3M™ Dermatac™ Drape and V.A.C.® Granufoam™ Dressing Kit – Small 10-Pk | 1 | DTGF10PKS |
| 3M™ Dermatac™ Drape and V.A.C.® Granufoam™ Dressing Kit – Medium 5-Pk | 1 | DTGF05PKM |
| 3M™ Dermatac™ Drape and V.A.C.® Granufoam™ Dressing Kit – Medium 10-Pk | 1 | DTGF10PKM |
| 3M™ Dermatac™ Drape and V.A.C.® Granufoam™ Dressing Kit – Large 5-Pk | 2 | DTGF05PKL |
| 3M™ Dermatac™ Drape and V.A.C.® Granufoam™ Dressing Kit – Large 10-Pk | 2 | DTGF10PKL |

For more information contact your local representative.

3M Company
2510 Conway Ave
St. Paul, MN 55144 USA

Phone 1-800-275-4524 (NPWT products)
1-800-228-3957
Web 3M.com/medical

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.

© 2021 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. 70-2011-8207-1 PRA-PM-US-03173 (05/21)