

A full line of restorative solutions to fit your unique needs.

With a rich legacy of oral care products and services, 3M is one of the world's most trusted companies. Filtek™ restoratives have been used in more than 1 billion restorations worldwide. Our family of Filtek™ dental restoratives lets you choose just the right option to fit any patient situation—from highly esthetic, to highly simple for an easier day.

Versatility & Esthetics

3M™ Filtek™ Z350XT Universal Restorative

- Exceptional esthetics
- Natural-looking fluorescence
- Excellent polish retention
- Exceptional handling
- Outstanding strength for anterior and posterior
- Wear resistance similar to enamel¹
- Can be safely warmed² and injection molded
- TRUE nanotechnology
- 36 shades, 4 opacities



3M™ Filtek™ Supreme Flowable Restorative

- Virtually no bubbles or material “run-on”
- Easy-to-bend cannula for improved access
- Ergonomic design is easy to hold and inject
- Excellent adaptation, polish retention and wear resistance
- Can be safely warmed³ and injection molded
- 7 shades



Efficiency & Simplicity

3M™ Filtek™ One Bulk Fill Restorative

- Increased opacity for improved esthetics
- One-step placement up to 5 mm
- Excellent handling
- Excellent adaptation
- Stress relief
- Excellent wear resistance
- Can be safely warmed² and injection molded
- TRUE nanotechnology
- 5 shades



3M™ Filtek™ Bulk Fill Flowable Restorative

- Virtually no bubbles or material “run-on”
- Easy-to-bend cannula for improved access
- Ergonomic design is easy to hold and inject
- Injectable, one-step placement up to 4mm
- Excellent adaptation, polish retention and low shrinkage
- Can be safely warmed³ and injection molded
- 4 shades






¹ In vitro data on file.

² Warm capsules to 70°C (158°F) for up to 1 hour.

³ Up to 70°C (158°F) for up to 1 hour, up to 25 times.

⁴ 3M internal data.

We have you covered for every class of direct restoration.

	Versatility & Esthetics		Efficiency & Simplicity	
	3M™ Filtek™ Supreme Ultra Universal Restorative (Incremental up to 2 mm)	3M™ Filtek™ Supreme Flowable Restorative (Incremental up to 2 mm)	3M™ Filtek™ One Bulk Fill Restorative (Bulk up to 5 mm)	3M™ Filtek™ Bulk Fill Flowable Restorative (Incremental up to 4 mm)
Class I 	✓	✓ (Liner/Base)	✓ (up to 4 mm)	✓ (Liner/Base)
Class II 	✓	✓ (Liner/Base)	✓ (up to 5 mm)	✓ (Liner/Base)
Class III 	✓	✓		✓
Class IV 	✓			
Class V 	✓	✓		✓
Minimally Invasive	✓	✓	✓	✓
Stress Bearing	✓		✓	

✓ = recommended use

Refer to Instructions for Use for full list of indications.



www.3M.com.my

3M Malaysia Sdn Bhd (7251 V)
 Level 8, Block F, Oasis Square
 No. 2, Jalan PJU 1A/7A, Ara Damansara
 47301 Petaling Jaya, Selangor, Malaysia
 Tel: +603 7884 2888
 Fax: +603 7884 2902
www.3M.com.my

3M and Filtek are trademarks of 3M.