

3M Science.
Applied to Life.™



VH202 3M™ Attest™ Vaporized
Hydrogen Peroxide Tri-Metric
Chemical Indicator 1348/1348E

Monitor three key variables with one simple innovation.

The first VH202 CI to independently monitor
all three variables – time, temperature,
and VH202 concentration.

For the first time, you can get a higher level of assurance with VH202 monitoring as you currently expect with steam monitoring. The innovative 3M™ Attest™ Tri-Metric Chemical Indicator, Type 4, streamlines your workflow by displaying clear sterilization results for heat and moisture sensitive instruments.



Detects failed VH202 sterilization cycles more accurately than any other CI tested.¹



Easy-to-read 'accept' and 'reject' zones. Clearly changes from blue to pink.²



The most colour stable of all VH202 CIs tested* – won't change colour before sterilization under ambient light.³



Proprietary protective film prevents ink smear on surgical instruments.⁴



Leading the way for better patient care

The first FDA cleared CI for VH202 sterilization categorised by the FDA as a multivariable CI for chemical vapour sterilization that independently monitors three critical variables.⁵


- time**
- temperature**
- VH202 concentration**

*Competitive Type 1 and Type 4 VH202 CIs.

Expect a higher level of performance.

VH2O2 sterilization is designed to perform at lower temperatures to process heat and moisture sensitive instruments and devices. With the 3M™ Attest™ VH2O2 Tri-Metric Chemical Indicator, Type 4, you get sterilization monitoring results you can count on.

- ▶ Easy-to-find among instruments in a pack
- ▶ Clear results at a glance. Moving front technology eliminates time spent interpreting colour changes.
- ▶ 2D barcode contains lot numbers and expiration dates to scan to your hospital’s current instrument tracking system
- ▶ Blue to pink colour change is consistent on all 3M chemical indicators for vaporized hydrogen peroxide



Product	Product name	Size	Items/bag	Per case
1348	3M™ Attest™ Vaporized Hydrogen Peroxide Tri-Metric Chemical Indicator	5.08 cm	250	4 bags
1348E	3M™ Attest™ Vaporized Hydrogen Peroxide Tri-Metric Chemical Indicator	10.16 cm	250	4 bags
1295	3M™ Attest™ Rapid Readout Biological Indicator (pink cap)	–	30	4 bags
1228	3M™ Comply™ Hydrogen Peroxide Indicator Tape	1.90 cm	–	24 rolls



For more information, or to request a sample, contact your sales representative or visit go.3M.com/TriMetric

3M United Kingdom PLC
 Charnwood Campus
 10 Bakewell Road
 Loughborough
 LE11 5RB
 01509 611611

3M Ireland
 The Iveagh Building
 Carrickmines Park
 Carrickmines
 Dublin 18
 00 353 (01) 280 3555

1 EM-05-692856, EM-05-692858
 2 EM-05-424823
 3 EM-05-692861
 4 EM-05-638321
 5 EM-05-673143

3M, Attest and Comply are trademarks of 3M Company. © 3M 2020. OMG57577.

