



Littmann®
Stethoscopes

Medical

A high degree of hearing.

Quality makes a difference when you're learning auscultation skills to help you diagnose patients. Start with a 3M™ Littmann® Stethoscope to hear the subtle sounds you don't want to miss. It's built to last and helps you perform your best.

The outstanding acoustics of 3M™ Littmann® Cardiology Stethoscopes can help you diagnose patients and make decisions confidently – even with hard-to-hear sounds like aortic stenosis murmurs, S3 gallop rhythms and faint pulmonary anomalies.

Designed to provide excellent sound quality, these stethoscopes feature tunable diaphragm technology, so you can hear low and high frequency sounds without repositioning or turning over the chestpiece. The dual-lumen tubing eliminates the rubbing noise of traditional double tubes, and the adjustable headsets and snap-tight soft-sealing eartips ensure an excellent acoustic seal and comfortable fit.



3M™ Littmann® Cardiology IV™ Stethoscope

Outstanding acoustics for detailed diagnostic auscultation with the versatility of adult and pediatric sides of the chestpiece.



3M™ Littmann® Master Cardiology™ Stethoscope

The best acoustic response in the 3M™ Littmann® Mechanical Stethoscope line.

3M™ Littmann® Cardiology IV™ Stethoscope



Standard stainless steel



6151 (22") Black 6152 Black 6154 Navy 6155 Hunter 6156 Plum



6158 Raspberry 6159 Rose Pink 6184 Burgundy

Black/black



6163 Black 6168 Navy 6200 Black 6201 Black 6203 Black

Mirror



6177 Black 6170 Burgundy

Rainbow/black



6165 Black 6205 Plum

High-polish smoke/smoke



6238 Grey

High-polish rainbow/smoke



6240 Black 6241 Raspberry

High-polish rainbow/black



6242 Navy 6239 Plum

High polish smoke/black



6202 Navy 6204 Black 6232 Black 6234 Caribbean

Champagne/black



6176 Burgundy 6190 Caribbean

Champagne/smoke



6179 Black

Product Features

Acoustic Performance	9 (1-10 scale)
Chestpiece Design/Technology	Machined stainless steel, two-sided
Headset Material	Lightweight aluminum alloy
Length/Total Weight	56 cm (22")/167 g, 69 cm (27")/177 g
Limited Warranty*	7 years
Diaphragm Details/Diameter	Two, single-piece tunable diaphragms 4.3 cm (1.7"), 3.4 cm (1.4")



Available in Canada from

3M Canada Company
Medical Solutions Division
300 Tartan Drive
London, Ontario N5V 4M9
Canada
1-800-364-3577
3M.ca

3M Company
2510 Conway Avenue
St. Paul, MN 55144
USA
1-800-228-3957

3M™ Littmann® Master Cardiology™ Stethoscope



Standard stainless steel



2160 Black 2163 Burgundy 2164 Navy 2167 Plum

Black/black



2161 Black

Smoke/smoke



2176 Black

Brass finish/brass finish



2175 Black

Product Features

Acoustic Performance	10 (1-10 scale)
Chestpiece Design/Technology	Polished stainless steel, single-sided
Headset Material	Lightweight aluminum alloy
Length/Total Weight	69 cm (27")/185 g
Limited Warranty*	7 years
Diaphragm Details/Diameter	Tunable 5.1 cm (2.0")

Colours shown are approximate.

*Please visit Littmann.ca for full warranty terms.

**Explore 3M™ Littmann® Stethoscopes at Littmann.ca
or call 1-800-364-3577.**

3M, Littmann, the L Littmann logo, Master Cardiology, the Littmann Quality Seal and Cardiology IV are trademarks of 3M. Used under license in Canada.
© 2021, 3M. All rights reserved. 2104-20330 E