

As dedicated and hardworking as you are.

Choose a classic

The compact 3M Littmann Classic III Stethoscope is perfect for monitoring and assessing children or adults in non-critical settings. The versatile stethoscope offers high acoustic sensitivity for exceptional performance, hours of comfortable use and comes in a variety of vibrant colors to reflect your personal style.

- Use in non-critical settings
- Detect normal and abnormal sounds and rhythms
- Fun and vibrant; match color and finish to your personality



Here's what makes the Classic III stethoscope stand out:



Tunable diaphragm

Tunable diaphragms on both the adult and pediatric sides of the chestpiece. The pediatric side is useful for small or thin patients, around bandages, and for carotid assessment.



Pediatric side

Pediatric side converts to an open bell. Just remove the single-piece diaphragm and replace it with a non-chill rim.



Single-piece diaphragm

Single-piece diaphragm is easy to attach and easy to clean because the surface is smooth without crevices.



Open bell

Open bell stays clear of dirt and debris by covering it with the small diaphragm.



Next-generation tubing

Not made with natural rubber latex or phthalate plasticizers in the tubing or any other component.



5 year promise


Turn to your trusted instrument for years to come with a five-year promise.



3M Littmann Stethoscopes

The personal diagnostic tool you depend on day after day.

Littmann Stethoscopes Acoustic Performance Levels



Lightweight II S.E. Stethoscope

Level 6

Ideal for taking blood pressure readings and making limited patient assessments.



Classic II Pediatric & Infant Stethoscope

Classic III Monitoring Stethoscope

Levels 7-8

Hard-working, dependable, clinical tools for physical assessment and monitoring of your patients, one long shift after another.



Cardiology IV Diagnostic Stethoscope

Master Cardiology Stethoscope

Levels 9-10

Provides outstanding acoustic sensitivity to help you perform detailed diagnostic auscultation. Useful in challenging and critical care environments.



CORE Digital Stethoscope

Up to 40x amplification (at peak frequency, vs. analog mode)

Level 10+

Offers high-performance acoustics with connectivity, amplification and noise reduction. Useful in challenging and critical care environments.