


Take listening to a whole new level.



A New Class of Clarity

The 3M™ Littmann® CORE Digital Stethoscope features the outstanding acoustics, comfort and quality of a Littmann cardiology-grade stethoscope enhanced by powerful digital technology. Features like up to 40x amplification* and active noise cancellation can help you auscultate with greater confidence. The Littmann CORE Digital Stethoscope connects with Eko software on a smart device** to visualize, record and share data.

*40X amplification based on comparison of analog versus electronically assisted amplified listening with maximum volume at the peak frequency (125 Hz).
**Smart device not included.



Connects to Eko software

Visualize, record and share heart sound waveforms by downloading the free Eko app on a smart device**



Up to 40x amplification

Up to 40x amplification (at peak frequency, vs. analog mode)



Active noise cancellation

Noise cancellation reduces unwanted background sounds



Analog & digital modes

Toggle between analog and amplified listening modes



Designed for comfort

Soft-sealing ear tips provide an excellent acoustic seal and comfortable fit



Tunable diaphragm

Tunable diaphragms on both the adult and pediatric sides of the chestpiece. The pediatric side is useful for small or thin patients, around bandages, and for carotid assessment.



Capture. Analyze. Share.

Capture, save and annotate 15, 30, 60 or 120-second recordings in the secure dashboard. The Littmann CORE Digital Stethoscope connects to Eko software on a smartphone, to create a library of heart sounds, so clinicians can monitor the progression of a disease or for education, and securely share recordings with colleagues for collaboration.