

3M™ Snap™ Therapy System – specifications sheet

3M™ Snap™ Therapy Cartridges

- Available in three preset pressure settings:
 - 75mmHg (Orange/Grey)
 - 100mmHg (Blue/Grey)
 - 125mmHg (Green/Grey)
- **Manages exudate & drainage in a 60mL reservoir**
- **Cartridge dimensions:**
 - 2.3"W x 5.2"H x 1"D (5.84cm x 12.97cm x 2.54cm)
- **Cartridge weight: 0.15lbs (68g)**
- **Material: Polycarbonate (PC)**



3M™ Snap™ Plus -125mmHg Therapy Cartridge

- Available in one preset pressure setting:
 - 125mmHg (Green/Grey)
- **Manages exudate & drainage in a 150mL reservoir**
- **Cartridge dimensions:**
 - 2.96"W x 7.3"H x 1.4"D (7.53cm x 18.54cm x 3.62cm)
- **Cartridge weight: 0.37lbs (167.4g)**
- **Material: Polycarbonate (PC)**




3M™ Snap™ Advanced Dressing Kit

- **All kits include:**
 - 1x proprietary hydrocolloid dressing
 - 1x hydrophobic blue foam
 - 1x one-way tube fitting (grey)
- **Hydrocolloid dressings available in two sizes:**
 - 10cm x 10cm (with 8cm x 8cm foam wound filler)
 - 15cm x 15cm (with 13cm x 13cm foam wound filler)

Dressing material:


- Proprietary hydrocolloid & polymer tubing

- **Storing information:**
 - Store dressings at room temperature
- **Foam material & porosity:**
 - Polyurethane foam with 400-600 micron pores



3M™ Snap™ Bridge dressing kit

- Purpose:
 - In addition to the standard dressing, this bridge dressing allows for improved application ease-of-use for planter surface wounds
- All kits include:
 - 1x foam bridge dressing (14cm x 11cm)
 - 1x hydrophobic blue foam (8cm x 8cm)
 - 1x one-way tube fitting (grey)
- Dressing material: Proprietary hydrocolloid & polymer tubing
- Storing information: Store dressings at room temperature
- Foam material & porosity: Polyurethane foam with 400-600 micron pores



3M™ Snap™ Therapy Strap (sizes & types)

- **Purpose:** Strap to attach the 3M™ Snap™ Therapy Cartridge to the body to help improve patient therapy discretion
- Available in the following sizes:
 - Small, length: 46cm
 - Medium, length: 53cm
 - Large, length: 61cm
- **Material:** Cotton
- **Color:** Black

for 3M™ Snap™ Therapy System



for 3M™ Snap™ Plus Therapy System



3M™ Snap™ SecurRing™ Hydrocolloid Skin Barrier (10 Pack)

- **Purpose:** Helps provide a seal in applications on uneven or challenging anatomy
- **Components:** 10 individually sealed pouches (sterile)
- **Type of material:** Proprietary Hydrocolloid
- **Size:** 5cm in diameter



	Product	Code	Ordering option
Therapy Cartridges	3M™ Snap™ -125mmHg Therapy Cartridge	SNPA125US	Eaches *
	3M™ Snap™ -100mmHg Therapy Cartridge	SNPA100US	Eaches *
	3M™ Snap™ -75mmHg Therapy Cartridge	SNPA075US	Eaches *
	3M™ Snap™ Plus -125mmHg Therapy Cartridge	SNPA125PLUS	Eaches *
Dressings	3M™ Snap™ 10 × 10cm Advanced Dressing Kit	SKTF10X10	Eaches *
	3M™ Snap™ 15 × 15cm Advanced Dressing Kit	SKTF15X15	Eaches *
	3M™ Snap™ 14 × 11cm Bridge Dressing Kit	BKTF14X11	Eaches *
	3M™ Snap™ 14 × 11cm Bridge Dressing Kit with SecurRing™ Hydrocolloid Skin Barrier	BKTF14X11S	Eaches *
Accessories	3M™ Snap™ SecurRing™ Hydrocolloid Skin Barrier	SRNG10	Eaches
	3M™ Snap™ Therapy Strap - Small	STPAS	Eaches
	3M™ Snap™ Therapy Strap - Medium	STPAM	Eaches
	3M™ Snap™ Therapy Strap - Large	STPAL	Eaches
	3M™ Snap™ Plus Therapy Strap - Small	STPASP	Eaches
	3M™ Snap™ Plus Therapy Strap - Medium	STPAMP	Eaches
	3M™ Snap™ Plus Therapy Strap - Large	STPALP	Eaches

* also available in cases of 10



3M Company
2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-275-4524 (NPWT products)
1-800-228-3957
Web 3M.com/medical

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals. Rx only.

© 2022 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.
70-2011-8275-8 PRA-PM-US-02122 (02/22)