



Prevena Restor™

Postoperative Recovery

Specification sheet

Item	Product description	Specifications	
3M™ Prevena™ Plus 125 Therapy Unit	<ul style="list-style-type: none"> Negative pressure options: Pre-set, continuous negative pressure wound therapy at -125 mmHg for up to 14 days Disposable, single patient use Rechargeable battery 	Weight with empty canister: 0.64 lbs (0.29 kg) Dimensions: Approx 8.9 cm x 16.3 cm x 5.49 cm	
3M™ Prevena Restor™ Arthro•Form™ Dressing Kits*		Polyurethane foam bolster	45.7 cm (A) x 31.3 cm (B) x 29.8 cm (C)
		Surface area of foam	768.9 cm ²
		Polyurethane drape	53.2 cm (A) x 38.9 cm (B) x 34.8 cm (C)
	3M™ SensaT.R.A.C.™ Pad/Technology	Tube length: 112 cm	
		Polyurethane foam bolster	33 cm (A) x 24.9 cm (B) x 29.8 cm (C)
		Surface area of foam	588.1 cm ²
Polyurethane drape		43.1 cm (A) x 35.1 cm (B) x 34.8 cm (C)	
SensaT.R.A.C. Pad/Technology	Tube length: 112 cm		
3M™ Prevena Restor™ Axio•Form™ Dressing Kits*		Polyurethane foam bolster	28 cm (A) x 7.5 cm (B) x 29 cm (C)
		Surface area of foam	635 cm ²
	Polyurethane drape	34.5 cm (A) x 15 cm (B) x 35.5 cm (C)	
	SensaT.R.A.C. Pad/Technology	Tube length: 112 cm	
3M™ Prevena Restor™ Bella•Form™ Dressing Kits*		Polyurethane foam bolster	29.0 cm (A) x 27.4 cm (B)
		Surface area of foam	576.8 cm ²
		Polyurethane drape	35.7 cm (A) x 33.7 cm (B)
	SensaT.R.A.C. Pad/Technology	Tube length: 112 cm	
		Polyurethane foam bolster	23.8 cm (A) x 22.8 cm (B)
		Surface area of foam	368.9 cm ²
		Polyurethane drape	30.5 cm (A) x 29.1 cm (B)
	SensaT.R.A.C. Pad/Technology	Tube length: 112 cm	
		Polyurethane foam bolster	21.6 cm (A) x 19.7 cm (B)
Surface area of foam		268.7 cm ²	
Polyurethane drape		28.2 cm (A) x 25.6 cm (B)	
SensaT.R.A.C. Pad/Technology	Tube length: 112 cm		
3M™ Prevena Restor™ Roto•Form™ Dressing Kits*		Polyurethane foam bolster	29 cm (A) x 31 cm (B) x 15 cm (C)
		Surface area of foam	548 cm ²
		Polyurethane drape	35 cm (A) x 36 cm (B) x 18 cm (C)
		SensaT.R.A.C. Pad/Technology	Tube length: 112 cm
3M™ Prevena™ Plus 150ml Canister	<ul style="list-style-type: none"> 150mL canister with connecting tubing Single use, sterile, disposable 	Canister	15.5 cm x 9 cm x 2.5 cm
		3M™ Prevena™ Plus Connector	Tube length: 45 cm
3M™ Prevena™ Y-Connector	Connects two Prevena Restor™ Bella•Form™ Dressings to a Prevena Plus 125 Therapy Unit	Dimensions:	3.4 cm x 1.3 cm x 2.6 cm
		Weight:	3.5 grams

*Dressings contain a skin friendly wicking interface layer that includes 0.019% ionic silver and polyurethane foam bolster, specially designed for application to the surgical area. Therapy should be continuous for a maximum of 7 or 14 days per dressing.

3M™ Prevena Restor™ Therapy System environmental conditions

Operating temperature range: 41°F (5°C) to 104°F (40°C)

Altitude range for optimum performance: -1,253 ft to 9,878 ft (-381.9 m to 3010 m)

Relative humidity range: 15-93%, non-condensing

Electrical specifications (3M™ Prevena™ Plus 125 Therapy Unit)

A rechargeable Lithium-ion battery pack with AC Mains disconnection IEC Classification: Medical Electrical Equipment

Electrical Safety: IEC 60601-1 (3.1 ed)

EMC: IEC 60601-1-2 (4th Ed)

IP22: Protection against solid objects greater than 12.5 mm and against liquid water falling for short periods of time

- Suitable for home use

Contains latex	Contains DEHP (Di(2-ethylhexyl)phthalate)	Contains mercury	Contains haz mat	Contains thimerosal	Radiopaque
No	No	No	No	No	No

Ordering information

SKU	Description	UOM
PRE5001	3M™ Prevena Restor™ Arthro●Form™ Incision Management System – 33 cm X 30 cm	1
PRE5055	3M™ Prevena Restor™ Arthro●Form™ Dressing Kits – 33 cm X 30 cm	5
PRE5101	3M™ Prevena Restor™ Arthro●Form™ Incision Management System – 46 cm X 30 cm	1
PRE5155	3M™ Prevena Restor™ Arthro●Form™ Dressing Kits – 46 cm X 30 cm	5
PRE5501	3M™ Prevena Restor™ Axio●Form™ Incision Management System – 29 cm X 28 cm	1
PRE5555	3M™ Prevena Restor™ Axio●Form™ Dressing Kits – 29 cm X 28 cm	5
PRE5221	3M™ Prevena Restor™ Bella●Form™ Incision Management System – 21 cm X 19 cm	1
PRE5255	3M™ Prevena Restor™ Bella●Form™ Dressing Kits – 21 cm X 19 cm	5
PRE5321	3M™ Prevena Restor™ Bella●Form™ Incision Management System – 24 cm X 22 cm	1
PRE5355	3M™ Prevena Restor™ Bella●Form™ Dressing Kits – 24 cm X 22 cm	5
PRE5421	3M™ Prevena Restor™ Bella●Form™ Incision Management System – 29 cm X 27 cm	1
PRE5455	3M™ Prevena Restor™ Bella●Form™ Dressing Kits – 29 cm X 27 cm	5
PRE5601	3M™ Prevena Restor™ Roto●Form™ Incision Management System 31 cm X 29 cm	1
PRE5655	3M™ Prevena Restor™ Roto●Form™ Dressing 31 cm X 29 cm	5

For more information, visit [3M.com/PrevenaRestor](https://www.3M.com/PrevenaRestor) or call **800-275-4524**



3M Company

2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-275-4524 (NPWT products)
1-800-228-3957

Web [3M.com/Medical](https://www.3M.com/Medical)

NOTE: Specific indications, limitations, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.

©2021 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks.
Unauthorized use prohibited. PRA-PM-US-01074 (07/21)