



Prevena Restor™

Postoperative Recovery

Specification sheet

Item	Product description	Specifications		
3M™ Prevena Restor™ ArthroForm™ Dressing Kits		Polyurethane foam bolster	45.7 cm (A) x 31.3 cm (B) x 29.8 cm (C)	
		Surface area of foam	768.9 cm ²	
		Polyurethane drape	53.2 cm (A) x 38.9 cm (B) x 34.8 cm (C)	
	3M™ SensaT.R.A.C.™ Pad/Technology	Tube length: 112 cm		
			Polyurethane foam bolster	33 cm (A) x 24.9 cm (B) x 29.8 cm (C)
			Surface area of foam	588.1 cm ²
Polyurethane drape	43.1 cm (A) x 35.1 cm (B) x 34.8 cm (C)			
3M™ SensaT.R.A.C.™ Pad/Technology	Tube length: 112 cm			
	3M™ Prevena Restor™ AxioForm™ Dressing Kits		Polyurethane foam bolster	28 cm (A) x 7.5 cm (B) x 29 cm (C)
			Surface area of foam	635 cm ²
Polyurethane drape		34.5 cm (A) x 15 cm (B) x 35.5 cm (C)		
3M™ SensaT.R.A.C.™ Pad/Technology		Tube length: 112 cm		
3M™ Prevena Restor™ BellaForm™ Dressing Kits		Polyurethane foam bolster	29.0 cm (A) x 27.4 cm (B)	
		Surface area of foam	576.8 cm ²	
		Polyurethane drape	35.7 cm (A) x 33.7 cm (B)	
	3M™ SensaT.R.A.C.™ Pad/Technology	Tube length: 112 cm		
			Polyurethane foam bolster	23.8 cm (A) x 22.8 cm (B)
			Surface area of foam	368.9 cm ²
	Polyurethane drape		30.5 cm (A) x 29.1 cm (B)	
	3M™ SensaT.R.A.C.™ Pad/Technology	Tube length: 112 cm		
			Polyurethane foam bolster	21.6 cm (A) x 19.7 cm (B)
Surface area of foam			268.7 cm ²	
Polyurethane drape	28.2 cm (A) x 25.6 cm (B)			
3M™ SensaT.R.A.C.™ Pad/Technology	Tube length: 112 cm			
	3M™ Prevena Restor™ RotoForm™ Dressing Kits		Polyurethane foam bolster	29 cm (A) x 31 cm (B) x 15 cm (C)
			Surface area of foam	548 cm ²
Polyurethane drape		35 cm (A) x 36 cm (B) x 18 cm (C)		
3M™ SensaT.R.A.C.™ Pad/Technology		Tube length: 112 cm		
3M™ Prevena Restor™ AdaptiForm™ Dressing Kits		Polyurethane foam bolster	49 cm (A) x 28 cm (B)	
		Surface area of foam	1,372 cm ²	
		3M™ V.A.C.® Drape (6)	38.7 cm x 26 cm (Each)	
		Hydrocolloid strips (10)	4.5 cm x 24 cm (with liner) 2 cm x 24 cm (without liner)	
		3M™ SensaT.R.A.C.™ Pad/Technology	Tube length: 112 cm	

- NOTE:**
- Dressings contain 1) polyurethane foam bolster specially designed for application to the surgical area, and 2) a skin-friendly wicking interface layer that includes 0.019% ionic silver
 - For maximum benefit, dressings should be applied immediately post surgery to clean, surgically closed incisions and continuously applied for up to a maximum of 7 days

Specification sheet (continued)

Item	Product description	Specifications
3M™ Prevena™ Plus 125 Therapy Unit (7 and 14 day therapy units available)	<ul style="list-style-type: none"> Negative pressure options: Pre-set, continuous negative pressure wound therapy at -125 mmHg Disposable, single patient use Rechargeable battery 	Weight with empty canister: 0.64 lbs (0.29 kg) Dimensions: Approx 8.9 cm x 16.3 cm x 5.49 cm
3M™ Prevena™ Plus 150ml Canister	<ul style="list-style-type: none"> 150mL canister with connecting tubing Single use, sterile, disposable 	Canister
		3M™ Prevena™ Plus Connector
3M™ Prevena™ Y-Connector	Connects two Prevena Restor™ Bella●Form™ Dressings to a Prevena Plus 125 Therapy Unit	Dimensions:
		Weight:

Environmental Conditions (3M™ Prevena Restor™ Dressings)

Operating temperature range: 41°F (5°C) to 104°F (40°C)

Altitude range for optimum performance: -1,253 ft to 9,878 ft (-381.9 m to 3010 m)

Relative humidity range: 15-93%, non-condensing

Electrical specifications (3M™ Prevena™ Plus 125 Therapy Unit)

A rechargeable Lithium-ion battery pack with AC Mains disconnection IEC Classification: Medical Electrical Equipment

Electrical Safety: IEC 60601-1 (3.1 ed)

EMC: IEC 60601-1-2 (4th Ed)

IP22: Protection against solid objects greater than 12.5 mm and against liquid water falling for short periods of time

- Suitable for home use

Ordering information

SKU	Description	UOM
PRE5001	3M™ Prevena Restor™ Arthro●Form™ Incision Management System – 33 cm X 30 cm	1
PRE5055	3M™ Prevena Restor™ Arthro●Form™ Dressing Kits – 33 cm X 30 cm	5
PRE5101	3M™ Prevena Restor™ Arthro●Form™ Incision Management System – 46 cm X 30 cm	1
PRE5155	3M™ Prevena Restor™ Arthro●Form™ Dressing Kits – 46 cm X 30 cm	5
PRE5501	3M™ Prevena Restor™ Axio●Form™ Incision Management System – 29 cm X 28 cm	1
PRE5555	3M™ Prevena Restor™ Axio●Form™ Dressing Kits – 29 cm X 28 cm	5
PRE5221	3M™ Prevena Restor™ Bella●Form™ Incision Management System – 21 cm X 19 cm	1
PRE5255	3M™ Prevena Restor™ Bella●Form™ Dressing Kits – 21 cm X 19 cm	5
PRE5321	3M™ Prevena Restor™ Bella●Form™ Incision Management System – 24 cm X 22 cm	1
PRE5355	3M™ Prevena Restor™ Bella●Form™ Dressing Kits – 24 cm X 22 cm	5
PRE5421	3M™ Prevena Restor™ Bella●Form™ Incision Management System – 29 cm X 27 cm	1
PRE5455	3M™ Prevena Restor™ Bella●Form™ Dressing Kits – 29 cm X 27 cm	5
PRE5601	3M™ Prevena Restor™ Roto●Form™ Incision Management System – 31 cm X 29 cm	1
PRE5655	3M™ Prevena Restor™ Roto●Form™ Dressing Kits – 31 cm X 29 cm	5
PRE6001	3M™ Prevena Restor™ Adapti●Form™ Incision Management System – 49 cm X 28 cm	1
PRE6055	3M™ Prevena Restor™ Adapti●Form™ Dressing Kits – 49 cm X 28 cm	5

For more information, visit [3M.com/PrevenaRestor](https://www.3M.com/PrevenaRestor) or call **800-275-4524**



3M Company

2510 Conway Ave.
St. Paul, MN 55144 USA

Phone 1-800-275-4524 (NPWT products)
1-800-228-3957

Web 3M.com/Medical

NOTE: Specific indications, limitations, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. Rx only.

©2022 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks.
Unauthorized use prohibited. 70-2013-1517-6