








# Prevena Restor™

Postoperative Recovery

	Item	Product Description	Specifications		
<b>Therapy Unit</b>	3M™ Prevena™ Plus 125 Therapy Unit (14 Day)	<ul style="list-style-type: none"> <li>Negative Pressure Setting: continuous NPWT at -125 mmHg for up to 14 days</li> <li>Disposable, single patient use</li> <li>Rechargeable battery</li> </ul>	<b>Weight with empty canister:</b> 0.29 kg (0.64 lbs) <b>Dimensions:</b> Approx. 8.9 x 16.3 x 5.49 cm		
<b>Dressings</b>	Specifications	<b>Arthro-Form™ Dressings</b>			
			Polyurethane foam bolster <span style="color: red;">■</span> 45.7 cm Polyurethane foam bolster <span style="color: green;">■</span> 31.3 cm Polyurethane foam bolster <span style="color: yellow;">■</span> 29.8 cm Surface area of foam 768.9 cm <sup>2</sup> Polyurethane drape <span style="color: red;">■</span> 53.2 cm Polyurethane drape <span style="color: green;">■</span> 38.9 cm Polyurethane drape <span style="color: yellow;">■</span> 34.8 cm SensaT.R.A.C.™ Dressing Tube length: 112 cm		
			Polyurethane foam bolster <span style="color: red;">■</span> 33 cm Polyurethane foam bolster <span style="color: green;">■</span> 24.9 cm Polyurethane foam bolster <span style="color: yellow;">■</span> 29.8 cm Surface area of foam 588.1 cm <sup>2</sup> Polyurethane drape <span style="color: red;">■</span> 43.1 cm Polyurethane drape <span style="color: green;">■</span> 35.1 cm Polyurethane drape <span style="color: yellow;">■</span> 34.8 cm SensaT.R.A.C.™ Dressing Tube length: 112 cm		
		<b>Bella-Form™ Dressings</b>			
			Polyurethane foam bolster <span style="color: red;">■</span> 29.0 cm Polyurethane foam bolster <span style="color: green;">■</span> 27.4 cm Surface area of foam 576.8 cm <sup>2</sup> Polyurethane drape <span style="color: red;">■</span> 35.7 cm Polyurethane drape <span style="color: green;">■</span> 33.7 cm SensaT.R.A.C.™ Dressing Tube length: 112 cm		
			Polyurethane foam bolster <span style="color: red;">■</span> 23.8 cm Polyurethane foam bolster <span style="color: green;">■</span> 22.8 cm Surface area of foam 368.9 cm <sup>2</sup> Polyurethane drape <span style="color: red;">■</span> 30.5 cm Polyurethane drape <span style="color: green;">■</span> 29.1 cm SensaT.R.A.C.™ Pad Tube length: 112 cm		
			Polyurethane foam bolster <span style="color: red;">■</span> 21,6 cm Polyurethane foam bolster <span style="color: green;">■</span> 19,7 cm Surface area of foam 268,7 cm <sup>2</sup> Polyurethane drape <span style="color: red;">■</span> 28,2 cm Polyurethane drape <span style="color: green;">■</span> 25,6 cm SensaT.R.A.C.™ Dressing Tube length: 112 cm		
		<b>Canister</b>	3M™ Prevena™ Plus 150 ml Canister	<ul style="list-style-type: none"> <li>150 ml canister with connecting tubing</li> <li>Single use, sterile, disposable</li> </ul>	Canister 15.5 x 9 x 2.5 cm 3M™ Prevena™ Plus Connector Tube length: 45 cm
		<b>Connector</b>	V.A.C.® Y-Connector	Connects two 3M™ Prevena Restor™ Bella-Form™ Dressings to a 3M™ Prevena™ Plus Therapy Unit or approved V.A.C.® Therapy Unit	Dimensions: 3.4 x 1.3 x 2.6 cm
					Weight: 3.5 grams

## 3M™ Prevena Restor™ Therapy System Environmental Conditions

- Operating Temperature Range: 5°C (41°F) to 40°C (104°F)
- Altitude Range for Optimum Performance: -381.9 m to 3010 m (-1,253 to 9,878 feet)
- Relative Humidity Range: 15-93%, non-condensing

## Electrical Specifications (3M™ Prevena™ Plus 125 Therapy Unit)

A rechargeable Lithium-ion battery pack with AC Mains disconnection IEC Classification: Medical Electrical Equipment

- Electrical Safety: IEC 60601-1 (3.1 ed)
- EMC: IEC 60601-1-2 (4th Ed)
- IP22: Protection against solid objects greater than 12.5 mm and against liquid water falling for short periods of time
- Suitable for home use

Contains Latex	Contains DEHP (Di(2-ethylhexyl) phthalate)	Contains Mercury	Contains Haz Mat	Contains Thimerosal	Radiopaque
No	Yes*	No	No	No	No

\*Contained within SensaT.R.A.C.™ Tubing

## Ordering Information

SKU	Product Name	Quantity
PRE5001	3M™ Prevena Restor™ Arthro•Form™ System Kit – 33 cm x 30 cm	1
PRE5055	3M™ Prevena Restor™ Arthro•Form™ Dressing – 33 cm x 30 cm	5
PRE5101	3M™ Prevena Restor™ Arthro•Form™ System Kit – 46 cm x 30 cm	1
PRE5155	3M™ Prevena Restor™ Arthro•Form™ Dressing – 46 cm x 30 cm	5
PRE5221	3M™ Prevena Restor™ Bella•Form™ System Kit – 21 cm x 19 cm	1
PRE5255	3M™ Prevena Restor™ Bella•Form™ Dressing – 21 cm x 19 cm	5
PRE5321	3M™ Prevena Restor™ Bella•Form™ System Kit – 24 cm x 22 cm	1
PRE5355	3M™ Prevena Restor™ Bella•Form™ Dressing – 24 cm x 22 cm	5
PRE5421	3M™ Prevena Restor™ Bella•Form™ System Kit – 29 cm x 27 cm	1
PRE5455	3M™ Prevena Restor™ Bella•Form™ Dressing – 29 cm x 27 cm	5

For more information about the 3M™ Prevena Restor™ Therapy System, contact your local representative.



Available in Canada from your authorized 3M-KCI distributors.  
KCI USA, Inc., a 3M Company  
KCI owned and operated by 3M Company

**KCI Medical Canada Inc.**  
75 Courtneypark Dr W, Unit 4  
Mississauga, ON  
L5W 0E3

**KCI USA, INC.**  
12930 IH 10 West  
San Antonio, TX  
78249

Follow local institutional protocols for infection control and waste disposal procedures. Local protocols should be based on the applicable local government environmental regulations.

**Note:** Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals.

© 2021, 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited. Used under license in Canada. PRA-PM-CA-00186 2105-20642 E