

Introducing the next generation of a clinical icon.



Open bell versatility.

The paediatric side converts to a traditional open bell by replacing the single-piece diaphragm with a non-chill rim.

Introducing the

3M™ Littmann® Classic III™ Stethoscope

Right out of the box, it lives up to its name.

The Littmann Classic III stethoscope is the latest version of the stethoscope that helps millions of medical professionals achieve their best. It brings new design, materials, and technology to the Littmann stethoscope series trusted by clinicians around the world. The new Classic III stethoscope offers you:



High acoustic sensitivity for performing general physical assessments.



Tunable diaphragms on both the adult and paediatric sides of the chestpiece.



Single-piece diaphragm that's easy to attach and easy to clean because the surface is smooth without crevices.



Next-generation tubing with improved resistance to skin oils and alcohol for longer tubing life. Less likely to pick up stains.

Using a Classic III stethoscope demonstrates what's important to you as a medical professional: quality, achievement, and providing the best possible care. Every patient. Every day.

Here's what makes the

3M™ Littmann® Classic III™ Stethoscope stand out.

Features

Tunable technology is a 3M invention that lets you shift between hearing high and low frequency sounds by simply adjusting how much you press on the chestpiece.

Soft-sealing eartips give you a comfortable fit and an excellent acoustic seal; they snap tightly onto the eartubes for safety.

Open bell stays clear of dirt and debris by keeping it covered with the small diaphragm.

Adjustable headset lets you easily set the tension for head size and comfort just by squeezing or pulling apart the eartubes.

No natural rubber latex or phthalate plasticisers are used in the tubing or any other component—to help protect human health and the environment.

Classic III stethoscope specifications	
Length	69cm / 27.2 inches
Weight	150g
Acoustic performance	7
CHESTPIECE	
Construction	Two-sided, adult and paediatric
Material	Stainless steel
Diaphragm type	Tunable, single-piece
Diaphragm diameter	Adult: 4.3cm / 1.7 in Paediatric: 3.3cm / 1.3 in
TUBING	
Construction	Single lumen (single acoustic channel within tube)
Material	Not made with natural rubber latex or phthalate plasticisers
HEADSET	
Material	Wide diameter aerospace alloy / anodized aluminium

Available tubing colours and finish combinations

TUBING COLOUR	STANDARD FINISH	SPECIAL FINISH	TUBING COLOUR	STANDARD FINISH	SPECIAL FINISH
Black	5620	Black 5803, Smoke 5811	Grey	5621	
Navy Blue	5622		Caribbean Blue	5623	Rainbow 5807
Raspberry	5626	Rainbow 5806	Hunter Green	5624	
Burgundy	5627		Orange	5629	
Chocolate		Copper 5809	Ceil Blue	5630	
Lavender	5832		Peal Pink	5633	
Turquoise	5835		Lime Green	5829	
Lemon-Lime	5839		Plum	5831	
Emerald	5840		Dark Olive Green		Smoke 5812
Rose Pink	5639				

Colours and Finishes

We periodically introduce new tubing colours and chestpiece finish combinations for the Classic III stethoscope.

Visit either www.littmann.com.au or www.littmann.co.nz for the latest product information.

Due to limitations of the printing process, the colors shown are approximations of the actual tubing colors.

For healthcare professional use only.
Registered under Act 737.
MDA Reg. No. : GMD29659566118A