

# A complete solution.

3M Sample Collection solutions aligned with  
3M™ Petrifilm™ Plates.

## Sample Collection Products



### 3M™ Swab-Samplers

- Self-contained swab
- Ideal for sampling hard-to-reach areas with long, pre-hydrated 9 cm swab
- Offered in multiple volumes for qualitative or quantitative use
- Offered in multiple buffers



### 3M™ Quick Swabs

- Ready and easy-to-use
- Quickly delivers 1 ml sample onto a 3M™ Petrifilm™ Plate
- Labour savings in dilution mixing and pipetting



### 3M™ Environmental Scrub Sampler with 10 ml Wide Spectrum Neutralizer

- Multi-layered non-woven composite material maintains integrity when sampling rough areas and scrub dot technology disrupts biofilm and enhances sample collection
- Stick is metal detectable and features an easy-to-use 1-step twist detachment mechanism
- 3M Wide Spectrum Neutralizer is formulated to effectively neutralize common sanitizers and is designed to be compatible with 3M™ Petrifilm™ Plates
- AOAC® Performance Tested Method™ Certificate #022104



### 3M™ Hydrated-Sponges and 3M™ Sponge-Sticks

- Ready-to-use, pre-hydrated sponges for easy and convenient sample collection
- Biocide-free cellulose sponge maintains viability for a wide range of organisms
- Durable sponge withstands scrubbing
- Offered in multiple buffers
- Allows sampling without direct handling of the sponge



### 3M™ Sample Bags

- Durable, quality construction
- Easy-to-open bag
- Wire ties for strong seals and reduced leakage

## Sample Preparation and Processing Products




### 3M™ Flip-Top Dilution Bottles

- Leak-resistance, strong, lightweight clear PET plastic
- Pre-filled, eliminate making and sterilizing media
- Hinged lid stays open for easy sample addition
- Offered in multiple buffers



### 3M™ Mini Flip-Top Dilution Vials

- Hinged lid stays open for easy sample addition
- Requires less sample material than regular Flip-Top dilution bottles
- Offered in multiple buffers



### 3M™ Homogenizer Bags

- Quality construction to survive your paddle blender



### 3M™ Dilution Pouches

- Flat bottom stand-up pouch ensures stability
- Wire tie holds top open for easy sample addition
- Multi-laminated pouch is puncture resistant



### 3M™ Dehydrated Media

- Convenient
- Consistent
- Cost-effective

# Create a seamless sample collection program with 3M's integrated product offerings.

Product	Description	Catalogue #
3M™ Quick Swab	1 ml Lethéen Broth	6432
	1 ml Lethéen Broth	6433
3M™ Swab-Sampler	1 ml Lethéen Broth	RS9601LET
	4 ml Lethéen Broth	RS9604LET
	5 ml Lethéen Broth	RS9605LET
	10 ml Lethéen Broth	RS96010LET
	10 ml Buffered Peptone Water	RS96010BPW
3M™ Environmental Scrub Sampler Stick	10 ml 3M™ Wide Spectrum Neutralizer	ESS10WSN
3M™ Environmental Scrub Sampler and Gloves	10 ml 3M™ Wide Spectrum Neutralizer	HES10WSN2G
3M™ Hydrated-Sponge	10 ml Buffered Peptone Water	HS10BPW
	10 ml Buffered Peptone Water, with 2 gloves	HS10BPW2G
	10 ml Lethéen Broth	HS10LET
	10mL Lethéen Broth, with 2 gloves	HS10LET2G
3M™ Sponge-Stick	Without Media	SSL100
	10 ml Buffered Peptone Water	SSL10BPW
	10 ml Lethéen Broth	SSL10LET
	10 ml Lethéen Broth, with 2 gloves	SSL10LETG
3M™ Dry-Sponge	Dry sponge, with 18 oz. sample bag	BP133ES
	Dry sponge with 18 oz. sample bag, with 2 gloves	BP237SPG
3M™ Lethéen Broth	500 g	BP0245500
3M™ Mini Flip-Top Vials with Lethéen	9 ml Lethéen Broth	BPPFV9LB
	10 ml Lethéen Broth	BPPFV10LB
3M™ Sample Bags, Plain	3" x 7", 4 oz., 2.5 mil	BP111S
	4.5" x 7.5", 12 oz., 2.5 mil	BP134S
	4.5" x 9", 15 oz., 2.5 mil	BP112S
	5.5" x 9", 22 oz., 3 mil	BP113S
	4.5" x 12", 27 oz., 2.5 mil	BP118S
	4.5" x 15", 36 oz., 2.5 mil	BP119S
	5.5" x 15", 50 oz., 3 mil	BP121S
	12" x 12", 120 oz., 4 mil	BP41212
	10" x 15", 135 oz., 4 mil	BP1015LP
12" x 18", 228 oz., 4 mil	BP41218	
3M™ Sample Bags, Printable	3" x 5", 2 oz., 2.5 mil	BP123S
	3" x 7", 4 oz., 2.5 mil	BP114S
	4.5" x 9", 15 oz., 2.5 mil	BP115S
	5.5" x 9", 22 oz., 3 mil	BP127S
	7" x 12", 55 oz., 3 mil	BP701
	10" x 15", 135 oz., 4 mil	BP1015

Product	Description	Catalogue #
3M™ Mini Flip-Top Vials	9 ml Butterfields Buffer	BPPFV9BFD
	9 ml Buffered Peptone Water	BPPFV9BPW
	10 ml Buffered Peptone Water	BPPFV10BPW
	9 ml MRD, Max Rec Diluent	BPPFV9MRD
	9 ml MRS Broth	BPPFV9MRS
3M™ Flip-Top Dilution Bottles	90 ml Buffered Peptone Water (ISO)	FTBPW19060
	90 ml Buffered Peptone Water	FTBPW9060
	90 ml Butterfields Buffer	FTBFD9060
	99 ml Butterfields Buffer	FTBFD9960
	90 ml Peptone Water	FTPW9060
	99 ml Peptone Water	FTPW9960
	90 ml Phosphate Buffer	FTPHB9060
3M™ Sample Bags, Blender	No wire tie, 7" x 12", 55 oz., 3 mil	BP912S
	No wire tie, 7" x 12", 55 oz., 4 mil	BP913S
3M™ Sample Bags, Blender with Filter	No wire tie, 7" x 12", 55 oz., 3 mil	BP914S
	No wire tie, 15" x 20", 405 oz., 3 mil	BP1520F
	No wire tie, 6.75" x 12", 55 oz., 3 mil	6469
3M™ Quick Dilution Pouch	99 ml Butterfields Buffer	QDBFD99
	225 ml Butterfields Buffer	QDBFD225
	225 ml Butterfields Buffer, with filter	QDBFD225F
	450mL Butterfields Buffer	QDBFD450
3M™ Bird Rinse	400 ml Butterfields Buffer	BRBFD40030
	400 ml Buffered Peptone Water	BRBPW40030
3M™ Buffered Peptone Water Broth	500 g	BP0115500
	2.5 kg	BP0115005
3M™ Buffered Peptone Water, ISO	500 g	BPW500
	2.5 kg	BPW025
	10 kg	BPW010
3M™ Salmonella Enrichment Base	500 g	SEB500
	2.5 kg	SEB025
3M™ Salmonella Enrichment Supplement	1 g/vial	SESUP001



**3M Food Safety Division**  
**3M Canada**  
 P.O. Box 5757  
 London, ON N6A 4T1  
 1-800-364-3577

3M, 3M Science. Applied to Life. and Petrifilm are trademarks of 3M.  
 Used under license in Canada. © 2021, 3M. All rights reserved.  
 All other trademarks are property of their respective owners. 2106-21003 E