

3M Science.
Applied to Life.™



3M™ Imprint™ 4
VPS Impression Material

Stress less over everyday impressions.

Say 'Yes' to the VPS impression material with the shortest intra-oral setting time. Move through your day with speed, control and precision.

Speed.


















Reduced chair time and stress for patients – Imprint 4 material has the shortest intra-oral setting time among leading VPS fast-setting impression materials.*

Control.

Long working time for you. Due to the unique chemistry, the setting time doesn't start until the end of the working time. That means you'll have a sufficient amount of time to load and seat the tray – no stressful race against the clock.
















Precision.

Imprint 4 material is hydrophilic, displacing moisture right from the start in the unset state – which helps prevent detail-disrupting voids.

	Product	Dispensing system	Viscosity Low  High	Setting version	Maximum working time (23°C/73°F) 	Maximum intra-oral syringe time (37°C/98°F) 	Intra-oral setting time (37°C/98°F) 
Tray materials	3M™ Imprint™ 4 Penta™ Putty		Putty	Regular set	1:30	—	2:30
	3M™ Imprint™ 4 Penta™ Heavy		Heavy body	Regular set	2:00	—	2:00
	3M™ Imprint™ 4 Penta™ Super Quick Heavy		Heavy body	Fast set	1:15	—	1:15
Wash materials	3M™ Imprint™ 4 Super Quick Ultra-Light			Fast set	—	0:35	1:15
	3M™ Imprint™ 4 Light			Regular set	—	1:00	2:00
	3M™ Imprint™ 4 Super Quick Light			Fast set	—	0:35	1:15
	3M™ Imprint™ 4 Regular			Regular set	—	1:00	2:00
	3M™ Imprint™ 4 Super Quick Regular			Fast set	—	0:35	1:15

Ordering information

Item	Product Information
Intro kits for 3M™ Pentamix™ 3	
71496	3M™ Imprint™ 4 Penta™ Putty Intro Kit
71494	3M™ Imprint™ 4 Penta™ Heavy Intro Kit
71495	3M™ Imprint™ 4 Penta™ Super Quick Heavy Intro Kit
71499	3M™ Imprint™ 4 Penta™ Putty Trial Kit
71497	3M™ Imprint™ 4 Penta™ Heavy Trial Kit
71498	3M™ Imprint™ 4 Penta™ Super Quick Heavy Trial Kit
Refills for 3M™ Pentamix™ Mixing Unit	
71486	3M™ Imprint™ 4 Penta™ Putty Refill
71484	3M™ Imprint™ 4 Penta™ Heavy Refill
71485	3M™ Imprint™ 4 Penta™ Super Quick Heavy Refill
Wash refills for 3M™ Garant™ Dispenser	
71487	3M™ Imprint™ 4 Regular
71488	3M™ Imprint™ 4 Light
71489	3M™ Imprint™ 4 Super Quick Regular
71490	3M™ Imprint™ 4 Super Quick Light
71491	3M™ Imprint™ 4 Super Quick Ultra-Light
Accessories	
77580	3M™ Garant™ Dispenser (1:1/2:1)
77949	3M™ Penta™ Mixing Tips – Red
71505	3M™ Intra-oral Syringes – Green (20pk)
71452	3M™ Garant™ Mixing Tips – Yellow (50pk)
71462	3M™ Intra-oral tips - Yellow (50pk)

Tray materials	Recommended wash materials
Intro kits for 3M™ Pentamix™ 3	
3M™ Imprint™ 4 Penta™ Putty 	3M™ Imprint™ 4 Light  3M™ Imprint™ 4 Regular 
3M™ Imprint™ 4 Penta™ Heavy 	3M™ Imprint™ 4 Light  3M™ Imprint™ 4 Regular 
3M™ Imprint™ 4 Penta™ Super Quick Heavy 	3M™ Imprint™ 4 Super Quick Light  3M™ Imprint™ 4 Super Quick Regular 
Two-step technique	
3M™ Imprint™ 4 Penta™ Putty 	3M™ Imprint™ 4 Super Quick Ultra-Light  3M™ Imprint™ 4 Light 
3M™ Imprint™ 4 Penta™ Heavy 	3M™ Imprint™ 4 Super Quick Light  3M™ Imprint™ 4 Light 



3M Oral Care
3M Centre, Cain Road,
Bracknell
Berkshire RG12 8HT
Tel: 0870 60 800 60

3M Ireland Limited
The Iveagh Building
Carrickmines Park, Carrickmines
Ireland, D18 X015
+353 (0)1 280 3555



 3M.co.uk/Dental
 **08705360036**
 **OralCareUK@mmm.com**

*3M Internal data

3M, Garant, Imprint, Penta and Pentamix are trademarks of 3M Company. Please recycle.
© 3M 2022. All rights reserved.
70-2013-0580-5 (Rev. 3/2022) OMG309055.