



Prevena™
Incision Therapy

One therapy unit compatible with all 3M™ Prevena™ Dressings

3M™ Prevena™ Plus 125 Therapy Unit - 14 day

Description:

Negative pressure options:

- Pre-set, continuous NPWT at -125mmHg for up to 14 days
- Disposable, single patient use
- Rechargeable battery

Specifications:

- Weight with empty canister:
0.29kg
- Dimensions:
Approx 8.9 × 16.3 × 5.49cm



* Actual size of Therapy unit

Help protect surgical incisions

Indications for use:


3M™ Prevena™ Plus 125 Therapy Unit - 14 day is intended to manage the environment of surgical incisions that continue to drain following sutured or stapled closure by maintaining a closed environment and removing exudate via the application of negative pressure wound therapy.

3M™ Prevena™ Therapy System clinical outcomes:

More than 70+ published clinical studies in a variety of surgical incisions and patient populations have demonstrated that the 3M™ Prevena™ Incision Management System is associated with positive clinical outcomes compared to standard post-operative dressings.

3M™ Prevena™ Dressings

- Contain a skin-friendly wicking interface layer, includes 0.019% ionic silver.
- Polyurethane foam bolster specially designed for application to the surgical area.
- Prevena dressings can be applied for a maximum of 7 days and Prevena Restor dressings for a maximum of 14 days.



Ordering information

Product #	Product Name	Qty	Product #	Product Name	Qty
PRE4010.S	3M™ Prevena™ Plus 125 Therapy Unit - 14 day	1	PRE5055.S	3M™ Prevena Restor™ Arthro•Form™ Dressing – 33cm x 30cm	5
PRE4095.S	3M™ Prevena™ Plus 150ml Canister	1	PRE5155.S	3M™ Prevena Restor™ Arthro•Form™ Dressing – 46cm x 30cm	5
PRE1155.S	3M™ Prevena™ Peel and Place Dressing – 13cm	5	PRE5255.S	3M™ Prevena Restor™ Bella•Form™ Dressing – 21cm x 19cm	5
PRE1055.S	3M™ Prevena™ Peel and Place Dressing – 20cm	5	PRE5355.S	3M™ Prevena Restor™ Bella•Form™ Dressing – 24cm x 22cm	5
PRE3255.S	3M™ Prevena™ Plus Peel and Place Dressing – 35cm	5	PRE5455.S	3M™ Prevena Restor™ Bella•Form™ Dressing – 29cm x 27cm	5
PRE4055.S	3M™ Prevena™ Plus Customizable Dressing	5	M6275066/10.S	3M™ V.A.C.® Y-Connector	10

To learn more about the 3M™ Prevena™ Plus Incision Management System, contact your local 3M representative.



KCI Medical Australia Pty Ltd
Building A, 1 Rivett Road
North Ryde NSW 2113
Ph:1300 524 822
3M.com.au/medical

[†] computer and bench top models.

1. Wilkes RP, Kilpadi DV, Zhao Y, Kazala R, McNulty A Closed incision management with negative pressure wound therapy (CIM): biomechanics Surg Innov 2012;19(1):67-75

3M™ Prevena™ Plus 125 Therapy Unit can be continuously applied for up to 14 days, and a 3M™ Prevena™ Restor™ Dressing has a maximum wear time of 14 days.

Note: Specific indications, limitations, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application.

© 2022 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorised use prohibited.
70-2013-1126-6