

3M Science.
Applied to Life.™



Small but mighty.

3M™ Victory Series™ Mini Tubes



Now available
pre-coated with
3M™ APC™ Flash-Free Adhesive!



Pre-coated
adhesive on base,
ready to bond



Time
savings



40% less bonding
time¹



No flash
to clean up



Ideal for
erupting molars

1. Compared with conventional bracket bonding technique, 3M internal data on file.

Seamless digital treatment design.

3M™ Victory Series™ Mini Tubes provide a well-designed, reliable and predictable orthodontic solution for partially erupted second molars. They have a mesially-extended, compound contoured, micro-etched bonding base that allows for superior fit and bond strength. The Mini Tube features a short, low profile hook that enhances patient comfort.

Bonding success and placement accuracy can be especially challenging on posterior teeth. The combination of the 3M™ Digital Bonding Tray and the pre-coated mini second molar tube offers a means to accurate and successful positioning.

Pair with 3M™ Oral Care Portal and 3M™ Digital Bonding Tray as part of your digital bonding solution.




Proven 3M technologies. Superior performance.

Equipped with all of the design features you've come to appreciate, Victory Series Mini Tubes now come pre-coated with 3M™ APC™ Flash-Free Adhesive.

Powered by 3M science, APC Flash-Free Adhesive leverages several of the highly advanced 3M Technology Platforms such as nanotechnology, adhesives and nonwoven materials to create the innovative design of APC Flash-Free Technology.

The APC Flash-Free Adhesive Coated Appliance System has shown a 40% reduction in the time required for bracket placement, positioning, and cleanup versus other bonding systems.¹



3M™ Victory Series™ Mini Tubes Part List

Description	Part Number
Upper Right 2nd Molar .018 (MBT)	4066-4082
Upper Right 2nd Molar .022 (MBT)	4066-5082
Upper Left 2nd Molar .018 (MBT)	4066-4071
Upper Left 2nd Molar .022 (MBT)	4066-5071
Lower Right 2nd Molar .018 (MBT)	4066-4044
Lower Right 2nd Molar .022 (MBT)	4066-5044
Lower Left 2nd Molar .018 (MBT)	4066-4033
Lower Left 2nd Molar .022 (MBT)	4066-5033



3M Oral Care

3M Australia Pty Limited
Building A, 1 Rivett Rd
North Ryde, NSW 2113
Ph: 1300 363 484
3M.com.au/ortho

3M New Zealand Limited
94 Apollo Drive
Rosedale, Auckland 0632
Ph: 0800 441 622
3M.co.nz/ortho

Visit 3M.com.au/APCFF or 3M.co.nz/APCFF or contact your 3M Representative.

3M, "3M Science. Applied to Life.", APC, and Victory Series are trademarks of 3M.
© 3M 2022. All rights reserved.