

In-box Accessories

	Cardiology		Classic			Lightweight
	Cardiology IV™ Stethoscope	Master Cardiology™ Stethoscope	Classic III™ Stethoscope	Classic II Pediatric	Classic II Infant	Lightweight II S.E.
						
ACCESSORIES						
Eartips, large (snap tight, soft-sealing)	✓	✓	✓	✓	✓	✓
Eartips, small (snap tight, soft-sealing)	✓	✓	✓	✓	✓	✓
Nonchill bell sleeve	✓		✓	✓	✓	✓
Special Procedures Adaptor (for difficult site or pediatric auscultation)		✓				

Why 3M™ Littmann® Stethoscopes are recognized as the gold standard.

Symbol of quality

Profession-leading innovations, precision engineering, advanced materials, and more, ensure consistency in every 3M™ Littmann® stethoscope.

Tunable technology

Isolate high- or low-frequency sounds simply by adjusting how much you press on the chestpiece.

Next-generation tubing

Improved resistance to alcohol and skin oils; hangs straight and flexes easily. Does not contain natural latex or phthalate plasticizers to protect the environment and human health.

Comfort for the physician

Adjustable headset, snap-tight soft sealing eartips ensure a comfortable fit and excellent acoustic seal.

A personal touch

Choose from a variety of 3M™ Littmann® Stethoscope models to reflect your style and meet your auscultation needs.

Visit our Product Finder for more information on choosing the right stethoscope.

See our Where to Buy page for a list of authorized 3M™ Littmann® distributors in Malaysia.

For healthcare professional use only.

Registered under Act 737.

MDA Reg. No. : GMD69864504618A ,GMD29659566118A, GMD65331541218A, GMD74701541318A, GMD67716462018A, GMD7716194717A, GMD59426461918A, GB32494283817 and GMD45321455618A

littmann.com.my