

3M™ V.A.C.® Therapy

3M™ ActiV.A.C.™ Therapy System



The 3M™ ActiV.A.C.™ Therapy Unit is portable Negative Pressure Wound Therapy (NPWT) designed for the ambulatory patient.

The features of this ergonomically designed therapy unit include:

For the patient:

- Lightweight. Weighs only 2.4 pounds
- Small size with a low profile that can be worn close to the body
- Easy-to-use, single-touch therapy on/off operation
- Alarm notifications that are easy to recognize and correct
- Easy, quick release 300 mL canister
- 14-hour battery for activities of daily living

For the clinician:

3M™ SensaT.R.A.C.™ Technology — only 3M™ V.A.C.® Therapy Units provide patented SensaT.R.A.C. Technology, a real time pressure feedback system.

- Adjusts pump output to the wound site, compensating for wound distance, wound position, exudate characteristics, and patient movement
- Delivers and maintains prescribed negative pressure for optimal healing outcomes

3M™ Seal Check™ Feature — designed to help clinicians and patients identify and troubleshoot negative pressure leaks.

- Instant feedback. Audio and visual cues are delivered in real-time, allowing easy location of dressing leaks
- Enhanced therapy confidence. A proper seal helps provide accurate delivery of prescribed therapy for optimal wound healing
- Possible reduction of unscheduled visits. Easier for clinicians to help troubleshoot problems over the phone

Settings guide — recommended therapy settings and pressure ranges are pre-programmed for indicated wound types:

- Easy setup. With pre-programmed settings by wound type, clinicians can initiate therapy quickly and efficiently
- Ease of use. Intuitive control functions make it easier to use, especially for clinicians who are not as familiar with V.A.C.® Therapy

Therapy history reporting — document patient compliance with V.A.C.® Therapy.

- Monitors therapy usage, settings, durations, and dressing and canister changes to help ensure patients are receiving the prescribed course of therapy
- USB port for convenient data downloads

Simplified touch screen — enhanced, simplified control functions:

- Easy navigation. With full-colour interface with intuitive touch screen controls
- Simplified patient mode. Designed to be easy to use and understand
- Easier training. With common control functions with the INFOV.A.C.™ Therapy System



3M™ V.A.C.® Dressing options with 3M™ SensaT.R.A.C.™ Technology.

Specially designed dressings, because no two wounds are alike.

3M™ SensaT.R.A.C.™ Pad

Only SensaT.R.A.C. Technology actively measures, monitors, and manages the accurate delivery of the physician's prescribed therapy to the wound site.

- Designed with patient comfort in mind
- Flexible pad material for easy application over body contours
- Discreet, low profile design



3M™ V.A.C.® Simplace™ Dressing

Designed to simplify the V.A.C.® Dressing placement process.

- Helps to promote healing by facilitating granulation tissue formation
- Fewer steps, faster application
- Spiral cut foam is simple to size; no scissors necessary
- Design allows for easier bridging



3M™ V.A.C.® Granufoam™ Bridge Dressing

Ideal for use with wounds in pressure sensitive areas that require bridging.

- Preassembled components specifically designed to simplify dressing application
- Integrated bridge allows for SensaT.R.A.C.™ Pad placement away from the wound site
- Wicking layers help ensure intact skin stays dry
- Allows V.A.C.® Therapy to be used in conjunction with off-loading therapies



3M™ V.A.C.® Granufoam™ Dressing

An advanced wound dressing to assist granulation tissue formation and enhance exudate removal.

- Can be easily trimmed to fit the contours of deep or irregularly shaped wounds
- Can be customized for bridging techniques when treating multiple wounds



3M™ V.A.C. Whitefoam™ Dressing

A comfortable and versatile wound dressing.

- May be used over split thickness skin grafts
- Helps promote flap and graft take
- Higher tensile strength (than V.A.C.® Granufoam™ Dressing) allows for easy removal from tunnels and undermining
- Recommended in situations where hypergranulation responses are likely



3M™ V.A.C.® Peel and Place Dressing Kit

All-in-one design makes application easy.

- The perforated, non-adherent layer allows for extended wear time up to 7-days¹ and fewer dressing changes
- Patented, all-in-one dressing allows for easy, time-saving application with less training
- The simple, versatile design of V.A.C.® Peel and Place Dressing is compatible with multiple V.A.C.® Therapy Units



Ordering Information

3M™ V.A.C.® Dressings with 3M™ SensaT.R.A.C.™ Pad

Part Number	Description	Case Quantity
Up to 7-day wear time dressing kits for V.A.C.® Therapy*		
EZ5SML/5	3M™ V.A.C.® Peel and Place Dressing Kit - Small	5
EZ5MED/5	3M™ V.A.C.® Peel and Place Dressing Kit - Medium	5
EZ5LRG/5	3M™ V.A.C.® Peel and Place Dressing Kit - Large	5
Up to 48-72 hours, or no less than 3 times per week dressing kits for V.A.C.® Therapy		
M8275051/10	3M™ V.A.C.® Granufoam™ Dressing Kit, Small	10
M8275051/5	3M™ V.A.C.® Granufoam™ Dressing Kit, Small	5
M8275052/10	3M™ V.A.C.® Granufoam™ Dressing Kit, Medium	10
M8275052/5	3M™ V.A.C.® Granufoam™ Dressing Kit, Medium	5
M8275053/10	3M™ V.A.C.® Granufoam™ Dressing Kit, Large	10
M8275053/5	3M™ V.A.C.® Granufoam™ Dressing Kit, Large	5
M8275065/5	3M™ V.A.C.® Granufoam™ Dressing Kit, X-Large	5
M8275045/5	3M™ V.A.C.® Simplace™ Dressing Kit, Medium	5
M8275046/5	3M™ V.A.C.® Simplace™ Dressing Kit, Small	5
M8275042/10	3M™ V.A.C.® Granufoam™ Bridge Dressing	10
M8275042/5	3M™ V.A.C.® Granufoam™ Bridge Dressing	5
M8275044/5	3M™ V.A.C.® Granufoam™ Bridge XG Dressing	5
M6275033/10	3M™ V.A.C. Whitefoam™ Dressing, Small, Foam Only	10
M6275034/10	3M™ V.A.C. Whitefoam™ Dressing, Large, Foam Only	10
M8275067/10	3M™ V.A.C. Whitefoam™ Dressing Kit, Large	10
M8275068/10	3M™ V.A.C. Whitefoam™ Dressing Kit, Small	10
Accessories for V.A.C.® Therapy		
M6275009/10	3M™ V.A.C.® Drape	10
M6275066/10	3M™ V.A.C.® Y-Connector	10
M6275069/10	3M™ V.A.C.® Tubing Cap	10
3M™ ActiV.A.C.™ Canisters		
M8275058/10	3M™ ActiV.A.C.™ Canister — 300 ml with Gel	10
M8275058/5	3M™ ActiV.A.C.™ Canister — 300 ml with Gel	5

*See V.A.C.® Peel and Place Dressing IFU for more information on wounds appropriate for these dressings.

For more information about 3M™ V.A.C.® Therapy, please call **1-800-668-5403** or visit: **3M.ca/Medical**

Note: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product Instructions for Use prior to application. This material is intended for healthcare professionals.

Note: Disposable components of the V.A.C.® (Vacuum Assisted Closure®) Therapy System, including the foam dressing (i.e., V.A.C.® Granufoam™, V.A.C. Granufoam Silver®, or V.A.C. Whitefoam™ and for Granufoam is V.A.C.® Granufoam™, tubing and drape are packaged sterile and are latex-free. V.A.C.® Therapy Unit canisters are packaged fluid path sterile and are latex free. All disposable components of the V.A.C.® Therapy System are for single use only. To help ensure safe and effective use, the V.A.C.® Granufoam®, V.A.C. Granufoam Silver® and V.A.C.® WhiteFoam Dressings are to be used only with V.A.C.® Therapy Units.

Reference

1. SAT-BSER-05-869347 VAC Peel and Place (Ganymede) BSER



Available in Canada from:

KCI Medical Canada, Inc.
75 Courtneypark Drive W., Unit #4
Mississauga
ON L5W 0E3
Canada
Phone: 1-800-668-5403
Web: 3M.ca/Medical

KCI USA, INC.
12930 IH 10 West
San Antonio, TX 78249

© Solventum 2024. Solventum, the S logo and other trademarks are trademarks of Solventum or its affiliates.
CAN_EN_70-2013-1391-6 2206-24256-E