

3M Science.
Applied to Life.™



**3M™ Bair Hugger™
Blanket Selection Guide**

Designed to be different

Effectiveness of forced-air warming systems depends on blanket design, particularly the evenness of heat distribution across the blanket^{1,2}

Moving air loses heat very quickly.^{1,2} Therefore air should be channeled efficiently to all parts of a blanket to provide even heat distribution.



Interconnected Air Channel System

Interconnected air channels run the length of the blanket providing even distribution of warm air to the patient.



Air Hole Perforation

A symmetrical pattern of small holes extends across the entire surface of the blankets to maximise recruitment of skin surface area.




Pre-sealed Hose Port

In selected blankets, two re-sealable hose ports enable versatile and convenient positioning of the warming unit.



Fluid Outlets


In selected blankets, fluid outlets minimise pooling of fluids on the surface of the blanket and cooling effect of excess fluids.



**Click on the buttons to view
a selection of 3M™ warming blankets
to suit your surgical needs.**



Click on the buttons to find out more



Head & Neck

**3M™ Bair Hugger™
Lower Body Warming Blanket**

52500


Pre-op
& Post-op

Warming
Gown

Paediatrics

Click on the buttons to find out more

Arm & Shoulder




**3M™ Bair Hugger™
Lower Body Warming Blanket**
52500



**3M™ Bair Hugger™ Multi-position
Upper Body Warming Blanket**
62200

- Pre-op & Post-op
- Warming Gown
- Paediatrics

Click on the buttons to find out more



- Pre-op & Post-op
- Warming Gown
- Paediatrics

Click on the buttons to find out more

Hand




3M™ Bair Hugger™ Multi-position Upper Body Warming Blanket
62200



3M™ Bair Hugger™ Lower Body Warming Blanket
52500



3M™ Bair Hugger™ Full Access Underbody Warming Blanket
63500



- Pre-op & Post-op
- Warming Gown
- Paediatrics

Click on the buttons to find out more

Torso



3M™ Bair Hugger™ Adult Underbody Warming Blanket
54500



3M™ Bair Hugger™ Full Access Underbody Warming Blanket
63500




3M™ Bair Hugger™ Lower Body Warming Blanket
52500



3M™ Bair Hugger™ Sterile Cardiac Warming Blanket
63000

Abdomen



- Pre-op
& Post-op
- Warming
Gown
- Paediatrics

Click on the buttons to find out more



3M™ Bair Hugger™ Lower Body Warming Blanket
52500



3M™ Bair Hugger™ Surgical Access Warming Blanket
57000




3M™ Bair Hugger™ Adult Underbody Warming Blanket
54500



3M™ Bair Hugger™ Full Access Underbody Warming Blanket
63500



3M™ Bair Hugger™ Multi-position Upper Body Warming Blanket
62200



- Pre-op & Post-op
- Warming Gown
- Paediatrics

Click on the buttons to find out more

Uro & Rectal



3M™ Bair Hugger™ Lithotomy Underbody Warming Blanket

58501




3M™ Bair Hugger™ Dual Port Torso Warming Blanket

54200



3M™ Bair Hugger™ Multi-position Upper Body Warming Blanket

62200



- Pre-op & Post-op
- Warming Gown
- Paediatrics

Click on the buttons to find out more

Hip



3M™ Bair Hugger™ Surgical Access Warming Blanket
57000




3M™ Bair Hugger™ Multi-position Upper Body Warming Blanket
62200



3M™ Bair Hugger™ Full Access Underbody Warming Blanket
63500



3M™ Bair Hugger™ Dual Port Torso Warming Blanket
54200



Foot & Knee

**3M™ Bair Hugger™ Multi-position
Upper Body Warming Blanket**

62200



**3M™ Bair Hugger™ Dual Port Torso
Warming Blanket**


54200

Pre-op
& Post-op

Warming
Gown

Paediatrics

Click on the buttons to find out more



- Pre-op & Post-op
- Warming Gown
- Paediatrics

Click on the buttons to find out more

Spine & Neuro



3M™ Bair Hugger™ Full Access Underbody Warming Blanket
63500




3M™ Bair Hugger™ Spinal Underbody Warming Blanket
57501



3M™ Bair Hugger™ Lower Body Warming Blanket
52500



3M™ Bair Hugger™ Surgical Access Warming Blanket
57000



Pre-op & Post-op



3M™ Bair Hugger™ Full Body Warming Blanket
30000



3M™ Bair Hugger™ Multi-Access Warming Blanket
31500


- Pre-op & Post-op**
- Warming Gown
- Paediatrics

Click on the buttons to find out more



Warming Gown


3M™ Bair Hugger™ Plus Gown
81002 & 81202



3M™ Bair Hugger™ Flex Gown
81003 & 81203



Click on the buttons to find out more



Click on the buttons to find out more

Paediatric Warming Blankets




3M™ Bair Hugger™ Paediatric Underbody Warming Blanket
55501



3M™ Bair Hugger™ Small Lower Body Warming Blanket
53700



3M™ Bair Hugger™ Large Paediatric Underbody Warming Blanket
55000



3M Australia Pty Limited
Building A, 1 Rivett Rd
North Ryde, NSW 2113
1300 363 454
www.3M.com.au

3M New Zealand Limited
94 Apollo Drive
Rosedale, Auckland 0632
0800 808 182
www.3M.com.nz

3M, Ranger, Bair Hugger and the other marks shown are marks and/or registered marks. Unauthorised use prohibited. All other marks are property of their respective owners. Please recycle. © 3M 2023. All rights reserved.