




| | |
|--|-------------|
| 3M™ Victory Series™ Superior Fit Buccal Tubes | 4.2 |
| 3M™ Victory Series™ Buccal Tubes | 4.4 |
| 3M™ Victory Series™ Molar Bonding Bases | 4.5 |
| 3M™ Victory Series™ First Molar Buccal Tubes..... | 4.6 |
| 3M™ Victory Series™ Second Molar Buccal Tubes... | 4.13 |
| 3M™ Victory Series™ Second Molar Mini Tubes | 4.15 |
| 3M™ Victory Series™ 2nd Bicuspid Tubes..... | 4.15 |
| Light Wire Buccal Tubes | 4.16 |

Victory Series™ Superior Fit Buccal Tubes

Superior fit.
Outstanding comfort.
Low profile tubes.



Victory Series Superior Fit Buccal Tubes available with:

APC™ Flash-Free Adhesive


Also available with APC™ II Adhesive pre-coating.

MBT™ Versatile+ Appliance System

Victory Series™ Superior Fit Buccal Tubes

Designed with extensive user input, Victory Series Superior Fit Buccal Tubes take advantage of modern computer design tools to offer performance enhancements you can both see and feel when you place them.

Well thought-out features improve gripping, handling and visual cues during placement. A large compound-contour base results in improved ease of tube positioning. Inside, Victory Series Superior Fit Buccal Tubes feature an exclusive “contoured funnel” design for easy wire insertion, designed to eliminate sharp corners and edges that contribute to binding and notching of the archwire. A low profile flush-mount hook and sleek tapered body add to patient comfort, and available APC™ Adhesive pre-coating make Victory Series Superior Fit Buccal Tubes an ideal choice.



MBT™ System Prescription

| Upper | | | Tube Length | | 0.018 in. | | 0.022 in. | |
|-----------|------------------------|----------|-------------|-----|-----------|----------|-----------|----------|
| Tooth | Type | T/O | inches | mm | L | R | L | R |
| 1st Molar | Single Convertible | -14T/10° | .170 | 4.3 | 069-5303 | 069-5314 | 069-6303 | 069-6314 |
| 1st Molar | Single Non-Convertible | -14T/10° | .187 | 4.7 | 069-5103 | 069-5114 | 069-6103 | 069-6114 |
| 1st Molar | Double Convertible | -14T/10° | .170 | 4.3 | 069-5153 | 069-5164 | 069-6153 | 069-6164 |
| 2nd Molar | Single Non-Convertible | -14T/10° | .147 | 3.7 | 069-3103 | 069-3114 | 069-4103 | 069-4114 |
| Lower | | | Tube Length | | 0.018 in. | | 0.022 in. | |
| Tooth | Type | T/O | inches | mm | L | R | L | R |
| 1st Molar | Single Convertible | -20T/0° | .170 | 4.3 | 069-5403 | 069-5414 | 069-6403 | 069-6414 |
| 1st Molar | Single Non-Convertible | -20T/0° | .187 | 4.7 | 069-5203 | 069-5214 | 069-6203 | 069-6214 |
| 1st Molar | Double Convertible | -20T/0° | .170 | 4.3 | 069-5253 | 069-5264 | 069-6253 | 069-6264 |
| 2nd Molar | Single Non-Convertible | -15T/0° | .147 | 3.7 | 069-3203 | 069-3214 | 069-4203 | 069-4214 |

Adding Prefix “4” is APC™ Flash-Free Adhesive System – Example: 4069-5303

Adding Prefix “3” is APC™ II Adhesive System – Example: 3069-5303

Check with 3M for availability of APC™ Adhesive Coating.

Roth* Prescription

| Upper | | | Tube Length | | 0.018 in. | | 0.022 in. | |
|-----------|------------------------|---------|-------------|-----|-----------|----------|-----------|----------|
| Tooth | Type | T/O | inches | mm | L | R | L | R |
| 1st Molar | Single Convertible | -10T/7° | .170 | 4.3 | 069-5313 | 069-5324 | 069-6313 | 069-6324 |
| 1st Molar | Single Non-Convertible | -10T/7° | .187 | 4.7 | 069-5113 | 069-5124 | 069-6113 | 069-6124 |
| 1st Molar | Double Convertible | -10T/7° | .170 | 4.3 | 069-5163 | 069-5174 | 069-6163 | 069-6174 |
| 2nd Molar | Single Non-Convertible | -10T/6° | .147 | 3.7 | 069-3113 | 069-3124 | 069-4113 | 069-4124 |
| Lower | | | Tube Length | | 0.018 in. | | 0.022 in. | |
| Tooth | Type | T/O | inches | mm | L | R | L | R |
| 1st Molar | Single Convertible | -25T/6° | .170 | 4.3 | 069-5413 | 069-5424 | 069-6413 | 069-6424 |
| 1st Molar | Single Non-Convertible | -25T/6° | .187 | 4.7 | 069-5213 | 069-5224 | 069-6213 | 069-6224 |
| 1st Molar | Double Convertible | -25T/6° | .170 | 4.3 | 069-5263 | 069-5274 | 069-6263 | 069-6274 |
| 2nd Molar | Single Non-Convertible | -30T/6° | .147 | 3.7 | 069-3213 | 069-3224 | 069-4213 | 069-4224 |

*3M version of this prescription. No endorsement by the Doctor is implied.

Adding Prefix “4” is APC™ Flash-Free Adhesive System – Example: 4069-5113

Adding Prefix “3” is APC™ II Adhesive System – Example: 3069-5113

Check with 3M on local availability.

Victory Series™ Buccal Tubes

Our compact, technically advanced Victory Series Buccal Tubes include a trouble-free integral convertible cap that converts easily and consistently, every time. All you have to do is lift and peel. This exclusive cap also offers solid, one piece construction, eliminating the possibility of braze or solder materials that can leave a residue in the archwire slot.

Victory Series Buccal Tubes are designed with funneled and radiused entries on both the headgear tube and archwire slot, resulting in easy archwire insertion. Smooth, patient-friendly hooks improve comfort for your patients. You'll also appreciate the strong construction, low-profile design and distal notch that secures tie backs. Increasing patient comfort and providing you with a well-designed, reliable and predictable orthodontic solution.



Victory Series™ Molar Bonding Bases

Now you can achieve a stronger, more consistent bond in every case with the Victory Series Molar Bonding Base. Look closely at its profile and you'll see a curved anatomical design for optimal tube-to-tooth fit and a buccal indent that naturally conforms to the surface of the molar, making it easier to correctly position the buccal tube in place.

The micro-etched 80 gauge mesh, along with the buccal indent and the correct anatomical fit, are designed for increased bond strength. The buccal indent assists in more accurate tube positioning; a key advantage that can eliminate time-consuming archwire adjustments. And with the Victory Series Molar Bonding base, you'll never have to worry about having to use tooth separators for band placement.

For added convenience, two different sizes accommodate a wide array of buccal attachments.



Victory Series Molar Bonding Bases available with:

APC™ II Adhesive

Also available with APC™ PLUS Adhesive pre-coating.

MBT™ Versatile+ Appliance System

Victory Series™ First Molar Buccal Tubes

Upper Convertible Double Combination



Compatible with:

**Forsus™ Fatigue Resistant
Device EZ2 Module**

| Slot Torque | | Distal Offset | Headgear | Length | L/R | .018 Slot | .022 Slot |
|-------------|--|------------------|----------|--------|-----|--------------|--------------|
| Arch | | | | | | | |
| 0° | | 0° | .045 | 3.6 mm | UL | 068-841 | 068-941 |
| 0° | | 0° | .045 | 3.6 mm | UR | 068-842 | 068-942 |
| 0° | | 0° | .051 | 3.6 mm | UL | 068-843 | 068-943 |
| 0° | | 0° | .051 | 3.6 mm | UR | 068-844 | 068-944 |
| 0° | | 7° | .045 | 4.3 mm | UL | 068-801 | 068-901 |
| 0° | | 7° | .045 | 4.3 mm | UR | 068-802 | 068-902 |
| 0° | | 7° | .051 | 4.3 mm | UL | 068-803 | 068-903 |
| 0° | | 7° | .051 | 4.3 mm | UR | 068-804 | 068-904 |
| 0° | | 7° | .045 | 3.6 mm | UL | 068-875 | 068-975 |
| 0° | | 7° | .045 | 3.6 mm | UR | 068-874 | 068-974 |
| -10° | | 7° | .045 | 4.3 mm | UL | 068-805 | 068-905 |
| -10° | | 7° | .045 | 4.3 mm | UR | 068-806 | 068-906 |
| -10° | | 7° | .051 | 4.3 mm | UL | 068-807 | 068-907 |
| -10° | | 7° | .051 | 4.3 mm | UR | 068-808 | 068-908 |
| -10° | | 7° | .045 | 3.6 mm | UL | 068-871 | 068-971 |
| -10° | | 7° | .045 | 3.6 mm | UR | 068-870 | 068-970 |
| -10° | | 7° | .051 | 3.6 mm | UL | 068-873 | 068-973 |
| -10° | | 7° | .051 | 3.6 mm | UR | 068-872 | 068-972 |
| -10° | | 14° | .045 | 4.3 mm | UL | 068-809 | 068-909 |
| -10° | | 14° | .045 | 4.3 mm | UR | 068-810 | 068-910 |
| -10° | | 14° | .051 | 4.3 mm | UL | 068-811 | 068-911 |
| -10° | | 14° | .051 | 4.3 mm | UR | 068-812 | 068-912 |
| -10° | | 14° | .045 | 3.6 mm | UL | 068-879 | 068-979 |
| -10° | | 14° | .045 | 3.6 mm | UR | 068-878 | 068-978 |

Upper Convertible Double Combination - MBT™ Versatile+ Appliance System



Compatible with:

**Forsus™ Fatigue Resistant
Device EZ2 Module**

| Slot Torque | | Distal Offset | Headgear | Length | L/R | .018 Slot | .022 Slot | MBT™ |
|-------------|--|------------------|----------|--------|-----|--------------|--------------|------|
| Arch | | | | | | | | |
| -14° | | 10° | .045 | 4.6 mm | UL | 068-893 | 068-993 | MBT™ |
| -14° | | 10° | .045 | 4.6 mm | UR | 068-894 | 068-994 | MBT™ |
| -14° | | 10° | .051 | 4.6 mm | UL | 068-895 | 068-995 | MBT™ |
| -14° | | 10° | .051 | 4.6 mm | UR | 068-896 | 068-996 | MBT™ |

Victory Series™ First Molar Buccal Tubes

Upper Convertible Double Combination With Gingival Headgear Tube



| Slot Torque | | Distal Offset | Headgear | Length | L/R | .018 Slot | .022 Slot | |
|-------------|--|---------------|----------|--------|-----|-----------|-----------|------|
| Arch | | | | | | | | |
| -10° | | 7° | .045 | 4.3 mm | UL | 068-821 | 068-921 | |
| -10° | | 7° | .045 | 4.3 mm | UR | 068-822 | 068-922 | |
| -14° | | 10° | .045 | 4.3 mm | UL | 068-885 | 068-985 | MBT™ |
| -14° | | 10° | .045 | 4.3 mm | UR | 068-886 | 068-986 | MBT™ |
| -14° | | 10° | .051 | 4.3 mm | UL | 068-887 | 068-987 | MBT™ |
| -14° | | 10° | .051 | 4.3 mm | UR | 068-888 | 068-988 | MBT™ |

Upper Convertible Double



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | |
|-------------|-----|---------------|--------|-----|-----------|---|-----------|---|
| Arch | Aux | | | | | | | |
| 0° | 0° | 0° | 3.6 mm | UL | 068-849 | 1 | 068-949 | 1 |
| 0° | 0° | 0° | 3.6 mm | UR | 068-850 | 2 | 068-950 | 2 |
| -10° | 0° | 7° | 3.6 mm | UL | 068-851 | 1 | 068-951 | 1 |
| -10° | 0° | 7° | 3.6 mm | UR | 068-852 | 2 | 068-952 | 2 |

Base

Adding Prefix "3" is APC™ II Adhesive System – Example: 3068-849

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5068-849

Upper Convertible Double Differentiated Slot



| Slot Torque | | Distal Offset | Main | Aux | Length | L/R | .018 Slot | .022 Slot |
|-------------|-----|---------------|------|------|--------|-------|-----------|-----------|
| Arch | Aux | | | | | | | |
| 0° | 0° | 0° | .022 | .018 | 3.6 mm | UL/LR | - | 068-953 |
| 0° | 0° | 0° | .022 | .018 | 3.6 mm | UR/LL | - | 068-954 |

Note

| | |
|--|--|
| 3 Use this prefix to order APC™ II Adhesive Coated Appliance System. | 5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive. |
| APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually. | Measurements are in millimeters unless otherwise noted. |
| 1 2 3 4 | Buccal Tube is available with a bondable base. |

Victory Series™ First Molar Buccal Tubes

Upper Convertible Triple Combination - Occlusal Headgear



Compatible with:

**Forsus™ Fatigue Resistant
Device EZ2 Module**

| Slot Torque | | Distal Offset | Headgear | Length | L/R | .018 Slot | .022 Slot | |
|-------------|-----|------------------|----------|--------|-----|--------------|--------------|------|
| Arch | Aux | | | | | | | |
| 0° | 0° | 0° | .045 | 3.6 mm | UL | 068-845 | 068-945 | |
| 0° | 0° | 0° | .045 | 3.6 mm | UR | 068-846 | 068-946 | |
| 0° | 0° | 0° | .051 | 3.6 mm | UL | 068-847 | 068-947 | |
| 0° | 0° | 0° | .051 | 3.6 mm | UR | 068-848 | 068-948 | |
| 0° | 0° | 15° | .045 | 4.0 mm | UL | 068-825 | - | |
| 0° | 0° | 15° | .045 | 4.0 mm | UR | 068-826 | - | |
| -10° | 0° | 7° | .045 | 3.6 mm | UL | 068-813 | 068-913 | |
| -10° | 0° | 7° | .045 | 3.6 mm | UR | 068-814 | 068-914 | |
| -10° | 0° | 7° | .051 | 3.6 mm | UL | 068-815 | 068-915 | |
| -10° | 0° | 7° | .051 | 3.6 mm | UR | 068-816 | 068-916 | |
| -10° | 0° | 14° | .045 | 3.6 mm | UL | 068-817 | 068-917 | |
| -10° | 0° | 14° | .045 | 3.6 mm | UR | 068-818 | 068-918 | |
| -10° | 0° | 14° | .051 | 3.6 mm | UL | 068-819 | 068-919 | |
| -10° | 0° | 14° | .051 | 3.6 mm | UR | 068-820 | 068-920 | |
| -14° | 0° | 10° | .045 | 4.3 mm | UL | 068-889 | 068-989 | MBT™ |
| -14° | 0° | 10° | .045 | 4.3 mm | UR | 068-890 | 068-990 | MBT™ |
| -14° | 0° | 10° | .051 | 4.3 mm | UL | 068-891 | 068-991 | MBT™ |
| -14° | 0° | 10° | .051 | 4.3 mm | UR | 068-892 | 068-992 | MBT™ |

MBT™
MBT™
MBT™
MBT™

Upper Convertible Triple Combination Differentiated Slot



Compatible with:

**Forsus™ Fatigue Resistant
Device EZ2 Module**

| Slot Torque | | Distal Offset | Headgear | Main | Aux | Length | L/R | .018 Slot | .022 Slot |
|-------------|-----|------------------|----------|------|------|--------|-----|--------------|--------------|
| Arch | Aux | | | | | | | | |
| -10° | 0° | 14° | .045 | .018 | .022 | 3.6 mm | UL | 068-853 | - |
| -10° | 0° | 14° | .045 | .018 | .022 | 3.6 mm | UR | 068-854 | - |
| -10° | 0° | 7° | .045 | .022 | .018 | 3.6 mm | UL | - | 068-959 |
| -10° | 0° | 7° | .045 | .022 | .018 | 3.6 mm | UR | - | 068-960 |
| 0° | 0° | 0° | .045 | .022 | .018 | 3.6 mm | UL | - | 068-957 |
| 0° | 0° | 0° | .045 | .022 | .018 | 3.6 mm | UR | - | 068-958 |

Victory Series™ First Molar Buccal Tubes

Upper Convertible Single



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | |
|-------------|--|---------------|--------|-----|-----------|------------|-----------|------------|
| Arch | | | | | Base | Jumbo Base | Base | Jumbo Base |
| 0° | | 0° | 3.6 mm | UL | 067-808 | 1 | 067-908 | 1 |
| 0° | | 0° | 3.6 mm | UR | 067-807 | 2 | 067-907 | 2 |
| 0° | | 6° | 3.6 mm | UL | 067-810 | 1 | 067-910 | 1 |
| 0° | | 6° | 3.6 mm | UR | 067-809 | 2 | 067-909 | 2 |
| -10° | | 7° | 4.3 mm | UL | 068-823 | 1 | 068-923 | 1 |
| -10° | | 7° | 4.3 mm | UR | 068-824 | 2 | 068-924 | 2 |

Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-808

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-808

Upper Convertible Single - MBT™ Versatile+ Appliance System



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | | MBT™ |
|-------------|--|---------------|--------|-----|-----------|------------|-----------|------------|------|
| Arch | | | | | Base | Jumbo Base | Base | Jumbo Base | |
| -14° | | 10° | 4.0 mm | UL | 068-897 | 1 | 068-997 | 1 | MBT™ |
| -14° | | 10° | 4.0 mm | UR | 068-898 | 2 | 068-998 | 2 | MBT™ |

Adding Prefix "3" is APC™ II Adhesive System – Example: 3068-897

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5068-897

Upper 1ST Non-Convertible Single - Funneled



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | | MBT™ |
|-------------|--|---------------|--------|-----|-----------|------------|-----------|------------|------|
| Arch | | | | | Base | Jumbo Base | Base | Jumbo Base | |
| -14° | | 10° | 4.6 mm | UL | 067-879 | 1 | 067-979 | 1 | MBT™ |
| -14° | | 10° | 4.6 mm | UR | 067-878 | 2 | 067-978 | 2 | MBT™ |

Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-879

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-879

Note

| | |
|--|--|
| 3 Use this prefix to order APC™ II Adhesive Coated Appliance System. | 5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive. |
| APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually. | Measurements are in millimeters unless otherwise noted. |
| 1 2 3 4 Buccal Tube is available with a bondable base. | 8 Use this prefix to order the large contoured Jumbo Base. |


Victory Series™ First Molar Buccal Tubes

Upper Non-Convertible Double Combination



| Slot Torque | | Distal Offset | Headgear | Length | L/R | .018 Slot | .022 Slot |
|-------------|--|---------------|----------|--------|-----|-----------|-----------|
| Arch | | | | | | | |
| 0° | | 0° | .045 | 5.1 mm | UL | 068-861 | 068-961 |
| 0° | | 0° | .045 | 5.1 mm | UR | 068-862 | 068-962 |
| 0° | | 0° | .051 | 5.1 mm | UL | 068-863 | 068-963 |
| 0° | | 0° | .051 | 5.1 mm | UR | 068-864 | 068-964 |
| 0° | | 6° | .045 | 5.1 mm | UL | - | 068-965 |
| 0° | | 6° | .045 | 5.1 mm | UR | - | 068-966 |
| 0° | | 6° | .051 | 5.1 mm | UL | 068-867 | 068-967 |
| 0° | | 6° | .051 | 5.1 mm | UR | 068-868 | 068-968 |

Upper Non-Convertible Triple Combination



Compatible with:

Forsus™ Fatigue Resistant Device EZ2 Module

| Slot Torque | | Distal Offset | Headgear | Hook | Length | L/R | .018 Slot | .022 Slot | |
|-------------|-----|---------------|----------|------|--------|------|-----------|-----------|----|
| Arch | Aux | | | | | | | | |
| 0° | 0° | 0° | .045 | CTR | 4.3 mm | UNIV | 068-882 | 068-982 | * |
| 0° | 0° | 0° | .051 | CTR | 4.3 mm | UNIV | 068-883 | 068-983 | * |
| 0° | 0° | 0° | .045 | MES | 4.3 mm | UL | 068-601 | - | ** |
| 0° | 0° | 0° | .045 | MES | 4.3 mm | UR | 068-602 | - | ** |

Lower Convertible Single



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | .022 Slot |
|-------------|--|---------------|--------|-----|-----------|-----------|
| Arch | | | | | | |
| 0° | | 0° | 3.6 mm | LL | 067-807 | 067-907 |
| 0° | | 0° | 3.6 mm | LR | 067-808 | 067-908 |
| 0° | | 6° | 3.6 mm | LL | 067-809 | 067-909 |
| 0° | | 6° | 3.6 mm | LR | 067-810 | 067-910 |
| -10° | | 0° | 4.3 mm | LR | 067-813 | 067-913 |
| -10° | | 0° | 4.3 mm | LL | 067-814 | 067-914 |
| -25° | | 6° | 3.6 mm | LL | 067-851 | 067-951 |
| -25° | | 6° | 3.6 mm | LR | 067-850 | 067-950 |
| -25° | | 6° | 4.3 mm | LL | 067-801 | 067-901 |
| -25° | | 6° | 4.3 mm | LR | 067-802 | 067-902 |
| -20° | | 0° | 4.3 mm | LL | 067-819 | 067-919 |
| -20° | | 0° | 4.3 mm | LR | 067-820 | 067-920 |
| -25° | | 0° | 4.3 mm | LL | 067-817 | 067-917 |
| -25° | | 0° | 4.3 mm | LR | 067-818 | 067-918 |

MBT™
MBT™

Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-807

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-807

Victory Series™ First Molar Buccal Tubes

Lower 1ST Convertible Single Low Profile



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | | .022 Slot | | | MBT™ |
|-------------|--|---------------|--------|-----|-----------|---|---|-----------|---|---|------|
| Arch | | | | | 3 | 8 | 4 | 8 | 3 | 8 | |
| -20° | | 0° | 4.3 mm | LL | 067-520 | 3 | 8 | 067-620 | 3 | 8 | MBT™ |
| -20° | | 0° | 4.3 mm | LR | 067-521 | 4 | 8 | 067-621 | 4 | 8 | MBT™ |




Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-520

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-520

Lower 1ST Non-Convertible Single - Funneled



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | | .022 Slot | | | MBT™ |
|-------------|--|---------------|--------|-----|-----------|---|---|-----------|---|---|------|
| Arch | | | | | 3 | 8 | 4 | 8 | 3 | 8 | |
| -20° | | 0° | 4.6 mm | LL | 067-877 | 3 | 8 | 067-977 | 3 | 8 | MBT™ |
| -20° | | 0° | 4.6 mm | LR | 067-876 | 4 | 8 | 067-976 | 4 | 8 | MBT™ |



Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-877

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-877

Lower Non-Convertible Single



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | | MBT™ |
|-------------|--|---------------|--------|-----|-----------|---|-----------|---|------|
| Arch | | | | | 3 | 4 | 3 | 4 | |
| -20° | | 0° | 4.3 mm | LL | 067-867 | 3 | 067-967 | 3 | MBT™ |
| -20° | | 0° | 4.3 mm | LR | 067-868 | 4 | 067-968 | 4 | MBT™ |




Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-867

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-867

Lower Convertible Double



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | | MBT™ |
|-------------|-----|---------------|--------|-----|-----------|---|-----------|---|------|
| Arch | Aux | | | | 3 | 4 | 3 | 4 | |
| 0° | 0° | 0° | 3.6 mm | LL | 068-850 | 3 | 068-950 | 3 | MBT™ |
| 0° | 0° | 0° | 3.6 mm | LR | 068-849 | 4 | 068-949 | 4 | |
| -25° | 0° | 7° | 3.6 mm | LL | 067-803 | 3 | 067-903 | 3 | MBT™ |
| -25° | 0° | 7° | 3.6 mm | LR | 067-804 | 4 | 067-904 | 4 | |
| -20° | 0° | 0° | 4.3 mm | LL | 067-861 | 3 | 067-961 | 3 | MBT™ |
| -20° | 0° | 0° | 4.3 mm | LR | 067-860 | 4 | 067-960 | 4 | |



Adding Prefix "3" is APC™ II Adhesive System – Example: 3068-850

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5068-850

Note


| | |
|--|--|
| 3 Use this prefix to order APC™ II Adhesive Coated Appliance System. | 5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive. |
| APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually. | Measurements are in millimeters unless otherwise noted. |
| 1 2 3 4 Buccal Tube is available with a bondable base. | 8 Use this prefix to order the large contoured Jumbo Base. |

Victory Series™ First Molar Buccal Tubes

Lower Convertible Double Differentiated Slot

| Slot Torque | | Distal Offset | Main | Aux | Length | L/R | .018 Slot | .022 Slot |
|-------------|-----|---------------|------|------|--------|-----|-----------|-----------|
| Arch | Aux | | | | | | | |
| -25° | 0° | 7° | .022 | .018 | 3.6 mm | LL | - | 067-955 |
| -25° | 0° | 7° | .022 | .018 | 3.6 mm | LR | - | 067-956 |

Lower Lip Bumper Convertible Double



| Slot Torque | | Distal Offset | Lip Bumper | Length | L/R | .018 Slot | .022 Slot |
|-------------|-----|---------------|------------|--------|-----|-----------|-----------|
| Arch | Aux | | | | | | |
| -25° | | 5° | .045 | 4.3 mm | LL | 061-251 | 061-351 |
| -25° | | 5° | .045 | 4.3 mm | LR | 061-252 | 061-352 |

Lower Lip Bumper Convertible Triple



| Slot Torque | | Distal Offset | Lip Bumper | Length | L/R | .018 Slot | .022 Slot |
|-------------|-----|---------------|------------|--------|-----|-----------|-----------|
| Arch | Aux | | | | | | |
| -25° | 0° | 5° | .045 | 4.3 mm | LL | 067-811 | 067-911 |
| -25° | 0° | 5° | .045 | 4.3 mm | LR | 067-812 | 067-912 |

Upper/Lower Non-Convertible Double



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | .022 Slot |
|-------------|-----|---------------|--------|------|-----------|-----------|
| Arch | Aux | | | | | |
| 0° | 0 | 0° | 4.3 mm | UNIV | 067-843 | 067-943 |

Victory Series™ Second Molar Buccal Tubes



Upper Second Non-Convertible Single - Funneled

| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | | MBT™/Sondhi |
|-------------|--|---------------|--------|-----|-----------|---|-----------|---|-------------|
| Arch | | | | | | | | | |
| -19° | | 8° | 3.9 mm | UL | 067-873 | 1 | 067-973 | 1 | MBT™/Sondhi |
| -19° | | 8° | 3.9 mm | UR | 067-872 | 2 | 067-972 | 2 | MBT™/Sondhi |

Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-873

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-873

Upper Second Non-Convertible Single with Distal Notch



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | | MBT™ |
|-------------|--|---------------|--------|-----|-----------|---|-----------|---|------|
| Arch | | | | | | | | | |
| -10° | | 6° | 3.6 mm | UL | 067-829 | 1 | 067-929 | 1 | MBT™ |
| -10° | | 6° | 3.6 mm | UR | 067-830 | 2 | 067-930 | 2 | |
| -10° | | 10° | 3.6 mm | UL | 067-831 | 1 | 067-931 | 1 | |
| -10° | | 10° | 3.6 mm | UR | 067-832 | 2 | 067-932 | 2 | |
| -14° | | 10° | 4.3 mm | UL | 067-857 | 1 | 067-957 | 1 | |
| -14° | | 10° | 4.3 mm | UR | 067-858 | 2 | 067-958 | 2 | |

Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-829

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-829

Upper/Lower Second Non-Convertible Single with Distal Notch

| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | | .022 Slot | |
|-------------|--|---------------|--------|-----|-----------|---|-----------|---|
| Arch | | | | | | | | |
| 0° | | 0° | 3.6 mm | UL | 067-825 | 1 | 067-925 | 1 |
| 0° | | 0° | 3.6 mm | UR | 067-826 | 2 | 067-926 | 2 |
| 0° | | 0° | 3.6 mm | LL | 067-826 | 3 | 067-926 | 3 |
| 0° | | 0° | 3.6 mm | LR | 067-825 | 4 | 067-925 | 4 |
| 0° | | 6° | 3.6 mm | UL | 067-827 | 1 | 067-927 | 1 |
| 0° | | 6° | 3.6 mm | UR | 067-828 | 2 | 067-928 | 2 |
| 0° | | 6° | 3.6 mm | LL | 067-828 | 3 | 067-928 | 3 |
| 0° | | 6° | 3.6 mm | LR | 067-827 | 4 | 067-927 | 4 |

Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-825

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-825

Note

| | |
|--|--|
| 3 Use this prefix to order APC™ II Adhesive Coated Appliance System. | 5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive. |
| APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually. | Measurements are in millimeters unless otherwise noted. |
| 1 2 3 4 Buccal Tube is available with a bondable base. | |

Victory Series™ Second Molar Buccal Tubes



Lower Second Non-Convertible Single - Funneled

| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | .022 Slot | MBT™ |
|-------------|--|---------------|--------|-----|-----------|-----------|------|
| Arch | | | | | | | |
| -10° | | 0° | 3.6 mm | LL | 067-536 3 | 067-636 3 | MBT™ |
| -10° | | 0° | 3.6 mm | LR | 067-537 4 | 067-637 4 | |
| -25° | | 6° | 3.6 mm | LL | 067-538 3 | 067-638 3 | MBT™ |
| -25° | | 6° | 3.6 mm | LR | 067-539 4 | 067-639 4 | |

Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-536

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-536

Lower Second Non-Convertible Single with Distal Notch and Hook

| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | .022 Slot |
|-------------|--|---------------|--------|-------|-----------|-----------|
| Arch | | | | | | |
| 0° | | 0° | 5.1 mm | UL/LR | 067-844 | 067-944 |
| 0° | | 0° | 5.1 mm | UR/LL | 067-845 | 067-945 |
| 0° | | 6° | 5.1 mm | UL/LR | 067-846 | 067-946 |
| 0° | | 6° | 5.1 mm | UR/LL | 067-847 | 067-947 |

Lower Second Non-Convertible Single with Distal Notch

| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | .022 Slot | MBT™ |
|-------------|--|---------------|--------|-----|-----------|-----------|------|
| Arch | | | | | | | |
| -10° | | 0° | 4.3 mm | LL | 067-816 3 | 067-916 3 | MBT™ |
| -10° | | 0° | 4.3 mm | LR | 067-815 4 | 067-915 4 | |
| -25° | | 6° | 3.6 mm | LL | 067-837 3 | 067-937 3 | MBT™ |
| -25° | | 6° | 3.6 mm | LR | 067-838 4 | 067-938 4 | |
| -25° | | 10° | 3.6 mm | LL | 067-839 3 | 067-939 3 | MBT™ |
| -25° | | 10° | 3.6 mm | LR | 067-840 4 | 067-940 4 | |
| -30° | | 6° | 3.6 mm | LL | 067-835 3 | 067-935 3 | MBT™ |
| -30° | | 6° | 3.6 mm | LR | 067-836 4 | 067-936 4 | |

Adding Prefix "3" is APC™ II Adhesive System – Example: 3067-816

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5067-816

Note

3 Use this prefix to order APC™ II Adhesive Coated Appliance System.

5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive.

APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually.

Measurements are in millimeters unless otherwise noted.

1 2 3 4 Buccal Tube is available with a bondable base.

Victory Series™ Second Molar Mini Tubes

Victory Series™ 2nd Bicuspid Tubes

Buccal Tubes For Partially Erupted Second Molars

Upper Victory Series™ Mini Tube



| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | .022 Slot |
|-------------|--|---------------|--------|-----|-----------|-----------|
| Arch | | | | | | |
| 0° | | 0° | 2.5 mm | UR | 066-4092 | 066-5092 |
| 0° | | 0° | 2.5 mm | UL | 066-4101 | 066-5101 |
| -10° | | 0° | 2.5 mm | UR | 066-4112 | 066-5112 |
| -10° | | 0° | 2.5 mm | UL | 066-4121 | 066-5121 |
| -14° | | 0° | 2.5 mm | UR | 066-4082 | 066-5082 |
| -14° | | 0° | 2.5 mm | UL | 066-4071 | 066-5071 |

MBT™
MBT™

Adding Prefix "3" is APC™ II Adhesive System – Example: 3066-4092

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5066-4092

Lower Victory Series™ Mini Tube

| Slot Torque | | Distal Offset | Length | L/R | .018 Slot | .022 Slot |
|-------------|--|---------------|--------|-----|-----------|-----------|
| Arch | | | | | | |
| 0° | | 0° | 2.5 mm | LR | 066-4024 | 066-5024 |
| 0° | | 0° | 2.5 mm | LL | 066-4013 | 066-5013 |
| -10° | | 0° | 2.5 mm | LR | 066-4044 | 066-5044 |
| -10° | | 0° | 2.5 mm | LL | 066-4033 | 066-5033 |
| -25° | | 0° | 2.5 mm | LR | 066-4064 | 066-5064 |
| -25° | | 0° | 2.5 mm | LL | 066-4053 | 066-5053 |

MBT™
MBT™

Adding Prefix "3" is APC™ II Adhesive System – Example: 3066-4024

Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5066-4024

These tubes are direct bond only, not available welded to bands.



Victory Series 2nd Bicuspid Tubes

| | | | |
|----------------|----------|--|------------------------------|
| 058-801 | 0 | Lower Left 2nd Bicuspid Tube .018 Slot | (-17° torque, 2° angulation) |
| 058-802 | 0 | Lower Right 2nd Bicuspid Tube .018 Slot | (-17° torque, 2° angulation) |
| 058-803 | 0 | Upper Left / Right 2nd Bicuspid Tube .018 Slot | (-7° torque, 0° angulation) |
| 058-901 | 0 | Lower Left 2nd Bicuspid Tube .022 Slot | (-17° torque, 2° angulation) |
| 058-902 | 0 | Lower Right 2nd Bicuspid Tube .022 Slot | (-17° torque, 2° angulation) |
| 058-903 | 0 | Upper Left / Right 2nd Bicuspid Tube .022 Slot | (-7° torque, 0° angulation) |

MBT™
MBT™
MBT™
MBT™
MBT™
MBT™

Offset Base

Adding Prefix "3" is APC™ II Adhesive System – Example: 3058-801


Adding Prefix "5" is APC™ PLUS Adhesive System – Example: 5058-801

Note

| | |
|--|--|
| 3 Use this prefix to order APC™ II Adhesive Coated Appliance System. | 5 Use this prefix to order APC™ PLUS Adhesive Coated Appliance System with color change adhesive. |
| APC™ II Buccal Tubes sold in units of five. APC™ PLUS Buccal Tubes sold in units of four. Non-Coated Buccal Tubes sold individually. | Measurements are in millimeters unless otherwise noted. |
| 1 2 3 4 | Buccal Tube is available with a bondable base. |


Light Wire Buccal Tubes

Upper First Molars




| V-Slot Size | Distal Offset | Light Wire Tube I.D. | Headgear I.D. | Headgear Length | Light Wire Length | L/R | Order Number |
|-------------|---------------|----------------------|---------------|-----------------|-------------------|-----|--------------|
| - | 6° | .036 | .045 | 4.5 mm | 6.5 mm | UL | 046-013 |
| - | 6° | .036 | .045 | 4.5 mm | 6.5 mm | UR | 046-014 |

Upper/Lower First Molars with Vertical Slot



| V-Slot Size | Distal Offset | Light Wire Tube I.D. | Headgear I.D. | Headgear Length | Light Wire Length | L/R | Order Number |
|-------------|---------------|----------------------|---------------|-----------------|-------------------|-------|--------------|
| .021 x .021 | 0° | .036 | | | 6.5 mm | UL/LR | 046-146 |
| .021 x .021 | 0° | .036 | | | 6.5 mm | LL/UR | 046-147 |
| .021 x .021 | 6° | .036 | | | 6.5 mm | UL/LR | 046-176 |
| .021 x .021 | 6° | .036 | | | 6.5 mm | LL/UR | 046-177 |

Upper/Lower First Molars



| Distal Offset | Light Wire Tube I.D. | Headgear I.D. | Headgear Length | Light Wire Length | L/R | Order Number |
|---------------|----------------------|---------------|-----------------|-------------------|-------|--------------|
| 0° | .036 | | | 4.5 mm | UL/LR | 046-114 |
| | | | | | LL/UR | 046-115 |

All Light Wire Buccal Tubes have hooks unless specified.