

Discharge with confidence.

3M™ Prevena™ Plus 125 Therapy Unit is indicated for use with select 3M™ V.A.C.® Dressings* for open wounds.

Modular approach allows you to select Prevena Plus 125 Therapy Unit 14 Days with compatible small- and medium-sized 3M™ V.A.C. Granufoam™ Dressings and 3M™ V.A.C. Whitefoam™ Dressings.



Simplified therapy setting with pre-set continuous negative pressure at -125mmHg for up to 14 days.*



Uses 3M™ SensaT.R.A.C.™ Technology to maintain and monitor pressure at the wound site.



Compatible with select 3M™ V.A.C.® Dressings you know and trust.



May help to facilitate transition from hospital to home.

* V.A.C.® Dressings should be changed by a medical professional every 48–72 hours but no less than three times per week.

Your care team is committed to sending your patient home with everything they need to continue NPWT.

Check to ensure patients discharged with Prevena Plus 125 Therapy Unit have the right amount of supplies depending on their condition and recommendations:

Supplies for 7 days of NPWT**



Prevena Plus 125 Therapy Unit, 14 days**
Includes: one Prevena Plus Therapy unit, one 150ml canister, therapy unit carrying case, lanyard, power supply and cord.



Up to three 3M™ Prevena™ Plus 150ml Canisters



Up to three packs of compatible small or medium-sized 3M™ V.A.C.® Dressing Kit

Supplies for 14 days NPWT**



Prevena Plus 125 Therapy Unit, 14-day**
Includes: one Prevena Plus Therapy unit, one 150ml canister, therapy unit carrying case, lanyard, power supply and cord.



Up to six Prevena Plus 150ml Canisters





Up to six packs of compatible small or medium-sized V.A.C.® Dressing Kit

**7 days of NPWT, the patient should be sent home with 2 to 3 dressing and enough canisters to support the level of exudate (each canister holds 150ml).

**14 days of NPWT, the patient should be sent home with 5 to 6 dressing and enough canisters to support the level of exudate (each canister holds 150ml).

Ordering information

Product name	Part number	
Therapy Units include: one 3M™ Prevena™ Plus 150 ml canister and 3M™ Prevena™ Plus connector, carrying case, lanyard, power supply and cord		
3M™ Prevena™ Plus 125 Therapy Unit – 14 day	PRE4010	
Replacement accessories		
3M™ Prevena™ Plus 150 ml Canister, 5-pack	PRE4095	

Simultaneously, you can continue working on getting the 3M™ ActiV.A.C.® Therapy order approved.

Product name	Part number	
3M™ V.A.C.® Granufoam™ Dressing Kit, Small, 10-Pack, 5-Pack	M8275051/10 M8275051/5	
3M™ V.A.C.® Granufoam™ Dressing Kit, Medium, 10-Pack, 5-Pack	M8275052/10 M8275052/5	
3M™ V.A.C.® Whitefoam™ Dressing Kit Small, 10-Pack, 5-Pack	M8275068/10 M8275068/5	
3M™ V.A.C.® Whitefoam™ Dressing Kit Large, 10-Pack, 5-Pack	M8275067/10 M8275067/5	
3M™ V.A.C.® Whitefoam™ Dressing, Small (foam only), 10-Pack	M6275033/10	
3M™ V.A.C.® Whitefoam™ Dressing, Large (foam only), 10-Pack	M6275034/10	
3M™ VAC® Drape	M6275009/10	

To learn more about 3M™ Prevena™ Plus 125 Therapy Unit and 3M NPWT dressings, contact your local 3M account representative.



KCI Medical Limited
Charnwood Campus
10 Bakewell Road
Loughborough
LE11 5RB
+44 (0)1509 611 611
www.3M.co.uk/Medical

KCI Medical Limited
The Iveagh Building
The Park, Carrickmines
D18 X015
Ireland
+353 (0)1 280 3555
www.3Mireland.ie

*Prevena Plus 125 Therapy Unit is compatible with small- and medium-sized V.A.C.® Dressings. V.A.C.® Dressings should be changed by a medical professional every 48–72 hours, but no less than 3 times per week.

Note: specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals.

© 2024 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorised use prohibited. OMG1073848.