

Delivering patient care wherever the healing is done

3M™ Prevena™ Plus 125 Therapy Unit
provides you with more options

Transitioning patients from hospital to home
with negative pressure wound therapy (NPWT)
is possible.



Simplified therapy setting with
pre-set continuous negative
pressure at -125mmHg for up
to 14 days*



Compatible with select 3M™
V.A.C.® Dressings* you know
and trust



Compact, lightweight and
portable, helping patients to
resume everyday activities



May help to facilitate transition
from hospital to home

For use with a range of select V.A.C.® Dressings
to help deliver the outcomes you need.



3M™ V.A.C.®
Granufoam™
Dressing Kit*



3M™ V.A.C.®
Whitefoam™
Dressing Kit

*Prevena Plus 125 Therapy Unit is compatible with 3M™ V.A.C.® Granufoam™ Dressings (small- to medium-sized)
and 3M™ V.A.C.® Whitefoam™ Dressings.

3M™ V.A.C.® Therapy you know and trust

V.A.C.® Therapy highlights:

- The delivery of V.A.C.® Therapy across care settings, brings benefits to patients and clinicians
- V.A.C.® Therapy is the most studied and used NPWT¹
- Engineered with proprietary 3M™ SensaT.R.A.C.™ Technology that delivers and adjusts pressure at the wound site



Backed by clinical evidence:



21%–34%

lower total costs across wound types at 12 months for patients receiving V.A.C.® Therapy² compared with the standard of care



When V.A.C.® Therapy was initiated early,* it **reduced healing time³** by

1/2 for acute wounds and
1/3 for chronic wounds



More than

75%

of published NPWT clinical evidence is based on V.A.C.® Therapy⁴



More than

10 million

wounds treated worldwide using V.A.C.® Therapy⁵

*Defined as treatment within the first 7 days for acute wounds and 30 days for chronic wounds from the first wound treatment date.

Navigating through a patient's care continuum

Discharge patients and confidently continue their wound healing at home with V.A.C.® Therapy

3M™ Prevena™ Plus 125 Therapy Unit, an easy-to-use, single-use disposable therapy unit is now compatible with select* V.A.C.® Dressings. Transition patients from hospital to home for up to 14 days of V.A.C.® Therapy with the V.A.C.® Dressings you know and trust.




*Prevena Plus 125 Therapy Unit is compatible with 3M™ V.A.C.® Granufoam™ Dressings (small- to medium-sized) and 3M™ V.A.C.® Whitefoam™ Dressings.

Transitioning patients with more options

If you're ready to discharge and your community 3M™ ActiV.A.C.™ Therapy Unit isn't ready, you now have the option to discharge with Prevena Plus 125 Therapy Unit for up to 14 days of V.A.C.® Therapy.

Discharge with the Prevena Plus 125 Therapy Unit until the ActiV.A.C. Therapy Unit becomes available or the wound is healed.










Prevena Plus 125 Therapy Unit is off-the-shelf available, disposable, and compatible with select V.A.C.® Therapy Dressings.



Check to ensure patients discharged with a 3M™ Prevena™ Plus 125 Therapy Unit have the right amount of supplies for their condition.

Prevena Plus 125 Therapy Unit can be continuously applied for up to 14 days. 3M™ V.A.C.® Dressings should be changed every 48–72 hours. For 14 days of Prevena Plus 125 Therapy, the patient will require five to six dressings and enough canisters to support the level of exudate (each canister holds 150 ml).

Ordering information

Product Name	Part number		Product Name	Part number	
Therapy Units include: one 3M™ Prevena™ Plus 150 ml canister and 3M™ Prevena™ Plus carrying case, lanyard, power supply and cord			3M™ V.A.C.® Granufoam™ Dressing Kit* Small, 10-Pack, 5-Pack	M8275051/10 M8275051/5	
3M™ Prevena™ Plus 125 Therapy Unit – 14 day	PRE4010		3M™ V.A.C.® Granufoam™ Dressing Kit* Medium, 10-Pack, 5-Pack	M8275052/10 M8275052/5	
3M™ Prevena™ Plus 150 ml Canister 5-pack	PRE4095		3M™ V.A.C.® Whitefoam™ Dressing Kit* Small, 10-Pack, 5-Pack	M8275068/10 M8275068/5	
			3M™ V.A.C.® Whitefoam™ Dressing Kit* Large, 10-Pack, 5-Pack	M8275067/10 M8275067/5	
			3M™ V.A.C.® Whitefoam™ Dressing Small (foam only), 10-Pack	M6275033/10	
			3M™ V.A.C.® Whitefoam™ Dressing Large (foam only), 10-Pack	M6275034/10	
			3M™ VAC® Drape	M6275009/10	

The 3M™ Prevena™ Plus Therapy Unit is also compatible with 3M™ Prevena™ Dressings and 3M™ Prevena Restor™ Dressings for comprehensive protection of the surgical site (incision and surrounding soft tissue).

Visit [3M.co.uk/Prevena](https://www.3M.co.uk/Prevena) for more information on Prevena™ Therapy for closed surgical incisions.



To learn more about 3M NPWT Solutions, contact your local 3M account representative.

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for these products and therapies. Please consult a clinician and product instructions for use prior to application. This material is intended for healthcare professionals.

References:

- 3M Data on File. Cumulative NPWT Wounds. 2021.
- Law AL, Krebs B, Karnik B, Griffin L. Comparison of healthcare costs associated with patients receiving traditional negative pressure wound therapies in the post acute setting. *Cureus*. 2020 Nov 30;12(11):e11790.
- Miller-Mikolajczyk C, Achi J, James R. Real world use: comparing early versus late initiation of negative pressure wound therapy on wound surface area reduction in patients at wound care clinics. Poster presented at The Wound Ostomy and Continence Nurses Society Annual Conference, June 22-26, 2013. Seattle, Washington.
- KCI. Percentage of V.A.C.® Therapy Articles vs. Comp Articles. May 7, 2020.
- 3M. Cumulative NPWT Wounds. 2018.

© 2024 3M. All rights reserved. 3M and the other marks shown are marks and/or registered marks. Unauthorized use prohibited.



KCI Medical Limited
Charnwood Campus
10 Bakewell Road
Loughborough
LE11 5RB
+44 (0)1509 611 611
www.3M.co.uk/Medical

KCI Medical Limited
The Iveagh Building
The Park, Carrickmines
D18 X015
Ireland
+353 (0)1 280 3555
www.3Mireland.ie