

3M™ Steri-Drape™ One-Piece Hip Drape with Side Pouches – 9055



Large Side Pockets

Secure Adhesion

This 3M Steri-Drape Hip Drape contains a reliable adhesive strip on the inside of the drape which provides secure adhesion to the skin throughout the procedure.

- Large side pockets with flaps for fluid collection and control
- Extra long and wide drape assures complete coverage, even during leg manipulation
- One-piece draping
- The fabric absorbs spills and splatters, preventing strikethrough and making multiple layers unnecessary



Lead
The Way

3M

3M™ Steri-Drape™ One-Piece Hip Drape with Side Pouches – 9055

3M Fabric

- Meets Flammability Class 1 (CFR Title 16 Part 1610)
- Meets AAMI Standard PB70:2012 *Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities, Selection of surgical gowns and drapes in health care facilities*
- Provides AAMI level 4 barrier performance throughout entire drape
- Low linting (INDA IST 160 1-83)
- Meets AORN 2014 Recommended Practices for materials, selection, and use of surgical drapes that represent an optimal level of practice



Ordering Information

3M™ Steri-Drape™ One-Piece Hip Drape with Side Pouches – 9055

Cat No.	Product Description	Packaging	Minimum Order	Pricing Unit
9055	Overall: 118" x 138" (300cm x 353cm) Aperture: 16" x 38" (43cm x 97cm)	1 Each 8 Each/Box 2 Box/Case	Case	Box

For more information about 3M™ Surgical Drapes or other 3M Health Care products:

- » Contact your local 3M representative
- » Call the 3M Health Care Customer Helpline at 800-228-3957
- » Visit www.3M.com/infectionprevention



Infection Prevention Division 3M Health Care

2510 Conway Avenue
St. Paul, MN 55144-1000
U.S.A.
1-800-228-3957
www.3M.com/infectionprevention

3M and Steri-Drape are trademarks of 3M.
©3M 2014. All rights reserved.
70-2009-7034-4