

## Express™ and Garant™ Dispenser Reliable. Easy to Use. Accurate.



Express

### Consultant's comments

- "Putty is great - easy to mix - love it"
- "Excellent system. Low viscosity stacks well, yet flows into the sulcus with ease, resulting in great impressions"
- "Great Wettability"
- "Good system for patients with a strong gag reflex . putty stays put"

## Express™ VPS Impression Material and Garant™ Dispenser – The System of Choice!

Express provides you with the reliability proven by millions of impressions made worldwide in over 20 years. Due to its various components, Express VPS impression material allows for maximum technique flexibility. You have the freedom to use either the two-step or simultaneous putty-wash technique. Due to its flow properties and hydrophilicity Express wash materials provide an excellent reproduction of details.

The new design of mixing tips allows faster and easier mixing with precise access thanks to the fine intraoral tip. Used in the putty-wash technique, the orange Express STD Putty mixes easily without streaking, while the thixotropic wash materials flow freely yet stays on your preparation without slumping.

The Garant™ Dispenser promises no more plugged cartridges because of reusable cap and an easier dispensing system. The additional advantages of this system are:

- Dual-port cartridge minimizes risk of cross contamination
- Garant™ dispenser works with all 3M ESPE syringeable impression materials
- Colour-coded cartridges and mixing tips for easy use.
- The Garant dispenser is completely steam autoclavable.
- Smaller mixing tip reduces waste up to 35%.

# New: Easier Application and Strain-Free Dispensing of Material with the Garant™ Dispenser.











**Express, the most popular VPS meets ADA/ISO Specifications and was awarded with 4-Stars from the Dental Advisor.**



*“3M ESPE Express is a highly recommended addition cured silicone designed for taking crown and bridge, denture and bite registration impressions. It received a 94% rating”*

Dental Advisor, Vol 16, No 6, July / August 1999 (rated + + + + 1/2)

## Ordering Information

Item Number	Product Information	
 <b>7300 H</b>	<b>Introductory Kit</b> · 2 x Putty Jar à 305ml (550g) · 2 x 50ml LB Regular Set cartridges · 1 x Garant Dispenser · 10 x Garant Mixing Tips (yellow)	· 10 x Garant Intra-oral Tips (yellow) · 1 x VPS Tray Adhesive · 2 x Putty spoons
 <b>7301</b> Light Body Fast Set	<b>Light Body, Fast Set Refill (Blue)</b> · 2 x 50ml cartridges · 10 x Garant Mixing Tips	
 <b>7302/7322</b> Light Body/Med Regular Set	<b>Light Body, Regular Set Refill (Green)</b> · 2 x 50ml cartridges · 10 x Garant Mixing Tips	<b>Medium Body, Regular Set Refill (Purple)</b> · 2 x 50ml cartridges · 10 x Garant Mixing Tips
 <b>7312</b>	<b>Putty Refill</b> · 2 x Putty Jar regular set, 305ml/550g each	
 <b>6160J</b>	<b>Bite Registration Putty Refill</b> · 2 x Putty Jar fast set, 305ml/550g each	
 <b>77580</b>	<b>Garant Dispenser</b> · 1 x Garant Dispenser	
 <b>71452</b> Mixing Tips	<b>Mixing Tips</b> · 50 x mix tips (yellow)	 <b>71462</b> <b>Intraoral Tips</b> 50 x Intraoral tips (yellow)
 <b>7307</b>	<b>VPS Tray Adhesive</b> · 1 x 24ml Tray Adhesive	



**Dual-port cartridge minimizes risk of cross contamination and plugged cartridges.**

# 3M ESPE

## Dental Products

3M Center  
Building 275-2SE-03  
St. Paul, MN 55144-1000  
USA

3M, ESPE, Express are trademarks of 3M or 3M ESPE AG. Used under license in Canada.

Printed in  
© 3M 2004  
00-0000-0000-0