

		5,4 YELLOW Bicuspid w/hk	3 RED Cuspid w/hk	Gingival				3 BLUE Cuspid w/hk	4,5 YELLOW Bicuspid w/hk		
	6 RED 1st Molar	5,4 YELLOW Univ BI	3 RED Cuspid	2 RED Lateral	1 RED Central	1 BLUE Central	2 BLUE Lateral	3 BLUE Cuspid	4,5 YELLOW Univ BI	6 RED 1st Molar	
<i>Patient's Right</i>				Occlusal				<i>Patient's Left</i>			
	5 BLACK 2nd Bicuspid	4 BLACK 1st Bicuspid	3 BLACK Cuspid	1,2 NO INK Anterior	1,2 NO INK Anterior	1,2 NO INK Anterior	1,2 NO INK Anterior	3 CREME Cuspid	4 CREME 1st Bicuspid	5 CREME 2nd Bicuspid	
6 BLACK 1st Molar	5 BLACK 2nd Bicuspid w/hk	4 BLACK 1st Bicuspid w/hk	3 BLACK Cuspid w/hk	Gingival				3 CREME Cuspid w/hk	4 CREME 1st Bicuspid w/hk	5 CREME 2nd Bicuspid w/hk	6 BLACK 1st Molar

Identification Information

<p>Labial Root Torque</p> <p>3</p>	<p>Regular Base</p> <p>4,5 4,5</p>	<p>4</p>
<p>180° Inversion Provides Lingual Root Torque</p> <p>3</p>	<p>Offset Base</p> <p>4,5 4,5</p>	<p>5</p>

3M
3M Unitek
Orthodontic Products
 2724 South Peck Road
 Monrovia, CA 91016 USA

EC REP **CE** 0086
3M Deutschland GmbH
Orthodontic Products
 Carl-Schurz-Str. 1
 D-41453 Neuss, Germany
 +49(0)8191/ 9474-0

6	5	4	3	2	1	1	2	3	4	5	6
6	5	4	3	2	1	1	2	3	4	5	6