

3M™ Steri-Drape™ Absorbent Prevention Fabric, Venous/Arterial Access Drape – 1192

The drape that fits. Accept nothing less.

The 3M™ Steri-Drape™ Venous/Arterial Access Drape – 1192 is designed to fit the basic needs for venous access procedures in the Cath Lab, Cardiovascular Surgery, and other specialty areas. It creates a sterile field for a variety of procedures:

- Pacemaker insertion and repairs
- Brachial and radial artery access
- Brachial and radial artery harvesting
- Central and PICC line insertions

The 1192 drape can be used as a stand alone drape, or in conjunction with other 3M surgical drapes.

The 3M™ Absorbent Prevention Fabric helps reduce the risk of exposure to blood and body fluids for both patients and healthcare workers by absorbing spills and splatters while preventing strike-through of fluids and reducing the risk of contamination of the operative site. No fluids get through, so double layering is unnecessary, saving time and money.


- **Large diamond shaped aperture makes landmarks easy to identify**
- **Reliable 3M adhesive around the aperture ensures that drape migration does NOT occur during the procedure**
- **Adhesive tabs enable easy attachment to CV drapes during peripheral artery harvesting**
- **Sized to provide body coverage for pacemaker and Central/PICC line insertions (neck to ankle)**
- **Low linting (INDA IST 160 1–83)**
- **3M Absorbent Prevention Fabric meets Flammability Class 1 (CFR Title 16 Part 1610)**
- **3M Absorbent Fabric meets AAMI Standard PB70: 2003 *Liquid barrier performance and classification of protective apparel and drapes intended for use in health care facilities*, and AAMI Technical Information Report TIR 11:2005, *Selection of surgical gowns and drapes in health care facilities***



1192 drape with 6661EZ 3M™ Ioban™ 2 Antimicrobial Incise Drape



1192 Brachial/Radial Access




1192 Central/PICC Line

3M™ Steri-Drape™ Absorbent Prevention Fabric, Venous/Arterial Access Drape – 1192

Ordering Information

Catalog Number Product Description Size	Packaging	UPN	Minimum Order	Pricing Unit
1192 3M™ Steri-Drape™ Absorbent Prevention Fabric, Venous/Arterial Access Drape	Each 12 Each/Box 4 Box/Case	0070738756797 3070738756797 5070738756797	1 Case	Each
1192NS Custom Packs Only	Each 68 Each/Case	0070387567989 50707387567984	1 Case	Each

Overall: 60" x 67" (152 cm x 170 cm)
Aperture: 3" x 8" (7.6 cm x 20 cm)



For more information

For more information about 3M surgical drapes or other 3M Health Care products:

- Contact your local 3M representative
- Call the 3M Health Care Customer Helpline: **1-800-228-3957**
- Visit the 3M Health Care products Web site: www.3M.com/infectionprevention

3M Infection Prevention Solutions. A vital step in infection prevention.



Health Care

3M Center, Building 275-4W-02
St. Paul, MN 55144-1000
U.S.A.
1 800 228-3957

www.3M.com/healthcare

3M Canada

Post Office Box 5757
London, Ontario N6A4T1
Canada
1 800-563-2921

Outside of USA, please contact your 3M Representative



Printed in U.S.A.
© 3M 2006
All Rights Reserved.
70-2009-7403-1 (106.4) DPI