

Express™ 2 VPS Impression Materials Portfolio

Product	Dispensing System	Viscosity Low  High	Setting Version	Working Time (23° C) 	Intra-oral Setting Time (37° C) 
TRAY MATERIALS					
Express™ 2 Penta™ Putty Soft 		Putty	Regular Set	1:30	3:30
Express™ 2 Penta™ Putty		Putty	Regular Set	1:30	3:00
Express™ 2 Penta™ H		Heavy Body	Regular Set	2:00	3:30
Express™ 2 Penta™ H Universal Quick		Heavy Body	Fast Set	1:30	2:30
Express™ 2 Penta™ H Quick		Heavy Body	Fast Set	1:15	2:30
WASH MATERIALS					
Express™ 2 Ultra-Light Body Quick			Fast Set	1:30	2:30
Express™ 2 Light Body <i>Flow</i>			Regular Set	2:00	3:30
Express™ 2 Light Body <i>Flow</i> Quick			Fast Set	1:30	2:30
Express™ 2 Light Body <i>Standard</i>			Regular Set	2:00	3:30
Express™ 2 Light Body <i>Standard</i> Quick			Fast Set	1:30	2:30
Express™ 2 Regular Body			Regular Set	2:00	3:30
Express™ 2 Regular Body Quick			Fast Set	1:30	2:30

Express™ 2 VPS Impression Materials per Technique

Impression Technique	Tray Material	Recommended Wash Materials	Wash Material
Express™ 2 Penta™ Putty Soft ▶ first hydrophilic auto-mix putty			
1-step technique 		Express™ 2 Light Body <i>Flow</i> ▶ offers great flow	
		Express™ 2 Regular Body ▶ stays very well at tooth	
Express™ 2 Penta™ Putty ▶ putty consistency			
2-step technique		Express™ 2 Ultra-Light Body Quick ▶ generates very thin layers	
		Express™ 2 Light Body <i>Flow</i> Quick ▶ offers great flow	
Express™ 2 Penta™ H ▶ hydrophilic heavy body			
1-step technique		Express™ 2 Light Body <i>Flow</i> ▶ offers great flow	
		Express™ 2 Light Body <i>Standard</i> ▶ stays very well at tooth	
Express™ 2 Penta™ H Universal Quick ▶ fast setting heavy body with excellent carving properties			
2-step technique		Express™ 2 Light Body <i>Flow</i> Quick ▶ offers great flow	
1-step technique		Express™ 2 Light Body <i>Standard</i> Quick ▶ stays very well at tooth	
Express™ 2 Penta™ H Quick ▶ fast setting hydrophilic heavy body			
1-step technique		Express™ 2 Light Body <i>Flow</i> Quick ▶ offers great flow	
		Express™ 2 Light Body <i>Standard</i> Quick ▶ stays very well at tooth	