

# 3M™ Scotchcast™ Conformable Roll Splint

## Commonly Asked Questions

**Question:**

What is the advantage of using fiberglass splints vs. plaster?

**Answer:**

Fiberglass is lightweight and has a lower profile; thus, patient is able to lead a more normal lifestyle. Fiberglass is stronger, there is less splint breakdown, and it is less messy for the clinician.

**Question:**

How does the Scotchcast splint set time compare to what I'm using now?

**Answer:**

You control the working time based on your application and needs.

**Question:**

What education and training is available for your splints?

**Answer:**

3M offers workshops on both splinting and casting.

**For More Information**

Contact your 3M Health Care Sales Representative, or call the 3M Health Care Customer Helpline at **1-800-228-3957**. These products can be ordered from your local distributor. Outside the United States, contact the local 3M subsidiary.



**Health Care**

3M Center, Building 275-4W-02  
St. Paul, MN 55144-1000  
USA  
1 800 228-3957  
[www.3M.com/healthcare](http://www.3M.com/healthcare)