



# Black is Back

**Our popular 3M™ Littmann® Black Edition Series returns with more choices.**

A Littmann Black Edition Stethoscope has always been a symbol of elegance and professionalism. That's why these models have been so popular, both as a personal choice and as gifts. It's also why we're now giving you more Black Edition models to choose from. In addition to Littmann Master Cardiology™ and Master Classic II Black Edition Stethoscopes, we now offer Black Editions of the Littmann Cardiology III and Classic II S.E. models. Pick the model that suits you or your gift recipient best. You can be sure anyone will look good in black.


To order, see your 3M Littmann authorized dealer.

Or, to learn more, visit [www.littmann.com](http://www.littmann.com).





**3M**

# 3M™ Littmann® Black Edition Stethoscopes

A trusted leader in innovative auscultation technology, Littmann stethoscopes are renowned worldwide for unsurpassed quality. The Littmann brand name is your assurance of precision, acoustical superiority, innovative design and exceptional performance.



A Classic II S.E. Stethoscope 2218BE    
 C Master Cardiology™ Stethoscope 2161  
B Cardiology III Stethoscope 3131BE    
 D Master Classic II Stethoscope 2141

| Catalog Number | Description   | Length           | Tube Color |   |
|----------------|---|------------------|------------|---|
| 2161           | <b>Littmann Master Cardiology™ Stethoscope</b> <ul style="list-style-type: none"> <li>Offers a handcrafted stainless steel chestpiece with the portable convenience of a single-sided stethoscope.</li> <li>The tunable diaphragm conveniently alternates between low- and high-frequency sounds with a simple pressure change on the chestpiece.</li> <li>Recommended for cardiology and high performance adult and pediatric use.</li> <li>Dual-lumen tubing.</li> </ul>  | 27 in. (69 cm)   | Black      |    |
| 3131BE         | <b>Littmann Cardiology III Stethoscope</b> <ul style="list-style-type: none"> <li>Offers two tunable diaphragms which allows the user to conveniently alternate between low- and high-frequency sounds without turning over the chestpiece.</li> <li>The large side can be used for adult patients, while the small side is especially useful for pediatric or thin patients, around bandages, and for carotid assessment.</li> <li>Recommended for cardiology and high performance adult and pediatric use.</li> <li>Dual-lumen tubing.</li> </ul> | 27 in. (69 cm)   | Black      |  |
| 2141           | <b>Littmann Master Classic II Stethoscope</b> <ul style="list-style-type: none"> <li>Traditional bell combined with a tunable diaphragm that allows the user to conveniently alternate between low- and high-frequency sounds without turning over the chestpiece.</li> <li>Single-lumen tubing.</li> </ul>   | 27 in. (69 cm)   | Black      |  |
| 2218BE         | <b>Littmann Classic II S.E. Stethoscope</b> <ul style="list-style-type: none"> <li>Features a two-sided chestpiece that offers high acoustic sensitivity for superior performance.</li> <li>The tunable diaphragm on one side augments the traditional bell function of the opposite side.</li> <li>Single-lumen tubing.</li> </ul>   | 28 in. (71,1 cm) | Black      |  |



## 3M Health Care Littmann Brand Products

3M Center, Building 275-4W-02  
St. Paul, MN, 55144-1000  
USA  
800 228 3957  
www.littmann.com

3M, Master Cardiology, and  
Littmann are trademarks of 3M.  
Please recycle. Printed in U.S.A.  
© 3M 2007, 2009. All rights reserved.  
70-2009-8642-3