









































Transcend™ 6000 Ceramic Brackets Bracket Identification System Roth Rx 011-458-4

<i>Patient's Right</i>					Gingival										<i>Patient's Left</i>				
6001-617, -717  -7° T, 0° A 1st & 2nd NAR BI W/HK	6001-615, -715  -7° T, 0° A 1st & 2nd BI W/HK	6001-606, -706  0° T, +11° A Cuspid W/HK											6001-605, -705  0° T, +11° A Cuspid W/HK	6001-615, -715  -7° T, 0° A 1st & 2nd BI W/HK	6001-617, -717  -7° T, 0° A 1st & 2nd NAR BI W/HK				
6001-618, -718  -7° T, 0° A 1st & 2nd NAR BI	6001-616, -716  -7° T, 0° A 1st & 2nd Bicuspid	6001-608, -708  0° T, +11° A Cuspid	6001-604, -704  +8° T, +9° A Lateral	6001-602, -702  +12° T, +5° A Central	6001-601, -701  +12° T, +5° A Central	6001-603, -703  +8° T, +9° A Lateral	6001-607, -707  0° T, +11° A Cuspid	6001-616, -716  -7° T, 0° A 1st & 2nd Bicuspid	6001-618, -718  -7° T, 0° A 1st & 2nd NAR BI	Occlusal									
6001-646, -746  -22° T, 0° A 2nd NAR BI	6001-644, -744  -22° T, 0° A 2nd BI	6001-638, -738  -17° T, 0° A 1st NAR BI	6001-636, -736  -17° T, 0° A 1st BI	6001-628, -728  -11° T, +7° A Cuspid	6001-621, -721  -1° T, 0° A Anterior	6001-621, -721  -1° T, 0° A Anterior	6001-621, -721  -1° T, 0° A Anterior	6001-621, -721  -1° T, 0° A Anterior	6001-627, -727  -11° T, +7° A Cuspid	6001-636, -736  -17° T, 0° A 1st BI	6001-638, -738  -17° T, 0° A 1st NAR BI	6001-644, -744  -22° T, 0° A 2nd BI	6001-646, -746 -22° T, 0° A 2nd NAR BI						
6001-645, -745  -22° T, 0° A 2nd NAR BI W/HK	6001-643, -743  -22° T, 0° A 2nd BI W/HK	6001-637, -737  -17° T, 0° A 1st NAR BI W/HK	6001-635, -735  -17° T, 0° A 1st BI W/HK	6001-626, -726  -11° T, +7° A Cuspid W/HK											6001-625, -725  -11° T, +7° A Cuspid W/HK	6001-635, -735  -17° T, 0° A 1st BI W/HK	6001-637, -737  -17° T, 0° A 1st NAR BI W/HK	6001-643, -743  -22° T, 0° A 2nd BI W/HK	6001-645, -745  -22° T, 0° A 2nd NAR BI W/HK
					Gingival														

How to use the Transcend™ 6000 Identification System

011-458-4



Made in U.S.A.



**3M Unitek
Orthodontic Products**

2724 South Peck Road
Monrovia, CA 91016 USA



European Representative:

3M Unitek GmbH

ESPE Platz
82229 Seefeld, Germany
0049/(0)8152-397-0

Caution: U.S. Federal Law restricts this device to sale by or on the order of an orthodontic professional.

© 2004 3M Company. All rights reserved.

REF 011-458-4 0405