



































Clarity™ Metal-Reinforced Ceramic Brackets Bracket Identification System

MBT™ Rx

Gingival									
6400-610 6400-710  YELLOW -7° T 0° A 1st & 2nd Bi w/hook	6400-808 6400-908  GREEN 0° T +8° A Cuspid w/hook	6400-682 6400-782  GREEN -7° T +8° A Cuspid w/hook					6400-683 6400-783  GREEN -7° T +8° A Cuspid w/hook	6400-807 6400-907  GREEN 0° T +8° A Cuspid w/hook	6400-609 6400-709  YELLOW -7° T 0° A 1st & 2nd Bi w/hook
6400-611 6400-711  YELLOW -7° T 0° A 1st & 2nd Bicuspid	6400-818 6400-918  GREEN 0° T +8° A Cuspid	6400-806 6400-906  GREEN -7° T +8° A Cuspid	6400-804 6400-904  RED +10° T +8° A Lateral	6400-802 6400-902  BLUE +17° T +4° A Central	6400-801 6400-901  BLUE +17° T +4° A Central	6400-803 6400-903  RED +10° T +8° A Lateral	6400-805 6400-905  GREEN -7° T +8° A Cuspid	6400-817 6400-917  GREEN 0° T +8° A Cuspid	6400-611 6400-711  YELLOW -7° T 0° A 1st & 2nd Bicuspid
Patient's Right					Occlusal		Patient's Left		
6400-690 6400-790  LT. GREY -17° T +2° A 2nd Bicuspid	6400-686 6400-786  LT. BLUE -12° T +2° A 1st Bicuspid	6400-698 6400-798  PURPLE 0° T +3° A Cuspid	6400-696 6400-796  PURPLE -6° T +3° A Cuspid	6400-810 6400-910  PINK -6° T 0° A Anterior	6400-810 6400-910  PINK -6° T 0° A Anterior	6400-695 6400-795  PURPLE -6° T +3° A Cuspid	6400-697 6400-797  PURPLE 0° T +3° A Cuspid	6400-687 6400-787  LT. BLUE -12° T +2° A 1st Bicuspid	6400-691 6400-791  LT. GREY -17° T +2° A 2nd Bicuspid
6400-692 6400-792  LT. GREY -17° T +2° A 2nd Bicuspid w/hook	6400-688 6400-788  LT. BLUE -12° T +2° A 1st Bicuspid w/hook	6400-700 6400-800  PURPLE 0° T +3° A Cuspid w/hook	6400-684 6400-784  PURPLE -6° T +3° A Cuspid w/hook	Gingival		6400-685 6400-785  PURPLE -6° T +3° A Cuspid w/hook	6400-699 6400-799  PURPLE 0° T +3° A Cuspid w/hook	6400-689 6400-789  LT. BLUE -12° T +2° A 1st Bicuspid w/hook	6400-693 6400-793  LT. GREY -17° T +2° A 2nd Bicuspid w/hook

Identification Information - Refer to Instruction Manual for details on Debonding, Cautions and Warnings.






Step 1



Step 2



Step 3



Made in U.S.A.



**3M Unitek
Orthodontic Products**
2724 South Peck Road
Monrovia, CA 91016 USA

CE
0086 93/42/EEC

European Representative:
3M Unitek GmbH
ESPE Platz
82229 Seefeld, Germany
0049/(0)8152-397-0

5	4	3	2	1	1	2	3	4	5
5	4	3	2	1	1	2	3	4	5

Caution: U.S. Federal Law restricts this device to sale by or on the order of an orthodontic professional.

© 2004 3M Company. All rights reserved.

REF 011-558-9 0405