




























Clarity™ Metal-Reinforced Ceramic Brackets Bracket Identification System **High Torque**

Patient's Right								Gingival				Patient's Left							
6400-674 6400-774  5,4 YELLOW 0° T 0° A 1st & 2nd Bicuspid w/hook		6400-656 6400-756  3 LT. GREEN +7° T +5° A Cuspid w/hook						6400-655 6400-755  3 LT. GREEN +7° T +5° A Cuspid w/hook		6400-675 6400-775  4,5 YELLOW 0° T 0° A 1st & 2nd Bicuspid w/hook									
6400-823 6400-923  5,4 YELLOW 0° T 0° A 1st & 2nd Bicuspid		6400-670 6400-770  3 LT. GREEN +7° T +5° A Cuspid		6400-654 6400-754  2 RED +14° T +8° A Lateral		6400-652 6400-752  1 LT. BLUE +22° T +5° A Central		6400-651 6400-751  1 LT. BLUE +22° T +5° A Central		6400-653 6400-753  2 RED +14° T +8° A Lateral		6400-671 6400-771  3 LT. GREEN +7° T +5° A Cuspid		6400-823 6400-923  4,5 YELLOW 0° T 0° A 1st & 2nd Bicuspid					
Patient's Right								Occlusal				Patient's Left							
6400-659 6400-759  5 BROWN -11° T 0° A 2nd Bicuspid		6400-660 6400-760  4 YELLOW -7° T 0° A 1st Bicuspid		6400-672 6400-772  3 PURPLE +7° T +5° A Cuspid		6400-824 6400-924  2,1 PINK 0° T 0° A Anterior		6400-824 6400-924  1,2 PINK 0° T 0° A Anterior		6400-673 6400-773  3 PURPLE +7° T +5° A Cuspid		6400-660 6400-760  4 YELLOW -7° T 0° A 1st Bicuspid		6400-659 6400-759  5 BROWN -11° T 0° A 2nd Bicuspid					
6400-676 6400-776  5 BROWN -11° T 0° A 2nd Bicuspid w/hook		6400-662 6400-762  4 YELLOW -7° T 0° A 1st Bicuspid w/hook		6400-658 6400-758  3 PURPLE +7° T +5° A Cuspid w/hook						6400-657 6400-757  3 PURPLE +7° T +5° A Cuspid w/hook		6400-661 6400-761  4 YELLOW -7° T 0° A 1st Bicuspid w/hook		6400-677 6400-777  5 BROWN -11° T 0° A 2nd Bicuspid w/hook					
Patient's Right								Gingival				Patient's Left							

Identification Information - Refer to Instruction Manual for details on Debonding, Cautions and Warnings.






Step 1



Step 2



Step 3



Made in U.S.A.



**3M Unitek
Orthodontic Products**
2724 South Peck Road
Monrovia, CA 91016 USA

CE
0086 93/42/EEC

European Representative:
3M Unitek GmbH
ESPE Platz
82229 Seefeld, Germany
0049/(0)8152-397-0

5	4	3	2	1	1	2	3	4	5
5	4	3	2	1	1	2	3	4	5

Caution: U.S. Federal Law restricts this device to sale by or on the order of an orthodontic professional.

© 2004 3M Company. All rights reserved.

REF 011-575-3 0405