



























Clarity™ Metal-Reinforced Ceramic Brackets Bracket Identification System

Roth Rx

6400-610 6400-710  5,4 YELLOW -7° T 0° A 1st & 2nd Bicuspid w/hook		6400-606 6400-706  3 LT. GREEN 0° T +11° A Cuspid w/hook		Gingival				6400-605 6400-705  3 LT. GREEN 0° T +11° A Cuspid w/hook		6400-609 6400-709  4,5 YELLOW -7° T 0° A 1st & 2nd Bicuspid w/hook					
6400-611 6400-711  5,4 YELLOW -7° T 0° A 1st & 2nd Bicuspid		6400-608 6400-708  3 LT. GREEN 0° T +11° A Cuspid		6400-604 6400-704  2 RED +8° T +9° A Lateral		6400-602 6400-702  1 LT. BLUE +12° T +5° A Central		6400-601 6400-701  1 LT. BLUE +12° T +5° A Central		6400-603 6400-703  2 RED +8° T +9° A Lateral		6400-607 6400-707  3 LT. GREEN 0° T +11° A Cuspid		6400-611 6400-711  4,5 YELLOW -7° T 0° A 1st & 2nd Bicuspid	
<i>Patient's Right</i>				Occlusal				<i>Patient's Left</i>							
6400-644 6400-744  5 BLACK -22° T 0° A 2nd Bicuspid		6400-636 6400-736  4 ORANGE -17° T 0° A 1st Bicuspid		6400-628 6400-728  3 PURPLE -11° T +7° A Cuspid		6400-621 6400-721  2,1 PINK -1° T 0° A Anterior		6400-621 6400-721  1,2 PINK -1° T 0° A Anterior		6400-627 6400-727  3 PURPLE -11° T +7° A Cuspid		6400-636 6400-736  4 ORANGE -17° T 0° A 1st Bicuspid		6400-644 6400-744  5 BLACK -22° T 0° A 2nd Bicuspid	
6400-638 6400-738  5 BLACK -22° T 0° A 2nd Bicuspid w/hook		6400-630 6400-730  4 ORANGE -17° T 0° A 1st Bicuspid w/hook		6400-626 6400-726  3 PURPLE -11° T +7° A Cuspid w/hook		Gingival				6400-625 6400-725  3 PURPLE -11° T +7° A Cuspid w/hook		6400-629 6400-729  4 ORANGE -17° T 0° A 1st Bicuspid w/hook		6400-637 6400-737  5 BLACK -22° T 0° A 2nd Bicuspid w/hook	

Identification Information - Refer to Instruction Manual for details on Debonding, Cautions and Warnings.







Step 1



Step 2



Step 3



Made in U.S.A.



**3M Unitek
Orthodontic Products**
2724 South Peck Road
Monrovia, CA 91016 USA

CE
0086 93/42/EEC

European Representative:
3M Unitek GmbH
ESPE Platz
82229 Seefeld, Germany
0049/(0)8152-397-0

5	4	3	2	1	1	2	3	4	5
5	4	3	2	1	1	2	3	4	5

Caution: U.S. Federal Law restricts this device to sale by or on the order of an orthodontic professional.

© 2004 3M Company. All rights reserved.

REF 011-545-4 0405