























**Transcend™ 6000 Ceramic Brackets** Bracket Identification System **Standard Edgewise**


<i>Patient's Right</i>				<i>Patient's Left</i>					
<b>Gingival</b>									
6001-805 6001-905	6001-805 6001-905	6001-805 6001-905	6001-801 6001-901	6001-801 6001-901	6001-801 6001-901	6001-801 6001-901	6001-805 6001-905	6001-805 6001-905	6001-805 6001-905
									
3,4,5	3,4,5	3,4,5	1,2	1,2	1,2	1,2	3,4,5	3,4,5	3,4,5
0° T 0° A 2nd Bicuspid	0° T 0° A 1st Bicuspid	0° T 0° A Cuspid	0° T 0° A Lateral	0° T 0° A Central	0° T 0° A Central	0° T 0° A Lateral	0° T 0° A Cuspid	0° T 0° A 1st Bicuspid	0° T 0° A 2nd Bicuspid
<b>Occlusal</b>									
6001-805 6001-905	6001-805 6001-905	6001-805 6001-905	6001-821 6001-921	6001-821 6001-921	6001-821 6001-921	6001-821 6001-921	6001-805 6001-905	6001-805 6001-905	6001-805 6001-905
									
3,4,5	3,4,5	3,4,5	1,2	1,2	1,2	1,2	3,4,5	3,4,5	3,4,5
0° T 0° A 2nd Bicuspid	0° T 0° A 1st Bicuspid	0° T 0° A Cuspid	0° T 0° A Anterior	0° T 0° A Anterior	0° T 0° A Anterior	0° T 0° A Anterior	0° T 0° A Cuspid	0° T 0° A 1st Bicuspid	0° T 0° A 2nd Bicuspid
<b>Gingival</b>									

## How to use the Transcend™ 6000 Identification System

011-461-2



Align Vertical Centerline of Indicator  
with Long Axis of the Tooth



Made in U.S.A.

**3M**

**3M Unitek  
Orthodontic Products**

2724 South Peck Road  
Monrovia, CA 91016 USA

**CE**

**0086 93/42/EEC**

European Representative:

**3M Unitek GmbH  
inter-Unitek GmbH**

Postfach 15 62  
82171 Puchheim, Germany  
49-(0)89-800-810

5	4	3	2	1	1	2	3	4	5
5	4	3	2	1	1	2	3	4	5

**Caution:** U.S. Federal Law restricts this device to sale  
by or on the order of an orthodontic professional.

Transcend is a trademark of 3M Unitek.

© 1998 3M Unitek REF 11-461-2 9806