




Science.
Applied to Life.™
























3M™ Impregum™ Polyether Impression Material

Take control of challenging cases.

You need the precision standard when it comes to complex cases and moist environments. Implants and large restorations demand unfailing accuracy. The long working time, excellent flowability and hydrophilic properties of Impregum materials give you control over your toughest challenges.



An implant treatment can only be as accurate as the impression.

	Product	Dispensing System	Viscosity	Setting Version	Working Time* at 23°C/73°F	Total Setting Time**
			Low  High			
Tray Materials	3M™ Impregum™ Penta™ Soft Heavy Body			Regular Set	2:30	6:00
	3M™ Impregum™ Penta™ Soft Medium Body			Regular Set	2:45	6:00
	3M™ Impregum™ Penta™ Medium Body			Regular Set	2:45	6:00
	3M™ Impregum™ Soft Medium Body (Tray)			Regular Set	1:45	6:00
	3M™ Impregum™ Penta™ Soft Quick Step Heavy Body			Fast Set	1:00	4:00
	3M™ Impregum™ Penta™ Soft Quick Step Medium Body			Fast Set	1:00	4:00
	3M™ Impregum™ Soft Quick Step Medium Body (Tray)			Fast Set	1:00	4:00
Wash Materials	3M™ Impregum™ Penta™ Soft Light Body			Regular Set	3:15	6:30
	3M™ Impregum™ Soft Light Body			Regular Set	2:00	5:30
	3M™ Impregum™ Soft Quick Step Light Body			Fast Set	1:00	4:00

* Working time includes mixing time. ** includes working time.

Ordering information

Item #	Product Information
3M™ Penta™ Intro Kits	
P31739	3M™ Impregum™ Penta™ Soft Heavy Body/Light Body 1 base paste HB – 300 ml, 1 catalyst HB – 60 ml, 1 3M™ Penta™ Cartridge HB, 2 cartridges LB – 50 ml each, 1 3M™ Penta™ Elastomer Syringe, 10 3M™ Penta™ Mixing Tips – Red, 5 polyether mixing tips, white; 5 polyether intraoral tips, white; 5 3M™ Intra-oral Syringes Purple; 1 polyether tray adhesive – 17 ml, 1 cartridge dispenser 1:1/2:1
P31779	3M™ Impregum™ Penta™ Soft Quick Step Heavy Body/Light Body 1 base paste HB – 300 ml, 1 catalyst HB – 60 ml; 1 3M™ Penta™ Cartridge HB, 2 cartridges LB – 50 ml each, 10 3M™ Penta™ Mixing Tips – Red, 5 polyether mixing tips, white, 5 polyether intraoral tips, white; 5 3M™ Intra-oral Syringes Purple, 1 polyether tray adhesive – 17 ml
P31735	3M™ Impregum™ Penta™ Soft Medium Body 1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Penta™ Cartridge MB, 1 3M™ Penta™ Elastomer Syringe, 10 3M™ Penta™ Mixing Tips – Red, 1 polyether tray adhesive – 17 ml
P31771	3M™ Impregum™ Penta™ Soft Quick Step Medium Body 1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Penta™ Cartridge MB, 1 3M™ Penta™ Elastomer Syringe, 10 3M™ Penta™ Mixing Tips – Red, 1 polyether tray adhesive – 17 ml
P31685	3M™ Impregum™ Penta™ Medium Body 1 base paste MB – 300 ml, 1 catalyst MB – 60 ml, 1 3M™ Penta™ Cartridge MB, 1 3M™ Penta™ Elastomer Syringe, 10 3M™ Penta™ Mixing Tips – Red, 1 polyether tray adhesive – 17 ml
Cartridge Intro Kits	
31786	3M™ Impregum™ Soft Tray Medium Body Intro Kit 3 tray cartridges – 50 ml each, 1 cartridge LB – 50 ml, 10 polyether tray mixing tips, purple, 5 polyether mixing tips, white, 5 polyether intraoral tips, white, 5 3M™ Intra-oral Syringes Purple, 1 polyether tray adhesive – 17 ml
31789	3M™ Impregum™ Soft Quick Step Tray Medium Body Intro Kit 3 tray cartridges – 50 ml each, 1 cartridge LB – 50 ml, 10 polyether tray mixing tips, purple; 5 polyether mixing tips, white; 5 polyether intraoral tips, white, 5 3M™ Intra-oral Syringes Purple, 1 polyether tray adhesive – 17 ml

Item #	Product Information
3M™ Penta™ Refills	
31741	3M™ Impregum™ Penta™ Soft Heavy Body Refill 2 base pastes – 300 ml each, 2 catalysts – 60 ml
31769	3M™ Impregum™ Penta™ Soft Quick Step Heavy Body Refill 2 base pastes – 300 ml each, 2 catalysts – 60 ml each
31731	3M™ Impregum™ Penta™ Soft Medium Body Refill 2 base pastes – 300 ml each, 2 catalysts – 60 ml each
31765	3M™ Impregum™ Penta™ Soft Quick Step Medium Body Refill 2 base pastes – 300 ml each, 2 catalysts – 60 ml each
31642	3M™ Impregum™ Penta™ Medium Body Refill 2 base pastes – 300 ml each, 2 catalysts – 60 ml
Cartridge Refills	
31785	3M™ Impregum™ Soft Tray Medium Body Refill 4 cartridges – 50 ml each, 10 polyether tray mixing tips, purple
31788	3M™ Impregum™ Soft Quick Step Tray Medium Body Refill 4 cartridges – 50 ml each, 10 polyether tray mixing tips, purple
31752	3M™ Impregum™ Soft Light Body Refill 4 cartridges – 50 ml each, 10 polyether mixing tips, white, 5 3M™ Intra-oral Syringes Purple
31767	3M™ Impregum™ Soft Quick Step Light Body Refill 4 cartridges – 50 ml each, 10 polyether mixing tips, white, 5 3M™ Intra-oral syringes Purple
Accessories	
P3794	3M™ Impregum™ Penta™ Soft Heavy Body Cartridge
P3774	3M™ Impregum™ Penta™ Soft Quick Step Heavy Body Cartridge
P3788	3M™ Impregum™ Penta™ Soft Medium Body Cartridge
P3776	3M™ Impregum™ Penta™ Soft Quick Step Medium Body Cartridge
P3796	3M™ Impregum™ Penta™ Soft Light Body Cartridge
P3792	3M™ Impregum™ Penta™ Medium Body Cartridge
77944	3M™ Penta™ Cartridge for 3M™ Pentamix™ Lite
71508	50 3M™ Intra-Oral Syringes Purple
77580	1 cartridge dispenser 1:1/2:1
77949	50 3M™ Penta™ Mixing Tips – Red
71210	Elastomer Syringe Set Kit includes: 1 syringe; 4 tips; 4 barrels; 1 brush



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3M.com/dental

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada

Phone 1-888-363-3685

Customer Care Center: 1-800-634-2249
3M.com/greatimpressions

3M, Impregum, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada. Please recycle. Printed in U.S.A. © 3M 2016. All rights reserved. 70-2013-0155-6 (Rev. 8/2016)