


**3M** Science.  
Applied to Life.™

3M™ Paradigm™ VPS Impression Material

## An exceptional value – from the company that stands for quality and leads in impressioning








- Excellent detail reproduction
- Hydrophilic material
- Excellent flow of the wash material
- Optimal working and setting times
- Good tear strength for both regular and low viscosity washes
- Excellent dimensional stability

Available in regular and fast set formulations in a range of viscosities. Suitable for all impression indications.



# Optimal working and setting times

## 3M™ Paradigm™ VPS Impression Materials Overview Sheet

				Ordering Information		
	Color	Working Time (including mixing)	Setting Time (from start of mix)	3M™ Penta™ Refills: 2 – 300 ml base pastes, 2 – 60 ml catalysts	Cartridge Refills: 2 – 50 ml cartridges	Cartridge Jumbo Packs: 30 – 50 ml cartridges
<b>Tray Materials</b>						
3M™ Paradigm™ VPS Heavy Body		2:00	4:00	<b>36716</b>	<b>5313</b>	<b>5313J</b>
3M™ Paradigm™ VPS Heavy Body Fast Set		1:15	2:30	<b>36717</b>	<b>5513</b>	<b>5513J</b>
3M™ Paradigm™ VPS Monophase		2:00	4:00	—	<b>5312</b>	—
<b>Wash Materials</b>						
3M™ Paradigm™ VPS Regular Body		1:00	4:00	—	<b>5315</b>	—
3M™ Paradigm™ VPS Regular Body Fast Set		0:40	2:30	—	<b>5515</b>	—
3M™ Paradigm™ VPS Light Body		1:00	4:00	—	<b>5311</b>	—
3M™ Paradigm™ VPS Light Body Fast Set		0:40	2:30	—	<b>5511</b>	—
<b>Material Delivery</b>				<b>Item No.</b>		
3M™ Pentamix™ Lite Automatic Mixing Unit				<b>77907</b>		
3M™ Pentamix™ 3 Automatic Mixing Unit				<b>77876</b>		
Universal Cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit				<b>71545</b>		
3M™ Garant™ VPS Mixing Tips – Yellow				<b>71452</b>		
3M™ Garant™ VPS Intraoral Mix Tip – Yellow				<b>71462</b>		
3M™ Garant™ VPS Mixing Tips – Green				<b>71450</b>		