


3M Science.
Applied to Life.™

3M™ Paradigm™ VPS Impression Material

An exceptional value – from the company that stands for quality and leads in impressioning








- Excellent detail reproduction
- Hydrophilic material
- Excellent flow of the wash material
- Optimal working and setting times
- Good tear strength for both regular and low viscosity washes
- Excellent dimensional stability

Available in regular and fast set formulations in a range of viscosities. Suitable for all impression indications.



Optimal working and setting times

3M™ Paradigm™ VPS Impression Materials Overview Sheet

				Ordering Information		
	Color	Working Time (including mixing)	Setting Time (from start of mix)	3M™ Penta™ Refills: 2 – 300 ml base pastes, 2 – 60 ml catalysts	Cartridge Refills: 2 – 50 ml cartridges	Cartridge Jumbo Packs: 30 – 50 ml cartridges
Tray Materials						
3M™ Paradigm™ VPS Heavy Body		2:00	4:00	36716	5313	5313J
3M™ Paradigm™ VPS Heavy Body Fast Set		1:15	2:30	36717	5513	5513J
3M™ Paradigm™ VPS Monophase		2:00	4:00	—	5312	—
Wash Materials						
3M™ Paradigm™ VPS Regular Body		1:00	4:00	—	5315	—
3M™ Paradigm™ VPS Regular Body Fast Set		0:40	2:30	—	5515	—
3M™ Paradigm™ VPS Light Body		1:00	4:00	—	5311	—
3M™ Paradigm™ VPS Light Body Fast Set		0:40	2:30	—	5511	—
Material Delivery				Item No.		
3M™ Pentamix™ Lite Automatic Mixing Unit				77907		
3M™ Pentamix™ 3 Automatic Mixing Unit				77876		
Universal Cartridge for 3M™ Pentamix™ 3 Automatic Mixing Unit				71545		
3M™ Garant™ VPS Mixing Tips – Yellow				71452		
3M™ Garant™ VPS Intraoral Mix Tip – Yellow				71462		
3M™ Garant™ VPS Mixing Tips – Green				71450		