



3M™ Bair Hugger™ Therapy
Spinal Underbody Blanket
Model 575

ACHIEVE
SCIP-Inf-10
GOALS



Warming
Made Simple



3M Infection Prevention Solutions

Innovation
On A Mission

3M


3M™ Bair Hugger™ Therapy Spinal Underbody Blanket Model 575



The spinal underbody blanket is designed to work with the open frame of the spinal surgery cradle, and the Wilson™ frame.



3M™ Bair Hugger™ spinal underbody blanket shown with the footboard and ProneView™ protective helmet.



This advanced forced-air warming blanket does not interfere with the adjustment of the patient support pads and provides clinicians with full patient visualization and unrestricted patient access.

www.bairhugger.com



3M™ Bair Hugger™ spinal underbody blanket shown with the sling and ProneView™ protective helmet.



3M™ Bair Hugger™ Therapy Spinal Underbody Blanket Model 575

Features

The under-the-patient design of the spinal underbody blanket provides full, unrestricted patient access

- Unique fluid outlets minimize pooling of fluids on the surface of the blanket
- Integrated tie strips secure the blanket to the table frame
- Removable section allows visualization of the patient's face throughout the procedure
- Soft, radiolucent materials accommodate imaging requirements
- Perforations on the side of the blanket allow either the table frame or Wilson™ frame to pass through the blanket
- One drape (included) helps maintain warm air around the head of the intubated patient

Ordering Information

Product	Part Number	Specifications
Spinal Underbody Blanket	57501	Size: 91" x 36" (213 x 91 cm) Weight: 8 oz (218 g) Drape: 24" x 48" (61 x 122 cm) Blankets/Case: 5

To order, please contact your 3M Patient Warming representative or call Customer Service at 1-800-733-7775.

Product Offerings

Underbody Blankets

Pediatric Underbody Model 555	Large Pediatric Underbody Model 550	Adult Underbody Model 545	Spinal Underbody Model 575	Lithotomy Underbody Model 585	Full Access Underbody Model 635	Sterile Full Access Underbody Model 637

Preoperative Blankets

Outpatient w/Boots Model 110	Outpatient Model 111

Postoperative Blankets

Full Body Model 300	Chest Access Model 305	Multi-Access Model 315	Pediatric Model 310

Specialty and Cardiac Blankets

Cath Lab Model 560	Surgical Access Model 570	Sterile Cardiac Model 630	Sterile Cardiac Access Model 645

*also see underbody blankets

Intraoperative Blankets

Upper Body Model 522	XL Upper Body Model 523	Lower Body Model 525	Torso Model 540	Dual Port Torso Model 542	Full Body Surgical Model 610

Pediatric Blankets

Small Lower Body Model 537	Pediatric Underbody Long Model 530

*also see underbody blankets



**Arizant Healthcare Inc.,
a 3M company**
10393 West 70th St.
Eden Prairie, MN 55344 USA
Phone 800-733-7775
Fax 800-775-0002
www.bairhugger.com

3M is a trademark of 3M Company, used under license in Canada. BAIR HUGGER, RANGER and the BAIR HUGGER logo are trademarks of Arizant Healthcare Inc., used under license in Canada. Please recycle. Printed in USA.
©2011 Arizant Healthcare Inc. All rights reserved. ProneView is a trademark of DUPACO, Inc. Wilson is a trademark of Getinge USA, Inc. 601903D 5/11