

RelyX™ Veneer Cement

Light Curing composite cement for veneers




Easy to use and provides
excellent aesthetic results

Indications:

Permanent cementation of:

- Porcelain
- Ceramic and
- Composite veneers

Clinical Technique:



Photos courtesy of Dr. C. Sabrosa, Brazil

The simplicity you want and reliability you need.

RelyX™ Veneer Cement is a veneer cement system with a customer-preferred delivery and colour match system that's easy to use and easy to clean up.

RelyX Veneer Cement provides simple and efficient bonding of porcelain, composite and ceramic veneers. Combined with the 3M™ ESPE™ RelyX™ Try-In Pastes, the cement provides a simple and efficient way to perform veneer procedures.

- Optimal colour match of the water-soluble try-in pastes to the cured cement shade
- Excellent colour stability
- Ideal viscosity of the cement for easy insertion and positioning of the veneer



3M ESPE

Features & Benefits:

Customer-preferred delivery and shading system is easy to use.

Try-in pastes are water soluble and residue-free and offer excellent colour match to the final cured cement for better patient satisfaction.

Easy to clean up.

Ideal, non-slump cement viscosity prevents veneer drift and is easy to dispense, apply and seat.

RelyX™ Veneer Cement's excellent adhesion to dentine and enamel reduces the potential for microleakage and marginal staining.

Light-cure only, resin-based cement provides exceptional colour stability. Additionally, the cement cures with a wide range of light sources.

Available in six commonly known and preferred shades: Translucent, B0.5/White, White Opaque, A1/Light Yellow, A3 Opaque/Yellow Opaque and A5 Opaque/Dark.



Fast and easy bonding with Scotchbond™ Universal.

Scotchbond Universal is a single-bottle adhesive solution by 3M ESPE. It reliably bonds to all surfaces, including enamel, dentine, porcelain, composites and ceramic. Scotchbond Universal also functions as a silane primer for glass ceramic and composite veneers. This eliminates the need for a separate primer and more importantly reduces your inventory.

Ordering Information

Item No.	Product Information
	RelyX™ Veneer Cement
7616	RelyX™ Veneer Cement Intro Kit: 1 syringe cement Translucent 3g; 1 syringe cement White Opaque 3g; 1 syringe cement A1/Light Yellow 3g; 1 syringe cement A3 Opaque/Yellow Opaque 3g; 1 syringe cement B0.5/White 3g; 1 syringe try-in paste Translucent 2g; 1 syringe try-in paste White Opaque 2g; 1 syringe try-in paste A1/Light Yellow 2g; 1 syringe try-in paste A3 Opaque/Yellow Opaque 2g; 1 syringe try-in paste B0.5/White 2g; 1 Scotchbond™ Universal 5 ml vial; 1 Scotchbond™ Universal Etchant 3 ml syringe; 25 Etchant Tips; 48 White Mix Wells; 50 Scotchbond™ Microbrushes; Step-by-Step Technique Guide
7616TK	RelyX™ Veneer Cement Trial Kit: 1 syringe cement Translucent 3g; 1 syringe try-in paste Translucent 2g; 20 Scotchbond™ Universal L-Pop™; 1 Scotchbond™ Universal Etchant 3 ml syringe; 25 Etchant Tips; Veneer Cement Shade card; Step-by-Step Technique Guide
	Scotchbond™ Universal Adhesive
41258	Scotchbond™ Universal Adhesive Refill: 1 vial Scotchbond™ Universal à 5 ml



RelyX™ Veneer Cement Intro Kit (7616)

RelyX™ Veneer Cement Trial Kit (7616TK)

3M ESPE

3M Deutschland GmbH
Location Seefeld
3M ESPE · ESPE Platz
82229 Seefeld · Germany
info3mespe@mmm.com
www.3MESPE.com



3M, ESPE, L-Pop, RelyX, and Scotchbond are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada.

©3M 2011. All rights reserved.

XX (12.2011)