



RelyX™ Ultimate Adhesive Resin Cement

Ultimate strength. It's that simple.



Ideal for:
CAD/CAM and glass ceramic restorations

Compare the simplicity of the 3M adhesive cementation system to other leading dental cements:

Only two components

Industry-leading bond strength with fewer working steps

When you need extra-strength bonding without the hassle of complex adhesive systems, choose RelyX™ Ultimate Adhesive Resin Cement in combination with Scotchbond™ Universal Adhesive. You'll get versatile, reliable adhesive cementation that's fast and easy to use – all with just two components.

- Proven, industry-leading bond strength
- Fast and easy procedure with only two components
- Total-etch, selective etch or self-etch procedures
- First-class esthetics with natural fluorescence

RelyX™ Ultimate + Scotchbond™ Universal Adhesive = 2

Multilink® Automix + Adhesive Primer A, Adhesive Primer B, Resin Primer = 4

Panavia™ F 2.0 + Adhesive Primer A, Adhesive Primer B, Primer, Ceramic Primer = 6

Variolink® II with Syntac + Adhesive Primer, Adhesive Primer, Primer, Resin Primer = 6

NX3 with OptiBond® Solo Plus™/All-In-One + Adhesive Primer or All-in-One Adhesive = 3

When to use.

RelyX™ Ultimate Adhesive Resin Cement in combination with Scotchbond™ Universal Adhesive is indicated for final cementation of lower strength materials.

- CAD/CAM and glass ceramic restorations
- Resin-nano ceramics (e.g. Lava™ Ultimate CAD/CAM Restorative)
- All-ceramic, composite or metal inlays, onlays, crowns and bridges
- 2 - 3-unit Maryland bridges and 3-unit inlay/onlay bridges*
- All-ceramic or composite veneers
- All-ceramic, composite, or metal restorations to implant abutments
- Posts and screws

* Excluded for patients with bruxism or periodontitis



Introductory Kit



Trial Kit



3M Oral Care
2510 Conway Avenue
St. Paul, MN 55144-1000 USA

Phone 1-800-634-2249
Web 3m.com/dentalcements

3M Canada
Post Office Box 5757
London, Ontario N6A 4T1
Canada

Phone 1-888-363-3685

Recommended Applications (++ Highly recommended / + Recommended)

Indication	Metal/ Metal Based	Glass Ceramics (including Li. Disilicate) (e.g. e.max®, VITA Mark II, IPS Empress 2)	Oxide Ceramics (e.g. Lava™ Plus, Brux Zir®, Procera®)	Resin Nano Ceramics (e.g. Lava™ Ultimate CAD/CAM Restorative)	Resin Composites (e.g. VITA ENAMIC®)
Inlays/ Onlays	+	++	++	++	++
Crowns/ Bridges	+	++	+		++
Endodontic Posts	+		+		+
Veneers		+		+	+
Maryland Bridges	++		++		
On Implant Abutments	+	+	+		++

Ordering Information

Item No.	Product	Contents
56891	RelyX™ Ultimate Introductory Kit (incl. 4 color shades)	(1) 5ml bottle Scotchbond™ Universal Adhesive, (1) 3ml etchant syringe, (25) dispensing tips for etchant Syringe, (50) mixing wells, (50) microbrushes, (4) 8.5g syringes of cement - 1 in each shade: A3O, B0.5, A1 and Translucent, (30) mixing tips - regular, (10) mixing tips - wide, (10) endo tips, (10) intra-oral tips, (1) Technique Guide, (1) Instructions for Use
56892	RelyX™ Ultimate Trial Kit (Translucent)	(1) 1.5ml bottle Scotchbond™ Universal Adhesive, (1) 3ml etchant syringe, (25) dispensing tips for etchant Syringe, (50) microbrushes, (1) 8.5g syringe cement - Translucent, (10) mixing tips - regular, (5) mixing tips - wide, (5) intra-oral tips, (1) Technique Guide, (1) Instructions for Use
56894	RelyX™ Ultimate Trial Kit (A1)	(1) 1.5ml bottle Scotchbond™ Universal Adhesive, (1) 3ml etchant syringe, (25) dispensing tips for etchant syringe, (50) microbrushes, (1) 8.5g syringe cement - shade A1, (10) mixing tips - regular, (5) mixing tips - wide, (5) intra-oral tips, (1) Technique Guide, (1) Instructions for Use

RelyX™ Ultimate Adhesive Resin Cement Refills

Item No.	Shade	Contents
56887	A1	(1) 8.5g syringe cement
56888	A3O Opaque	(10) mixing tips - regular
56889	B0.5 Bleach	(5) mixing tips - wide
56890	Translucent	(5) intra-oral tips

Scotchbond™ Universal Adhesive

Item No.	Product	Size
41258	Scotchbond™ Universal Adhesive Refill	(1) 5ml vial
41263	Scotchbond™ Universal Etchant – Etching Gel	(2) 3ml syringes

3M, “3M Science. Applied to Life.”, Lava, RelyX and Scotchbond are trademarks of 3M Company or 3M Deutschland GmbH. Used under license in Canada. All other trademarks are owned by their respective owners.

Please recycle. © 3M 2016.
All rights reserved. 70-2013-0505-2 12596HB REV 05/2017