

3M™ Avagard™ Hand Hygiene Solutions Floor Stand 9244



The First Line of Defense In High-Traffic Areas.

Stand up to germs anywhere.

Now the brand that started the battle gives you a new way to gain the upper hand. A sleek, solidly constructed floor stand that carries 3M™ Avagard™ Foaming Instant Hand Antiseptic (70% v/v ethyl alcohol), 3M™ Avagard™ D Instant Hand Antiseptic with Moisturizers (61% w/w ethyl alcohol) or 3M™ Avagard™ Gel Instant Hand Antiseptic with Moisturizers (66% w/w ethyl alcohol). The perfect way to help fight HAIs in high-traffic zones.

- > No tools required to assemble stand*
- > Made of powder coated steel
- > Easy to clean drip tray
- > Self-righting base, up to 40 degrees past center (without dispenser)

*Tools required to mount dispenser.

Lead
The Way

3M

3M™ Avagard™ Hand Hygiene Solutions

Floor Stand 9244



AVAGARD™ Ordering Information

Catalog Number	Product Name	Size/Description	Item/Case
9244	3M™ Avagard™ Floor Stand	18 in. W x 54¼ in. H x 18 in. D (16.2 lbs.)	1
9240	3M™ Avagard™ Universal Hands Free Wall Dispenser	Use with 9322A, 9230 and 9338 1000 mL bottle	4
9322A	3M™ Avagard™ Foaming Instant Hand Antiseptic (70% v/v ethyl alcohol)	33.8 fl. oz. (1000 mL) (Dispenser bottle)	5
9230	3M™ Avagard™ D Instant Hand Antiseptic with Moisturizers (61% w/w ethyl alcohol)	33.8 fl. oz. (1000 mL) (Dispenser bottle)	5
9338	3M™ Avagard™ Gel Instant Hand Antiseptic with Moisturizers (66% w/w ethyl alcohol)	33.8 fl. oz. (1000 mL) (Dispenser bottle)	5

To learn more about 3M Avagard Hand Hygiene solutions, call the 3M Health Care Helpline at 1-800-228-3957 or visit www.3M.com/Avagard



3M Health Care
Infection Prevention Division
2510 Conway Avenue
St. Paul, MN 55144-1000
U.S.A.
1 800 228-3957
www.3m.com/infectionprevention

Please recycle. Printed in U.S.A.
3M and Avagard are trademarks of 3M.
© 3M 2014. All rights reserved.
70-2010-8459-0