


3M™ Bair Hugger™ Therapy
Surgical Access Blanket
Model 570

ACHIEVE
SCIP-Inf-10
GOALS



Warm from the Start

This unique blanket design accommodates supine, lateral, and prone positioning for surgical procedures such as spinal, abdominal, hip replacement and transplant surgery.





Features

- Large surgical access window
- Adaptable to supine, lateral, and prone positions
- Adhesive strips on patient side of surgical access window secure blanket position
- Non-inflated foot drape designed to protect feet from thermal injury
- Two resealable hose ports
- Attached head drape


Ordering Information

Product	Part Number	Specifications
Surgical Access Blanket Model 570	57000	Size: 84" x 36" (213 x 91 cm) Window Size: 13.5" x 22" (34 x 56 cm) Weight: 6 oz (181 g) Head Drape: 24" x 16" (61 x 41 cm) Units/Case: 10


To order, please contact your local 3M Patient Warming sales representative or call **Customer Service at 1-800-733-7775**.

Product Offerings


Preoperative Blankets




Postoperative Blankets




Specialty and Cardiac Blankets




Intraoperative Blankets



Pediatric Blankets



Underbody Series Blankets



**Arizant Healthcare Inc.,
a 3M company**
10393 West 70th St.
Eden Prairie, MN 55344 USA
Phone 800-733-7775
Fax 800-775-0002
www.bairhugger.com

3M is a trademark of 3M Company, used under license in Canada. BAIR HUGGER and the BAIR HUGGER logo are trademarks of Arizant Healthcare Inc., used under license in Canada. ©2012 Arizant Healthcare Inc. All rights reserved. 601215F 7/12