



3M™ Bair Hugger™ Therapy XL Upper Body Blanket Model 523



XL Upper Body Blanket
 Model 523



Upper Body Blanket
 Model 522

An extra large upper body blanket for larger patients

All the features you expect from a 3M™ Bair Hugger™ upper body blanket plus some changes specifically designed for larger patients including:

- 10" longer to accommodate a higher rise through the chest.
- 12" wider to maximize warming in the chest, shoulders and upper arms, but still narrow enough to provide surgical access to the abdomen.
- One additional air channel across the torso for a larger chest.
- Wider cut out at the chest to facilitate proper placement on a larger torso.


Ordering Information

XL Upper Body Blanket	Size	Weight	Drape
Model 52301	84" x 36" (213 x 91 cm)	5.4 oz (153.1 g)	24" x 24" (61 x 61 cm)


Product Offerings

Preoperative Blankets

Outpatient
w/Booties
Model 110




Outpatient
Model 111



Postoperative Blankets

Full Body
Model 300




Chest
Access
Model 305



Pediatric
Model 310



Multi-
Access
Model 315




Specialty and Cardiac Blankets


Cath Lab
Model 560



Surgical
Access
Model 570



Sterile
Cardiac
Model 630




Sterile Cardiac
Access
Model 645




Intraoperative Blankets

Upper Body
Model 522




XL Upper Body
Model 523




Lower Body
Model 525



Torso
Model 540



Dual Port
Torso
Model 542




Full Body
Surgical
Model 610



Pediatric Blankets

Pediatric
Long
Model 530



Small
Lower Body
Model 537



Underbody Series Blankets


Pediatric
Underbody
Model 555




Large Pediatric
Underbody
Model 550




Adult
Underbody
Model 545




Spinal
Underbody
Model 575




Lithotomy
Underbody
Model 585



Full Access
Underbody
Model 635



Sterile Full Access
Underbody
Model 637



To order, please contact your 3M Patient Warming representative or call Customer Service at 1-800-733-7775.



**Arizant Healthcare Inc.,
a 3M company**
10393 West 70th St.
Eden Prairie, MN 55344 USA
Phone 800-733-7775
Fax 800-775-0002
www.bairhugger.com

3M is a trademark of 3M Company, used under license in Canada. BAIR HUGGER and the BAIR HUGGER logo are trademarks of Arizant Healthcare Inc., used under license in Canada. ©2012 Arizant Healthcare Inc. All rights reserved. 601860E 8/12