



# 3M™ Bair Hugger™ Therapy XL Upper Body Blanket Model 523



**XL Upper Body Blanket**  
 Model 523



**Upper Body Blanket**  
 Model 522

## An extra large upper body blanket for larger patients

All the features you expect from a 3M™ Bair Hugger™ upper body blanket plus some changes specifically designed for larger patients including:

- 10" longer to accommodate a higher rise through the chest.
- 12" wider to maximize warming in the chest, shoulders and upper arms, but still narrow enough to provide surgical access to the abdomen.
- One additional air channel across the torso for a larger chest.
- Wider cut out at the chest to facilitate proper placement on a larger torso.


### Ordering Information

XL Upper Body Blanket	Size	Weight	Drape
Model 52301	84" x 36" (213 x 91 cm)	5.4 oz (153.1 g)	24" x 24" (61 x 61 cm)


## Product Offerings

### Preoperative Blankets

Outpatient  
w/Booties  
Model 110




Outpatient  
Model 111




### Postoperative Blankets


Full Body  
Model 300




Chest  
Access  
Model 305



Pediatric  
Model 310



Multi-  
Access  
Model 315




### Specialty and Cardiac Blankets


Cath Lab  
Model 560



Surgical  
Access  
Model 570



Sterile  
Cardiac  
Model 630




Sterile Cardiac  
Access  
Model 645




### Intraoperative Blankets

Upper Body  
Model 522




XL Upper Body  
Model 523




Lower Body  
Model 525



Torso  
Model 540



Dual Port  
Torso  
Model 542




Full Body  
Surgical  
Model 610



### Pediatric Blankets

Pediatric  
Long  
Model 530



Small  
Lower Body  
Model 537



### Underbody Series Blankets


Pediatric  
Underbody  
Model 555



Large Pediatric  
Underbody  
Model 550




Adult  
Underbody  
Model 545




Spinal  
Underbody  
Model 575




Lithotomy  
Underbody  
Model 585



Full Access  
Underbody  
Model 635



Sterile Full Access  
Underbody  
Model 637



To order, please contact your 3M Patient Warming representative or call Customer Service at 1-800-733-7775.



**Arizant Healthcare Inc.,  
a 3M company**  
10393 West 70th St.  
Eden Prairie, MN 55344 USA  
Phone 800-733-7775  
Fax 800-775-0002  
www.bairhugger.com

3M is a trademark of 3M Company, used under license in Canada. BAIR HUGGER and the BAIR HUGGER logo are trademarks of Arizant Healthcare Inc., used under license in Canada. ©2012 Arizant Healthcare Inc. All rights reserved. 601860E 8/12